



endolite
get busy living

Female Shifting Pyramid Adapter Part No. **189128**



Instructions for Use

1 Description and purpose

These instructions are for use by the Clinician/ Practitioner.

Application:

The Female Shifting Pyramid Adaptor is for use exclusively as part of a lower limb prosthesis. The device connects an open ended Endolite slotted chassis or top housing to a male pyramid receiver. This allows for shifting and angular & rotational adjustment between the two limb parts.

- Intended for single user only
- There are no known contra-indications if used in accordance with these instructions


Ensure that the user has understood any Instructions for use, drawing particular attention to the safety information.

Part No. 189128

2 Construction

Principal Parts:

Female Pyramid	Titanium
Collar	Aluminium Alloy
Pyramid Screws	Stainless Steel
Clamp Screw and inserts	Stainless Steel



3 Function

The Female Shifting Pyramid Adapter comprises two main components.

The movement between these two components prior to tightening the pyramid screws and clamp screw allows for alignment of the limb to suit the individual user.

When tightened securely using Loctite and torque settings shown, the device offers a secure way to connect two parts of the limb.

4 Maintenance

Maintenance is to be carried out by competent personnel only.

An annual visual inspection is recommended

- Ensure all screws are secure. If not remove and clean screws, reapply Loctite and tighten all screws to the correct torque settings.
- Check for visual defects that may affect proper function.

The wearer should be advised:

The user should be advised to contact their Clinician/Practitioner if their condition changes.

Any changes in performance of the device must be reported to the Clinician/Practitioner.

Changes in performance may include:

- Any unusual noises
- Excessive play or loss of alignment

5 Limitations on use

Intended Life:

The service life of the product is covered by the warranty period.

A local risk assessment should be carried out based upon activity and usage.

This product is recommended for use with other Endolite products.

Lifting Loads:

User weight and activity is governed by the stated limits.



Load carrying by the user should be based on a local risk assessment.

Environment:

Avoid exposing the product to corrosive elements such as water, acids and other liquids. Avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Exclusively for use between temperatures of -10 and 50°C

6 Assembly and alignment



Adjustment:

1. To adjust loosen clamp screw
2. Hold collar to slide/rotate as required
3. Tighten clamp screw to 8Nm torque.


7 Technical Specification:

Part No.	189128
Operating and Storage Temperature Range:	-10°C to 50°C 14°F to 122°F
Product Weight:	114g
Maximum User Weight:	150kg
Recommended Activity:	up to K4
Range of adjustment:	360 degrees rotation 7 degrees angular
AP Shift (dependent upon mating housing/chassis)	typically 15mm
Build Height:	28mm
Connections:	Distal Endolite Slotted Chassis or top housing Proximal Female Pyramid

Materials:

Titanium, Aluminium Alloy and Stainless Steel

Key Dimensions:



8 Warranty and Service:

Warranted for 3 years.

The user should be aware that changes or modifications not expressly approved could void the warranty.

9 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Blatchford Products Limited with sole responsibility according to appendix VII of the guidelines.

Distributor

Head Office

Chas A Blatchford & Sons Ltd
Lister Road
Basingstoke
Hampshire
RG22 4AH
United Kingdom
Tel: +44 (0) 1256 316600
Fax: +44 (0) 1256 316617
Email: sales@blatchford.co.uk
www.endolite.co.uk

Customer Services UK

Prosthetic and Orthotic Products
11 Atlas Way
Atlas North
Sheffield
S4 7QQ
United Kingdom
Tel: +44 (0) 114 263 7900
Fax: +44 (0) 114 263 7901
Email: sales@blatchford.co.uk
www.endolite.co.uk

endolite North America

1031 Byers Road
Miamisburg
Ohio 45342
USA
Tel: 800.548.3534
Fax: 800.929.3636
Email: info@endolite.com
www.endolite.com

endolite Germany

Endolite Deutschland GmbH
Holzstr. 5
95336 Mainleus
GERMANY
Tel: +49 9229 9737 001
Fax: +49 9229 9737 006
Email: info@endolite.de
www.endolite.de

endolite France

Parc d'Activités de l'Aéroport, 125 Impasse
Jean-Baptiste Say
34470 PEROLS
FRANCE
Tel: 00 33 (0) 467 820 820
Fax: 00 33 (0) 467 073 630
Email: contact@endolite.fr
www.endolite.fr

endolite India Ltd

A4 Naraina Industrial Area
Phase - 1
New Delhi
INDIA – 110028
Tel: 91 11 45689955
Fax: 91 11 25891543
Email: endolite@vsnl.com
www.endoliteindia.com

