Clinician's Manual

EN	Instructions for Use	page 2
FR	Instructions D'Utilisation	page 12
DE	Gebrauchsanweisung	seite 22
IT	Istruzioni per L'uso	pagina 32
ES	Instrucciones de Uso	pagina 42

Avalon^{K2}

endolite

1 Description and purpose

EN

Application:

These instructions are for use by the Clinician/Practitioner.

The Avalon^{K2} foot is to be used exclusively as part of a lower limb prosthesis. The Avalon^{K2} dorsi flexes after mid stance and remains dorsi flexed at toe off through the swing phase to give greater toe clearance for improved confidence and security. The controlled plantar flexion provided by the Avalon^{K2} at heel contact improves time to foot flat which enhances stability on uneven surfaces and slopes.

This device is recommended for users that have the potential to achieve Activity Level 2 who may benefit from enhanced stability and an increase in confidence on uneven surfaces and slopes.

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users in Activity Level 1 who would benefit from the enhanced stability offered by the Avalon^{K2}, but this decision should be made with sound and thorough justification.

Contra-indications:

Due to the cushioned heel and limited energy return from the foot keel this device is not suitable for Activity Level 3 or 4 individuals. This type of user will be better served by a specially designed prosthesis optimized for their needs.

This device is not suitable for users with poor balance.

Intended for single user.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs , stairs or uneven surfaces. Typical of the limited community ambulator.

Order Example:

AV	25L
----	-----

e.g. AV25L

Size
Side

Available from size 22 to size 30:

AV22L to AV30R


AV22LD to AV30RD

(add 'D' for a dark tone foot shell)

2 Construction

Principal parts:

- Hydraulic Body Assembly including pyramid (Aluminum/St. Stl./Titanium)
- Carrier Assembly (Aluminum/St. Stl.)
- Keel (Thermoplastic Composite)
- Keel attachment screws (St. Stl.)
- Glide Sock (UHM PE)
- Foot shell (PU)

3 Function

The Avalon^{K2} comprises a hydraulic body with a valve that can be adjusted to increase or reduce hydraulic resistance to plantar flexion and dorsi flexion simultaneously. The hydraulic body is connected to a carrier assembly via two pivot pins. A keel is attached to the carrier assembly using stainless steel screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

4 Maintenance

Maintenance must be carried out by competent personnel.

It is recommended that the following maintenance is carried out annually:

- Visually inspect the foot shell and glide sock, check for damage or wear and replace if necessary.

There are no other serviceable parts in the foot assembly.

The wearer should be advised:

Any changes in performance of this device must be reported to the practitioner

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise

Cleaning:

Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers.

5 Limitations on use:

Intended life:

A local risk assessment should be carried out based upon activity and usage if usage is beyond recommended activity level.

Lifting loads:

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment:

The Avalon^{K2} is water proof to a maximum depth of 1 metre.

Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.


Foot products must be adequately finished to prevent water ingress into the foot shell where possible. If water enters the foot shell, the limb should be inverted and dried before further use.

Exclusively for use between -15°C and 50°C.

It is recommended that only Endolite products be used in conjunction with the Avalon^{K2}.


6 Bench Alignment

Keep the build line between pivots as shown, using shift and/or tilt devices as necessary.

7 Biomimetic Alignment

The aim of alignment is to achieve a *balance point* while standing and set the hydraulically damped range of motion. The aim of damping adjustment is to fine tune the ankle-foot roll-over stiffness characteristics until a comfortable gait is achieved. Due to the increased range of motion provided by the ankle the user may experience the need for more voluntary control and initially find the ankle disconcerting during setup. This should quickly pass upon completion of satisfactory setup.


- * Ensure the user is not resting on the DF (Dorsiflexion) limit when standing in an upright position.

8 Biomimetic Adjustment

Note: Carry out static alignment while ensuring the user has some means of support such as parallel bars. This is standing alignment only.

Check the Bench Alignment and the heel height



The device should encourage some degree of self adjustment to achieve a sense of balance for the user during standing.

938327/5-0116

9 Dynamic Adjustment

The user should experience the ankle moving with the body through the gait cycle. There should be no effort exerted by the user to overcome the hydraulic resistance of the ankle.

Resistance adjustment procedure:


Guidance:
Following dynamic adjustment, trial the foot/ankle on ramps and stairs. Ensure the user is comfortable with the kind of terrain he/she may normally be expected to encounter. If the user reports any issues with comfort, usability or range of movement of the ankle, adjust accordingly. Caution: decreasing plantarflexion resistance may cause knee instability for Transfemoral users

10 Technical Data

Principal Materials	Aluminum, Stainless Steel, Titanium, PU, Thermoplastic Composite, UHM-PE
Component Weight [size 26]	780g (1lb 11ozs)
Maximum User Weight	150kg (330lb)
Recommended Activity Level	2
Size Range	22 to 30cm
Build Height [See diagram below]	115mm
Heel Height	10mm
Range of ankle movement	6 degrees plantarflexion to 6 degrees dorsiflexion
Operating and Storage Temperature Range	-15°C to 50°C (5°F to 122°F)
Proximal connection	Male Pyramid (Endolite)


Fitting Length

12 Assembly Instructions

1


Place foot in glide sock and slide assembly into foot shell.

The diagram shows a foot being inserted into a sock. An arrow points to the foot being inserted. Another arrow points to the sock being slid onto the foot. A third arrow points to the foot being inserted into a foot shell.

2

Ensure the heel of the keel is fully engaged in the location feature in the foot shell.

The diagram shows a cutaway section through the foot shell. An arrow points to the cutaway section. Another arrow points to the heel engagement.

3


Abraid the top surface of the foot shell. Bond 12mm Pelite interface to foot shell and shape as appropriate. Use Evostik 528 (926213) or equivalent.

The diagram shows a foot shell with a 12mm Pelite interface bonded to the top surface. An arrow points to the 12mm Pelite interface.

4

Bond foam cosmesis to top surface of the Pelite using Thixofix adhesive (926204), Evostik 528 (926213) or equivalent and shape to suit.

The diagram shows a foot shell with a foam cosmesis bonded to the top surface. An arrow points to the foam cosmesis.

13 Replacement Parts

Foot shell (for dark add 'D')			
Small	Medium	Large	Extra Large
22L 539020	24L 539024	26L 539028	28L 539032
22R 539021	24R 539025	26R 539029	28R 539033
23L 539022	25L 539026	27L 539030	29L 539034
23R 539023	25R 539027	27R 539031	29R 539035
			30L 539036
			30R 539037

Glide Sock, Sizes 22-30

Part No. 405815

Allen 2.5 (Hex key, valve adjuster tool)

Part No. 940234

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

Warranty:

The Avalon^{K2} foot is warranted for 2 years from original date of purchase unless otherwise stated.

Warranty for the Glide sock and Footshell is 12 months.

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Endolite Deutschland GmbH
endolite India Ltd