



SPECIFICATIONS

Part No.	Alloy	Adapter Style	Offset	Center Hole Size	Product Weight	Build Height	Max Patient Weight	Screw Torque (M8x12mm)
3T09	Titanium	Pyramid	None	10 mm.	46 gr.	12 mm.	136 Kg. / 300 lbs.	
3T10	Titanium	Pyramid	None	8.2 mm.	50 gr.	12 mm.	136 Kg. / 300 lbs.	
3T11	Titanium	Pyramid	Axial	8.2 mm.	60 gr.	12 mm.	136 Kg. / 300 lbs.	
3T12	Titanium	Pyramid	Diagonal	8.2 mm.	60 gr.	12 mm.	136 Kg. / 300 lbs.	
3T20	Titanium	Receiver	None	8.2 mm.	65 gr.	17 mm.	136 Kg. / 300 lbs.	15 Nm.
3S10	S. Steel	Pyramid	None	8.2 mm.	84 gr.	12 mm.	136 Kg. / 300 lbs.	
3S11	S. Steel	Pyramid	Axial	8.2 mm.	108 gr.	12 mm.	136 Kg. / 300 lbs.	
3S12	S. Steel	Pyramid	Diagonal	8.2 mm.	108 gr.	12 mm.	136 Kg. / 300 lbs.	
3S20	S. Steel	Receiver	None	8.2 mm.	102 gr.	17 mm.	136 Kg. / 300 lbs.	15 Nm.
3A10	Aluminum	Pyramid	None	8.2 mm.	51 gr.	12 mm.	100 Kg. / 220 lbs.	
3A20	Aluminum	Receiver	None	8.2 mm.	40 gr.	17 mm.	100 Kg. / 220 lbs.	15 Nm.

USAGE

This product is designed to be used for lower limb external prosthetic applications. Before fitting this product onto the patient, all parts must be joined properly and all set screws must be secured enough for the patient to walk. In order to make final adjustments, there should be no loose screws.

ATTENTION

After the final adjustments, apply single drop of removable Loctite® 242 Blue® to each of the set/cap screws and then torque the product. Improper installation or uneven, insufficient or excessive torque may cause bolt or adapter failure. Always use quality metric tools for installation.

WARRANTY

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.