



SPECIFICATIONS

Part No.	Alloy	Adapter Style	Product Height	Product Weight	Max Patient Weight	Torque (M8x12mm)
3T50	Titanium	Receiver/Receiver	32 mm	88 gr.	136 Kg. / 300 lbs.	15 Nm
3T51	Titanium	Receiver/Receiver	45 mm	98 gr.	136 Kg. / 300 lbs.	15 Nm
3T52	Titanium	Receiver/Receiver	60 mm	110 gr.	136 Kg. / 300 lbs.	15 Nm
3T53	Titanium	Receiver/Receiver	75 mm	120 gr.	136 Kg. / 300 lbs.	15 Nm
3T60	Titanium	Pyramid/Receiver	43 mm	76 gr.	136 Kg. / 300 lbs.	15 Nm
3T70	Titanium	Pyramid/Pyramid	35 mm	60 gr.	136 Kg. / 300 lbs.	N/A
3S50	S. Steel	Receiver/Receiver	32 mm	135 gr.	136 Kg. / 300 lbs.	15 Nm
3S51	S. Steel	Receiver/Receiver	45 mm	156 gr.	136 Kg. / 300 lbs.	15 Nm
3S52	S. Steel	Receiver/Receiver	60 mm	179 gr.	136 Kg. / 300 lbs.	15 Nm
3S53	S. Steel	Receiver/Receiver	75 mm	197 gr.	136 Kg. / 300 lbs.	15 Nm
3S60	S. Steel	Pyramid/Receiver	43 mm	124 gr.	136 Kg. / 300 lbs.	15 Nm
3S70	S. Steel	Pyramid/Pyramid	35 mm	109 gr.	136 Kg. / 300 lbs.	N/A
3A50	Aluminum	Receiver/Receiver	32 mm	65 gr.	100 Kg. / 220 lbs.	15 Nm
3A51	Aluminum	Receiver/Receiver	45 mm	80 gr.	100 Kg. / 220 lbs.	15 Nm
3A52	Aluminum	Receiver/Receiver	60 mm	94 gr.	100 Kg. / 220 lbs.	15 Nm
3A53	Aluminum	Receiver/Receiver	75 mm	112 gr.	100 Kg. / 220 lbs.	15 Nm

USAGE

This product is designed to be used for lower limb external prosthetic length and alignment corrections. Before fitting this product onto the patient, all parts must be joined properly and all set screws must be secured enough for the patient to walk. In order to make final adjustments, there should be no loose screws.

ATTENTION

After the final adjustments, apply single drop of removable Loctite® 242 Blue® to each of the set/cap screws and then torque the product. Improper installation or uneven, insufficient or excessive torque may cause bolt or adapter failure. Always use quality metric tools for installation.

WARRANTY

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.