EN	Instructions for Use	1
FR	Instructions D'Utilisation	20
DE	Gebrauchsanweisung	39
IT	Istruzioni per L'uso	59
ES	Instrucciones de Uso	79
NO	Bruksanvisning	98
RU	Инструкция протезиста	117
TR	Kullanım Talimatları	136

echelon  
VAC

endolite

# 1 Description and Purpose

EN

## Application:

These instructions are for use by the practitioner.

The EchelonVAC foot is to be used exclusively as part of a lower limb prosthesis.

The EchelonVAC provides limited self-alignment of the prosthesis on varied terrain and following footwear changes. It is intended to improve on postural sway and symmetry while easing abnormal pressures at the socket interface. A moderate energy return foot with multi-axial ankle movement. Independent heel and toe springs provides some axial deflection. The split toe provides good ground compliance.

In addition to the visco-elastic, self-aligning hydraulic ankle it generates an elevated vacuum in the range 12-17in Hg.

This device is recommended for users that have the potential to achieve Activity Level 3 who may benefit from enhanced stability and an increase in confidence on uneven surfaces.

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users in the Activity Levels 2 and 4\* who would benefit from the enhanced stability offered by the EchelonVAC, but this decision should be made with sound and thorough justification.

\*(maximum user weight 100 kg (220 lb) and always use one higher spring rate category than shown in the Spring Set Selection table).

 **To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.**

 **Be aware of finger trap hazard at all times.**

## Contra-indications:

This device may not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

It may not be suitable for use on individuals with poor balance, especially for bilateral use. If the user has any pertinent circulation condition seek medical advice if there is a possible risk of adverse reactions.

It is NOT recommended for use for:

- Wearers with poor cognitive function
- Users on dialysis
- Users with neuromas preventing weight bearing
- Use where a large range of heel height is required without re-alignment

EchelonVAC should only be fitted by suitably trained practitioners and should only be used with suitable, well fitting total contact sockets. There should be no reliefs or voids into which tissue may be drawn by the vacuum.

- If multiple walled sockets are used there should be no voids in their construction
- There should be no excessive flares to socket brim or trim lines

Intended for single user.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding Maintenance.

## Spring Set Selection

		User Weight								
Activity		44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125	kg
		(100-115)	(116-130)	(131-150)	(151-170)	(171-195)	(196-220)	(221-255)	(256-275)	(lbs)
3		1	2	3	4	5	6	7	8	Foot Spring Set

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

#### Note:

If in doubt choosing between two categories, choose the higher rate spring set.

Foot Spring set recommendations shown are for transtibial users.

For transfemoral users we suggest selecting a spring set one category lower, refer to fitting advice Section 7 to ensure satisfactory function and range of movement.

### Order Example:

<b>EVAC</b>	<b>25</b>	<b>L</b>	<b>N</b>	<b>3</b>	<b>S</b>
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Sandal Toe

Available from size 22 to size 30:

EVAC22L1S to EVAC30R8S

EVAC22L1SD to EVAC30R8SD

(add 'D' for a dark tone foot shell)

\*Sizes 25-27 only. For all other sizes, omit the Width field.

e.g. EVAC25LN3S, EVAC22R4S, EVAC27RW4SD


After continuous use the ankle casing may become hot to the touch.

# 2 Construction


## Principal Parts:

- Hydraulic Body Assembly including pyramid (Aluminium/St. Stl./Titanium)
- Carrier Assembly (Aluminium/St. Stl.)
- Heel & Toe Springs (e-Carbon)
- Spring Attachment Screws (Titanium/St.Stl.)
- Glide Sock (UHM PE)
- Foot Shell (PU)
- Vacuum Parts (PU, Nylon, Aluminium)


## Vacuum System Parts

- 1 Vac Barb 90deg
- 2a Auto-Expulsion Valve
- 2b Threaded Housing
- 3 Check Valve
- 4 In-line Filter
- 5 Vacuum Tubing

Vacuum Connector

---

## 3 Function

The EchelonVAC comprises a hydraulic body assembly containing adjustable hydraulic valves. The valves can be independently adjusted to increase and reduce hydraulic resistance of plantar flexion and dorsiflexion.

The hydraulic body also houses a pneumatic chamber and piston which, via one-way valves and a filter, creates a vacuum that can be passed via tubing to a prosthetic socket. For maximum vacuum effect the check valve should be positioned close to the ankle. The number of steps necessary to create an elevated vacuum will vary depending on the free space/air in the system. The use of multiple socks may require an increased number of steps to reach an elevated vacuum.

NB. if high hydraulic resistances are used such that they restrict ankle movement, the ability to generate vacuum may be compromised.

The hydraulic body assembly is connected to a carrier assembly via two pivot pins. Heel and toe springs are attached to the carrier assembly using titanium and stainless steel screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

---

## 4 Maintenance

Maintenance must be carried out by competent personnel.

It is recommended that the following maintenance is carried out annually:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use, this does not affect the function or strength of the foot.

The wearer should be advised:

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise
- Lack of vacuum

The practitioner must also be informed of:

- Any changes in body weight and/or activity level.
- Discoloration of the residual limb.

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

### **Cleaning:**

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.



## 4.1 Vacuum System Maintenance Guide

### 1 Visual Inspection


Make a visual inspection of the system parts paying attention to connections, these must be air tight to ensure integrity of the vacuum. Inspect the tubes and ensure they are firmly connected and are not kinked or split. The socket arrangement should also be inspected to check the integrity of the vacuum seals.

### 2 Check Valve

The check valve retains the vacuum created in the socket. It must be connected with the direction arrow pointing towards the EchelonVAC ankle.


## 4.2 Vacuum System Checklist


### 1 Cleaning the Check Valve, vacuum inlet

Disconnect the check valve and connect a syringe to the proximal tube with the flow direction arrow pointing away from the syringe. If the valve is working correctly the syringe should only push inwards. If the valve is blocked use the syringe to clear the valve with a 'blast of air' (Do not use compressed air). If it is still blocked clean it with distilled water using the syringe. If the valve is still not working replace it (409663 or 409863).

### 2 Cleaning the Check Valve, exhaust port

Check the exhaust valve is working correctly by connecting a syringe to the distal tube and clamp the proximal tube. Use a 'blast of air' to clear it through (Do not use compressed air). If the exhaust valve is working correctly and retaining the vacuum it should not be possible to draw the syringe plunger back out again.

---

## 5 Limitations on Use

For use only by appropriately trained practitioners.

Should only be used with well fitting total surface bearing sockets with no reliefs or voids, which have been constructed with air tight sockets and a suspension sleeve to create an air tight seal proximally.

### **Intended life:**

A local risk assessment should be carried out based upon activity and usage.

### **Lifting loads:**

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### **Environment:**

Avoid exposing the EchelonVAC to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Suitable for outdoor use


Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

It is recommended that only Endolite products be used in conjunction with the EchelonVAC.

# 6 Bench Alignment

## 6.1 Static Alignment


Keep the build line between pivots as shown, using shift and/or tilt devices as necessary.

Align transfemoral devices according to the instructions for use that are supplied with the knee.

### Tilt setting

Align limb to achieve range of motion shown.


For the Vacuum System Assembly see section 8.0  
Wrap the vacuum tubing around the pylon as shown and position the Check Valve close to the ankle for the best vacuum performance.

\* Allow for users own footwear

## 6.2 Biomimetic Alignment


The aim of alignment is to achieve a “balance point” while standing and set the hydraulically damped range of motion. The aim of damping adjustment is to fine tune the ankle-foot roll-over stiffness characteristics until a comfortable gait is achieved. Due to the increased range of motion provided by the ankle the user may experience the need for more voluntary control and initially find the ankle disconcerting during setup. This should quickly pass upon completion of satisfactory setup.

\* Ensure that the user is relaxed and not resting on the dorsiflexion limit.

## 6.3 Biomimetic Adjustment

NB: Carry out static alignment while ensuring the user has some means of support such as parallel bars. This is standing alignment only.

Use shift for static alignment and standing.


The device should encourage some degree of self adjustment to achieve a sense of balance for the user during standing.

## 6.4 Dynamic Adjustment

### Adjustment of the hydraulic valves.

The user should experience the ankle moving smoothly with the body through the gait cycle with no additional effort required by the user to overcome the hydraulic resistance of the ankle. During this procedure the user should walk at normal speed, in a straight line on a level surface.

#### Guidance:

Following dynamic adjustment, trial the foot/ankle on ramps and stairs. Ensure the user is comfortable with the kind of terrain he/she may normally be expected to encounter. If the user reports any issues with comfort, usability or range of movement of the ankle, adjust accordingly.

## 7 Fitting Advice

The correct alignment (A-P position), range of motion (distribution of plantar flexion to dorsiflexion) and adjustment of the hydraulic settings are critical in achieving a smooth roll over and correct slope adaptation (see 6.3).

The user should feel the vacuum effect after taking approximately 15-20 steps depending on the initial socket fit.

The springs for the EchelonVAC foot will be supplied assembled with heel and toe springs of the same category. If after following the instructions below you still have problems with the function please contact the sales team in your area for advice.

Any of the following:

- Incorrect spring selection
- Incorrect A-P shift alignment
- Incorrect distribution of plantar flexion and dorsiflexion range

will have a negative effect on function and stability.

	<i>Symptoms:</i>	<i>Remedy</i>
1.	<p>Sinking at heel strike</p> <p>Difficulty in achieving a smooth progression to mid stance</p> <p>User feels they are walking up hill or forefoot feels excessively long</p>	<ol style="list-style-type: none"> <li>1. Increase plantar flexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too anteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that the plantar flexion range is not excessive</li> <li>4. Check spring category is not too soft, if so fit a higher rate spring</li> </ol>
2.	<p>Progression from heel strike to mid stance is too rapid</p> <p>Difficulty in controlling the energy return from the foot at the heel strike (reduced knee stability)</p> <p>User feels heel is too hard, fore foot is too short</p>	<ol style="list-style-type: none"> <li>1. Reduce plantar flexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too posteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that there is adequate plantarflexion range</li> <li>4. Check the spring category is not too high for the weight and activity of the patient, if so fit lower rate spring</li> </ol>
3.	<p>Heel contact and progression feel OK but:</p> <p>Forefoot feels too soft</p> <p>Forefoot feels too short</p> <p>User feels they are walking down hill, possibly with reduced knee stability</p> <p>Lack of energy return</p>	<ol style="list-style-type: none"> <li>1. Increase dorsiflexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too posteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that there is not excessive dorsiflexion range</li> <li>4. Check the spring category is not too soft for the weight and activity of the patient, if so fit higher rate spring</li> </ol>

	<i>Symptoms:</i>	<i>Remedy</i>
4.	Forefoot feels too rigid Forefoot feels too long Feels like walking up hill	<ol style="list-style-type: none"> <li>1. Reduce dorsiflexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too anteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that there is sufficient dorsiflexion range</li> <li>4. Check the spring category is not too rigid for the weight and activity of the user, if so fit a lower rate spring</li> </ol>

### Vacuum System

	<i>Symptoms:</i>	<i>Cause/Remedy</i>
1.	Unable to generate a vacuum	Vacuum tube(s) split or disconnected Inspect and repair/replace as necessary
		Check and clean/replace the check valve
		Filter blocked, replace Filter
		Limited ankle movement creating insufficient vacuum due to: <ol style="list-style-type: none"> <li>1. Excessive PF/DF setting</li> <li>2. Footwear</li> </ol>
2.	Unable to maintain a vacuum	Vacuum tube(s) split or disconnected Inspect and repair/replace as necessary
		Check and clean/replace the check valve
		Leakage at socket valve/barbs Reseal valve/barb
		Porous socket Seal with lacquer/re-make
		Check the integrity of the vacuum seal at the socket/ residuum interface

# 8 Assembly Instructions


**Be aware of finger trap hazard at all times.**

## Footshell removal


Insert shoe horn behind Heel Spring.


Rotate shoe horn as shown to remove shell.


## Spring Replacement


Remove heel spring and screws.

Remove toe spring screw, replace toe. Upon reassembly, use Loctite 243 (926012) and torque to 35 Nm. Ensure toe spring is central to the carrier.

Re-assemble with replacement heel spring. Use Loctite 243 (926012) and torque to 15 Nm.

## 8 Assembly Instructions *(continued)*

**7** If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

**8**

928017

Lubricate toe and heel if required. (Foot shell is pre-lubricated).  
Fit sock as shown.

**9**

Slide carrier/heel spring assembly into the foot shell.

Toe spring location in foot shell.

**10**

Use a suitable lever to encourage the heel spring into location in the foot shell.

**11**

Ensure heel spring is engaged into slot.

heel spring location slot

**12**


Ensure glide sock does not get trapped when assembling to female pyramid part.

If a cosmetic finish is required please contact a member of the Endolite Sales Team.

## 8 Assembly Instructions *(continued)*

### Vacuum System Assembly


- 1 Push the In-line Filter on to a short piece of vacuum tubing and connect it to the Vac Barb 90deg on the socket. The barb outlet should be positioned pointing downwards by releasing the centre screw and rotating the body of the unit. Tighten the screw when the barb outlet is in the required position.
- 2 Attach a length of vacuum tubing to the In-line Filter and wrap it around the pylon. Connect the other end of the tube to the Check Valve ensuring that the flow arrow points towards the ankle. For maximum vacuum, position the Check Valve close to the inlet on the EchelonVAC. Connect a short length of vacuum tubing from the Check Valve to the inlet on the ankle to complete the vacuum system.

## 9 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight [Size 26N]:	930 g (2 lb 1 oz)
Recommended Activity Level:	2, 3, 4
Maximum User Weight:	125 kg (275 lb)
Proximal Alignment Attachment:	Male Pyramid (Endolite)
Range of Hydraulic Ankle Motion: (excludes additional range of motion provided by heel and toe springs)	6 degrees plantar flexion to 3 degrees dorsiflexion
Build Height:	[Sizes 22-24] 120 mm
[See diagram below]	[Sizes 25-26] 125 mm
	[Sizes 27-30] 130 mm
Heel Height:	10 mm
Maximum Vacuum:	17in Hg

### Fitting length

Size	A
22-24	120 mm
25-26	125 mm
27-30	130 mm

Size	B
22-26	65 mm
27-28	70 mm
29-30	75 mm

# 10 Ordering Information

Spring Kits				
Rate	Foot Sizes			
	Small (S)	Medium (M)	Large (L)	Extra Large (XL)
	22–24	25–26	27–28	29–30
Set 1	539801S	539810S	539819S	539828S
Set 2	539802S	539811S	539820S	539829S
Set 3	539803S	539812S	539821S	539830S
Set 4	539804S	539813S	539822S	539831S
Set 5	539805S	539814S	539823S	539832S
Set 6	539806S	539815S	539824S	539833S
Set 7	539807S	539816S	539825S	539834S
Set 8	539808S	539817S	539826S	539835S

Foot Shell (for dark add 'D')		
Size/Side	Narrow	Wide
22L	539038S	-
22R	539039S	-
23L	539040S	-
23R	539041S	-
24L	539042S	-
24R	539043S	-
25L	539044SN	539044SW
25R	539045SN	539045SW
26L	539046SN	539046SW
26R	539047SN	539047SW
27L	539048SN	539048SW
27R	539049SN	539049SW
28L	-	539050S
28R	-	539051S
29L	-	539052S
29R	-	539053S
30L	-	539054S
30R	-	539055S

Item	Part. No.
Glide Sock (Sizes 22–26)	531011
Glide Sock (Sizes 27–30)	532811
DF/PF Adjuster Key, 4.0 A/F Allen	940236
Vacuum System Parts:	
Socket Connection Kit	409663
Check Valve Service Kit	409863

## **Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

## **CE Conformity**

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

## **Warranty**

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months.

This warranty does not apply to:

Consumable parts including the vacuum tubing, in-line filter and valves unless a failure has occurred due to a defect in materials or workmanship.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See Endolite website for the current full warranty statement.

## **Environmental Aspects**

Where possible the components should be recycled in accordance with local waste handling regulations.

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