# Avalon<sup>K2-VAC</sup>

## Instructions for Use

AVAC24L-AVAC30R

|    |                            |     |
|----|----------------------------|-----|
| EN | Instructions for Use       | 3   |
| FR | Instructions D'Utilisation | 21  |
| DE | Gebrauchsanweisung         | 39  |
| IT | Istruzioni per L'Uso       | 57  |
| ES | Instrucciones de Uso       | 75  |
| NO | Bruksanvisning             | 93  |
| RU | Инструкция протезиста      | 111 |
| TR | Kullanım Talimatları       | 129 |

**Blatchford:**

# 1 Description and Purpose

EN

## Application

These instructions are for the practitioner.

The term device is used throughout these instructions for use to refer to AvalonK2VAC.

This device is to be used exclusively as part of a lower limb prosthesis.

It provides limited self-alignment of the prosthesis on varied terrain and following footwear changes. It is intended to improve on postural sway and symmetry while easing abnormal pressures at the socket interface.

This device dorsiflexes after mid-stance and remains dorsiflexed at toe-off through the swing phase to give greater toe clearance for improved confidence and security.

The controlled plantar flexion provided by this device at heel contact improves time to foot flat which enhances stability on uneven surfaces and slopes.

In addition to it functioning as a self-aligning hydraulic ankle it generates an elevated vacuum in the range 406-575 mbar/12-17 inHg.

## Activity Level

This device is recommended for users that have the potential to achieve Activity Level 2 who may benefit from enhanced stability and an increase in confidence on uneven surfaces and slopes.

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users in Activity Level 1 who would benefit from the enhanced stability offered by this device, but this decision should be made with sound and thorough justification.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs , stairs or uneven surfaces.  
Typical of the limited community ambulator.

## Contra-indications

Due to the cushioned heel and limited energy return from the foot keel this device is not suitable for Activity Level 3 or 4 individuals. This type of user will be better served by a specially designed prosthesis optimized for their needs.

It may not be suitable for use on individuals with poor balance , especially for bilateral use.

If the user has any pertinent circulation condition seek medical advice if there is a possible risk of adverse reactions.

It is NOT recommended for use for:

- Wearers with poor cognitive function
- Users on dialysis
- Users with neuromas preventing weight bearing
- Use where a large range of heel height is required without re-alignment

The device should only be fitted by suitably trained practitioners and should only be used with suitable, well fitting total contact sockets. There should be no reliefs or voids into which tissue may be drawn by the vacuum.

- If multiple walled sockets are used there should be no voids in their construction
- There should be no excessive flares to socket brim or trim lines

## 1.1 Safety Information

**The Warning symbol highlights safety information which must be followed carefully.**

**Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.**

**After continuous use the ankle casing may become hot to the touch.**

**To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.**


**Be aware of finger trap hazard at all times.**

Intended for a single user.

## 2 Construction


### Principal Parts:

- Hydraulic Body Assembly including pyramid (aluminum/stainless steel/titanium)
- Carrier Assembly (aluminum/stainless steel)
- Keel (thermoplastic composite)
- Keel Attachment Screws (stainless steel)
- Glide Sock (UHM PE)
- Foot Shell (PU)
- Vacuum Parts (PU, nylon, aluminum)


### Vacuum System Parts

- 1a Auto-Expulsion Valve
- 1b Threaded Housing
- 2 Check Valve
- 3 In-line Filter
- 4 Vacuum Tubing

### Vacuum Connector

---

## 3 Function

The device comprises a hydraulic body assembly containing adjustable hydraulic valves. The valves can be independently adjusted to increase and reduce hydraulic resistance of plantar flexion and dorsiflexion.

The hydraulic body also houses a pneumatic chamber and piston which, via one-way valves and a filter, creates a vacuum that can be passed via tubing to a prosthetic socket. For maximum vacuum effect the check valve should be positioned close to the ankle. The number of steps necessary to create an elevated vacuum will vary depending on the free space/air in the system. The use of multiple socks may require an increased number of steps to reach an elevated vacuum.

Note... If high hydraulic resistances are used such that they restrict ankle movement, the ability to generate vacuum may be compromised.

The hydraulic body is connected to a carrier assembly via two pivot pins. A keel is attached to the carrier assembly using stainless steel screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

---

## 4 Maintenance

Maintenance must be carried out by competent personnel.

We recommend the following annual maintenance:

- Visually inspect the foot shell and glide sock, check for damage or wear and replace if necessary
- Carry out vacuum system maintenance according to Section 4.1 and 4.2

There are no other serviceable parts in the foot assembly.

The wearer should be advised that any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise
- Lack of vacuum

The practitioner must also be informed of:

- Any changes in body weight and/or activity level
- Discoloration of the residual limb

### **Cleaning**

Use a damp cloth and mild soap to clean outside surfaces.

Do not use aggressive cleansers.



## 4.1 Vacuum System Maintenance Guide

### 1. Visual Inspection


Make a visual inspection of the system parts paying attention to connections, these must be air tight to ensure integrity of the vacuum. Inspect the tubes and ensure they are firmly connected and are not kinked or split. The socket arrangement should also be inspected to check the integrity of the vacuum seals.

### 2. Check Valve

The check valve retains the vacuum created in the socket. It must be connected with the direction arrow pointing towards the device.


## 4.2 Vacuum System Checklist


### 1. Cleaning the Check Valve, vacuum inlet

Disconnect the check valve and connect a syringe to the proximal tube with the flow direction arrow pointing away from the syringe. If the valve is working correctly the syringe should only push inwards. If the valve is blocked use the syringe to clear the valve with a 'blast of air' (Do not use compressed air). If it is still blocked clean it with distilled water using the syringe. If the valve is still not working replace it (409663 or 409863).

### 2. Cleaning the Check Valve, exhaust port

Check the exhaust valve is working correctly by connecting a syringe to the distal tube and clamp the proximal tube. Use a 'blast of air' to clear it through (Do not use compressed air). If the exhaust valve is working correctly and retaining the vacuum it should not be possible to draw the syringe plunger back out again.

---

## 5 Limitations on Use

Should only be used with well fitting total surface bearing sockets with no reliefs or voids, which have been constructed with air tight sockets and a suspension sleeve to create an air tight seal proximally.

### **Intended Life:**

A local risk assessment should be carried out based upon activity and usage.

### **Lifting Loads:**

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### **Environment:**

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Suitable for outdoor use

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

It is recommended that only Blatchford products be used in conjunction with this device.

# 6 Bench Alignment

## 6.1 Static Alignment



Keep the build line between pivots as shown, using shift and/or tilt devices as necessary.

**Align transfemoral devices according to the instructions for use that are supplied with the knee.**


### Tilt Setting

Align limb to achieve range of motion shown.

For the Vacuum System Assembly see Section 8.

Wrap the vacuum tubing around the pylon as shown and position the Check Valve close to the ankle for the best vacuum performance.


\* Allow for user's own footwear

## 6.2 Biomimetic Alignment

The aim of alignment is to achieve a *balance* point while standing and set the hydraulically damped range of motion. The aim of damping adjustment is to fine tune the ankle-foot roll-over stiffness characteristics until a comfortable gait is achieved.

Due to the range of motion provided by the ankle the user may experience the need for more voluntary control and initially find the ankle disconcerting during setup. This should quickly pass upon completion of satisfactory setup.


- \* Ensure the user is not resting on the DF (Dorsiflexion) limit when standing in an upright position.

## 6.3 Biomimetic Adjustment

Carry out static alignment while ensuring the user has some means of support such as parallel bars. This is standing alignment only.


Use shift for static alignment and standing.

The device should encourage some degree of self adjustment to achieve a sense of balance for the user during standing.

## 6.4 Dynamic Adjustment

### Adjustment of the hydraulic valves

The user should experience the ankle moving smoothly with the body through the gait cycle with no additional effort required by the user to overcome the hydraulic resistance of the ankle. During this procedure the user should walk at normal speed, in a straight line on a level surface.

### Guidance

Following dynamic adjustment, trial the foot/ankle on ramps and stairs. Ensure the user is comfortable with the kind of terrain he/she may normally be expected to encounter. If the user reports any issues with comfort, usability or range of movement of the ankle, adjust accordingly.

## 7 Fitting Advice

The correct alignment (A-P position), range of motion (distribution of plantar flexion to dorsiflexion) and adjustment of the hydraulic settings are critical in achieving a smooth roll over and correct slope adaptation (see 6.3).

The user should feel the vacuum effect after taking approximately 15-20 steps depending on the initial socket fit.

The following conditions will have a negative effect on the function and stability of the device:

- Incorrect A-P shift alignment
- Incorrect distribution of plantar flexion and dorsiflexion range

### Ankle Assembly

|    | Symptoms  | Remedy   |
|----|---|--|
| 1. | Sinking at heel strike<br>Difficulty in achieving a smooth progression to mid stance<br>User feels they are walking up hill or forefoot feels excessively long                  | <ol style="list-style-type: none"> <li>1. Increase plantar flexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too anteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that the plantar flexion range is not excessive</li> </ol> |
| 2. | Progression from heel strike to mid stance is too rapid<br>User feels heel is too hard, fore foot is too short, or reduced knee stability.                                      | <ol style="list-style-type: none"> <li>1. Reduce plantar flexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too posteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that there is adequate plantar flexion range</li> </ol>     |
| 3. | Heel contact and progression feel OK but:<br>Forefoot feels too soft<br>Forefoot feels too short<br>User feels they are walking down hill, possibly with reduced knee stability | <ol style="list-style-type: none"> <li>1. Increase dorsiflexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too posteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that there is not excessive dorsiflexion range</li> </ol>    |
| 4. | Forefoot feels too rigid<br>Forefoot feels too long<br>Feels like walking up hill   | <ol style="list-style-type: none"> <li>1. Reduce dorsiflexion resistance</li> <li>2. Check A-P shift alignment; ensure foot is not too anteriorly positioned</li> <li>3. Check distribution of plantar flexion and dorsiflexion movement; ensure that there is sufficient dorsiflexion range</li> </ol>          |

## Vacuum System

|   | Symptoms                    | Cause/Remedy   |
|---|-----------------------------|--|
| 1.  | Unable to generate a vacuum | Vacuum tube(s) split or disconnected   |
|   |                             | Inspect and repair/replace as necessary  |
|   |                             | Check and clean/replace the check valve  |
|   |                             | Filter blocked, replace Filter   |
|   |                             | Limited ankle movement creating insufficient vacuum due to:<br>1. Excessive PF/DF setting<br>2. Footwear |
| 2.  | Unable to maintain a vacuum | Vacuum tube(s) split or disconnected   |
|   |                             | Inspect and repair/replace as necessary  |
|   |                             | Check and clean/replace the check valve  |
|   |                             | Leakage at socket valve/barbs  |
|   |                             | Reseal valve/barb  |
|   |                             | Porous socket  |
|   |                             | Seal with lacquer/re-make  |
| Check the integrity of the vacuum seal at the socket/residuum interface |                             |  |

# 8 Assembly Instructions

## Hydraulic Ankle and Foot Shell Assembly


Be aware of finger trap hazard at all times.

- 1** Place foot in glide sock and slide assembly into foot shell.


- 2** Ensure the heel of the keel is fully engaged in the location feature in the foot shell.

- 3** Abraid the top surface of the foot shell. Bond 12 mm Pelite interface to foot shell and shape as appropriate. Use Evostik 528 (926213) or equivalent.

- 4** Bond foam cosmesis to top surface of the Pelite using Thixofix adhesive (926204), Evostik 528 (926213) or equivalent and shape to suit.


## 8 Assembly Instructions (continued)

### Vacuum System Assembly


1. Push the In-line Filter on to a short piece of vacuum tubing and connect it to the Auto-Expulsion Valve.
2. Attach a length of vacuum tubing to the In-line Filter and wrap it around the pylon. Connect the other end of the tube to the Check Valve ensuring that the flow arrow points towards the ankle. For maximum vacuum, position the Check Valve close to the inlet on the device. Connect a short length of vacuum tubing from the Check Valve to the inlet on the ankle to complete the vacuum system.

## 9 Technical Data

|   |   |
|---|---|
| Operating and Storage Temperature Range | -15°C to 50°C (5°F to 122°F)                        |
| Component Weight (size 26)              | 583 g (1 lb 5 oz)                                   |
| Maximum User Weight                     | 150 kg (330 lb)                                     |
| Activity Level                          | 2   |
| Size Range                              | 24 to 30 cm   |
| Build Height<br>(See diagram below)     | 122 mm  |
| Heel Height                             | 10 mm   |
| Range of Hydraulic Ankle Motion         | 6 degrees plantar flexion to 3 degrees dorsiflexion |
| Proximal connection                     | Male Pyramid (Blatchford)                           |
| Maximum Vacuum:                         | 575 mbar (17 inHg)                                  |

### Fitting Length

# 10 Ordering Information

## Order Example

|             |           |               |
|-------------|-----------|---------------|
| <b>AVAC</b> | <b>25</b> | <b>L</b>      |
|             | Size      | Side<br>(L/R) |

Available from size 24 to size 30:

AVAC24L to AVAC30R

AVAC24LD to AVAC30RD

*(add 'D' for a dark tone foot shell)*

e.g. AVAC25L

## Foot Shell

| Medium     | Large      | Extra Large |
|------------|------------|-------------|
| 24L 539024 | 26L 539028 | 28L 539032  |
| 24R 539025 | 26R 539029 | 28R 539033  |
| 25L 539026 | 27L 539030 | 29L 539034  |
| 25R 539027 | 27R 539031 | 29R 539035  |
|            |            | 30L 539036  |
|            |            | 30R 539037  |

(For dark add 'D' to Part No.)

## Glide Sock

| Sizes: | Part. No. |
|--------|-----------|
| 24-30  | 405815    |

## Other Parts

| Item                              | Part. No. |
|-----------------------------------|-----------|
| DF/PF Adjuster Key, 4.0 A/F Allen | 940236    |
| <b>Vacuum System Parts</b>        |           |
| Socket Connection Kit             | 409663    |
| Check Valve Service Kit           | 409863    |

## **Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

## **CE Conformity**

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

## **Warranty**

This device is warranted for 24 months—glide sock 3 months—foot shell 12 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

This warranty does not apply to:

Consumable parts including the vacuum tubing, in-line filter and valves unless a failure has occurred due to a defect in materials or workmanship.

## **Environmental Aspects**

Where possible the components should be recycled in accordance with local waste handling regulations.

## **Manufacturer's Registered Address**

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© Blatchford Products Limited 2019. All rights reserved.
Endolite India Ltd