Epirus

Instructions for Use

EP22L1S-EP30R8S
EP22L1SD-EP30R8SD

EN	Instructions for Use	3
FR	Instructions D'Utilisation	13
DE	Gebrauchsanweisung	23
IT	Istruzioni per L'Uso	34
ES	Instrucciones de Uso	45

1 Description and purpose

Application

These instructions are for the practitioner.

The term *device* is used throughout these instructions for use to refer to Epirus.

This device is to be used exclusively as part of a lower limb prosthesis.

A moderate-energy-return foot with resilient multi-axial ankle movement. Independent heel and toe spring provides some axial deflection. The split toe provides good ground compliance.

This device is recommended for users that have the potential to achieve Activity Level 3.

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users in the Activity Levels 2 and 4* who would benefit from the smoother transition from heel to toe offered by this device, but this decision should be made with sound and thorough justification.

*Maximum user weight 100 kg and always use one higher spring rate category than shown in the Spring set selection table.

To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

Contra-indications

This device might not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Intended for single user.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

This is a guide for initial spring & plantar flexion buffer selection, refer to Section 7 *Fitting Advice* for final choice.

Activity	User weight								kg (lb)
	44-52 (100-115)	53-59 (116-130)	60-68 (131-150)	69-77 (151-170)	78-88 (171-195)	89-100 (196-220)	101-116 (221-255)	117-125 (256-275)	
3	1	2	3	4	5	6	7	8	Foot Spring set
	None		Medium			Firm			<i>Suggested Plantar Flexion Buffer</i>

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Note:

If in doubt choosing between two categories, choose the higher rate spring set.


Foot Spring set recommendations shown are for transtibial users.

For transfemoral users or lower activity user we suggest selecting a spring set one category lower, refer to fitting advice Section 7 to ensure satisfactory function and range of movement.

2 Construction

Principal parts:

- Carrier Assembly (aluminum/St. Stl.)
- Buffer (natural rubber)
- Heel & Toe Springs (e-carbon)
- Spring Attachment Screws (titanium/ St. Stl.)
- Glide Sock (UHM PE)
- Foot Shell (PU)

3 Function

This device comprises an e-carbon toe and independent heel spring combined with resilient multi-axial ankle movement, ankle function can be customised to the individuals requirements via the use of interchangeable buffers.

Heel and toe springs are attached to the carrier using titanium and stainless steel screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

4 Maintenance

Maintenance must be carried out by competent personnel.

It is recommended that the following maintenance is carried out annually:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use, this does not affect the function or strength of the foot.
- Remove the dome, check the condition of the plantar flexion buffer and replace if necessary

The user should be advised:

Any changes in performance of this device must be reported to the practitioner,

changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noises

The practitioner must also be informed of any changes in body weight and/or activity level.

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

5 Limitations on Use

Intended life

A local risk assessment should be carried out based upon activity and usage.

Lifting loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing this device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.


Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for outdoor use

We recommend using only Blatchford products with this device.

6 Bench Alignment

Static alignment

Setup length

With flexion, adduction and abduction properly accommodated, set the limb length 5 mm longer than the sound side to allow for compression and deflection of the foot springs during gait. This should be re-assessed once the dynamic trial commences and the length adjusted accordingly.

Build line

This should fall between the centre line of the pyramid and 10 mm anterior (with heel height properly accommodated). The socket should be positioned accordingly.

Dynamic alignment

Coronal plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal Plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are **evenly loaded** and that both are **touching** the floor.

* Allow for users own footwear

7 Fitting Advice

Spring sets are supplied as matched pairs i.e. the heel and toe spring are designed to work together. Progression throughout the stance phase should be smooth; the combined heel and ankle function is key to this process.

If a smooth progression is not achieved, the following should be considered:

	Symptoms	Remedy
Ankle too soft	<ul style="list-style-type: none"> • Rapid plantar flexion and possible foot slap • Lack of energy return from heel spring 	Fit plantar flexion buffer. A choice of medium and firm buffer is available (the medium version is castellated), see figs 1 & 2 below. For trial purposes the buffer does not need to be bonded. Final selection should be bonded in position using Loctite 424.
Heel spring too soft	<ul style="list-style-type: none"> • Sinking at heel strike • Difficulty climbing over the toe • (toe feels too hard) 	<ol style="list-style-type: none"> 1. Move socket anteriorly in relation to the foot (excess movement may result in drop-off) 2. Fit a stiffer spring set
Heel spring too hard	<ul style="list-style-type: none"> • Rapid transition from heel strike through stance phase • Difficulty in controlling heel action, foot jars into mid-stance • Foot feels too rigid 	<ol style="list-style-type: none"> 1. Move socket posteriorly in relation to foot 2. Fit softer spring set

Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

fig 1. Showing Dome removed ready for fitting buffer

fig 2. Medium Buffer fitted

8 Foot Shell Removal


Be aware of finger trap hazard at all times.

- 1** Insert shoe horn behind heel spring.


- 2** Rotate shoe horn as shown to remove shell.

9 Assembly Instructions

- 1** Assemble toe spring onto carrier, use Loctite 243 (926012) and torque to 35 Nm.

- 2** Assemble heel spring onto carrier, use Loctite 243 (926012) and torque to 15 Nm.

- 3** Cover appropriate lines on carrier with permanent black marker to leave spring set number showing.

- 4**
- 
- 928017**
Lubricate toe and heel if required.
(Foot shell is pre-lubricated)


Fit glide sock over the carrier assembly as shown.

9 Assembly Instructions *(continued)*

5

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

6

Slide carrier/heel spring assembly into the foot shell.

7

Toe spring location in foot shell.

8

Use a suitable lever to encourage the heel spring into the location in the foot shell.

9

Ensure heel spring is engaged into slot.

Heel spring location slot

10

Ensure glide sock does not get trapped when assembling to female pyramid part.


10 Disassembly Instructions

Dis-assembly is the reverse of assembly.

11 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (Size 26N):	610 g (1 lb 6 oz)
Activity Level:	2–3
Maximum User Weight:	125 kg (275 lb)
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Adjustment	+/- 7° angular
Build Height: (See diagram below)	85 mm
Heel Height	10 mm

Fitting length

Size	A
22–26	65 mm
27–28	70 mm
29–30	75 mm

12 Ordering Information

Order example

EP	25	L	N	3	S
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Sandal Toe

Available from size 22 to size 30:

EP22L1S to EP30R8S

EP22L1SD to EP30R8SD

(add 'D' for a dark tone foot shell)

*Sizes 25-27 only. For all other sizes, omit the Width field.

e.g. EP25LN3S, EP22R4S, EP27RW4SD

Spring Kits				
Rate	Foot Sizes			
	Small (S)	Medium (M)	Large (L)	Extra Large (XL)
	22-24	25-26	27-28	29-30
Set 1	539801S	539810S	539819S	539828S
Set 2	539802S	539811S	539820S	539829S
Set 3	539803S	539812S	539821S	539830S
Set 4	539804S	539813S	539822S	539831S
Set 5	539805S	539814S	539823S	539832S
Set 6	539806S	539815S	539824S	539833S
Set 7	539807S	539816S	539825S	539834S
Set 8	539808S	539817S	539826S	539835S

Foot Shell (for dark add 'D')		
Size/Side	Narrow	Wide
22L	539038S	-
22R	539039S	-
23L	539040S	-
23R	539041S	-
24L	539042S	-
24R	539043S	-
25L	539044SN	539044SW
25R	539045SN	539045SW
26L	539046SN	539046SW
26R	539047SN	539047SW
27L	539048SN	539048SW
27R	539049SN	539049SW
28L	-	539050S
28R	-	539051S
29L	-	539052S
29R	-	539053S
30L	-	539054S
30R	-	539055S

Item	Part. No.
Glide Sock (Sizes 22-26)	531011
Glide Sock (Sizes 27-30)	532811
Plantar flexion buffer: Medium	533710
Plantar flexion buffer: Firm	533711

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

Warranty

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Manufacturer's Registered Address

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Blatchford SAS
Endolite India Ltd