Navigator

Instructions for Use

NAV22LWW-NAV30RBB
NAV22LWWD-NAV30RBBD

EN	Instructions for Use	2
FR	Instructions D'Utilisation	12
DE	Gebrauchsanweisung	22
IT	Istruzioni per L'Uso	32
ES	Instrucciones de Uso	42
NO	Bruksanvisning	52
RU	Инструкция протезиста	62
TR	Kullanım Talimatları	72

1 Description and Purpose

Application

These instructions are for the practitioner.

The term *device* is used throughout these instructions for use to refer to Navigator.

This device is an ankle-foot unit that is for use as a component of a lower limb prosthesis.

Make sure that the user understands all the information in these instructions for use, particularly the maintenance and safety information.

This device comprises an integrated multi-axial ankle and a resilient keel which offer the user enhanced ground compliance.

Activity Level

This device is recommended for Activity Level 1–3 users.

As an exception to our recommendation, Activity Level 4* users may only use this device based upon the outcome of a prosthetic assessment.

*Maximum user weight: 100 kg (220 lb).

Contra-indications

This device is not suitable for use in competitive sports events.

This device is for a single user.

Ankle Ball and Snubber Selection


The ankle ball and the snubber are color coded to indicate their hardness.

Use the table to choose an ankle ball and snubber combination that suits the weight and Activity Level of the user. The letter on the left refers to the ankle ball. For example, WO refers to a white ankle ball and an orange snubber.

- White (W) = Soft
- Orange (O) = Medium
- Blue (B) = Hard

Notes:

- Bilateral users might need stiffer snubbers than unilateral users.
- Trans-femoral users might need softer ankle balls than trans-tibial users.

Location of color codes

Weight		Activity Level		
kg	lb	1	2	3
44–52	100–115	WW	WW	WW
53–59	116–130	WW	WW	WO
60–68	131–150	WW	WO	OO
69–77	151–170	WO	OO	OB
78–88	171–195	OO	OB	BB
89–100	196–220	OB	BB	BB
101–125	221–275	BB	BB	BB

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.


Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as, curbs, stairs or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

2 Construction

Principal parts:

- Dome (aluminum alloy)
- Snubber (rubber)
- Top Housing (aluminum alloy)
- Ankle Ball and Pyramid (aluminum alloy/ rubber/stainless steel/titanium)
- Keel (fiber-reinforced polymer)
- M6x27 Socket Cap Screw (plated carbon steel)
- Glide Sock (UHM PE)
- Foot Shell (polyurethane)

3 Function

This device is an ankle-foot unit that comprises both an integrated ankle and a resilient keel inside a sandal-toe foot shell. The ankle center is positioned anatomically to promote a natural gait from heel strike to toe-off. The plantar flexion and dorsiflexion resistances can be adjusted by changing the ankle ball and the snubber to suit the weight and the Activity Level of the user.

4 Maintenance

Only competent personnel should carry out maintenance.

Carry out the following maintenance annually:

- Check the foot shell and glide sock for damage or wear, and replace if necessary;
- Check the ankle ball and the snubber for damage or wear, and replace if necessary;
- Make sure the socket cap screws are tight. If they are loose, tighten them to the correct torque setting;
- Relubricate the snubber.

Advise the user to report the following to the practitioner:

- Changes in body weight or activity level
- Changes in the performance of this device, for example, unusual noises or play.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces. **Do not** use aggressive cleansers.

5 Limitations on Use

Intended life

A local risk assessment should be carried out based upon activity and usage.

Lifting loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Avoid abrasive environments such as those containing sand, as these may promote premature wear.


Suitable for outdoor use

Exclusively for use between -15 °C and 50 °C (5 °F and 122 °F).

We recommend using Blatchford products with this device.

6 Alignment

Additional components

Follow the alignment instructions that are supplied with any additional components.

Static alignment

Sagittal plane

With the flexion accommodated, the load line should pass through the center of the socket and through the center of the pyramid.

Make sure that, when the user is standing, the heel and the toe are evenly loaded and the foot is in full contact with the floor.

Dynamic alignment

Coronal plane

Minimise M-L thrust by adjusting the relative positions of the socket and the foot .

Sagittal plane

Make sure that the user transitions smoothly from heel strike to toe-off.

* Accommodate the heel height of the user's footwear


7 Fitting Advice

During dynamic alignment the user can experience gait abnormalities. Use this table to correct some common gait abnormalities.


Gait Abnormality	Action
<ul style="list-style-type: none">Sinking at heel strikeFoot slap	<ol style="list-style-type: none">Confirm A-P alignment and heel height accommodation.Fit a harder ankle ball.
<ul style="list-style-type: none">Rapid transition from heel strike through stance phaseDifficulty controlling heel action (e.g. the foot jars into mid-stance)Foot feels too rigid	<ol style="list-style-type: none">Confirm A-P alignment and heel height accommodation.Fit a softer ankle ball.
<ul style="list-style-type: none">Drop-offFeels like walking downhill	<ol style="list-style-type: none">Confirm A-P alignment and heel height accommodation.Fit a harder snubber.
<ul style="list-style-type: none">Feels like walking uphill	<ol style="list-style-type: none">Confirm A-P alignment and heel height accommodation.Fit a softer snubber.

8 Changing the Ankle Ball

The main function of the ankle ball is to provide plantar flexion resistance. Harder ankle balls provide more resistance and softer ankle balls provide less resistance.


- Clean the ankle ball and the keel pocket with acetone or a similar solvent.
- Align the flats of the male pyramid along the axis of the keel and rest the ankle ball in the keel pocket.
- Insert the alignment screw into the ankle ball through the hole in the bottom of the keel.

4. Clean inside of the top housing with acetone or a similar solvent.
5. Place the top housing over the ankle ball.

Note: Make sure the opening of the top housing faces posteriorly.

6. Apply pressure to the top housing while maintaining pyramid alignment. Then, engage the socket cap screws into the housing.


7. Turn each screw, in sequence, ½ to 1 turn at a time until the top housing comes into contact with keel.

Note: Make sure that the alignment screw stays in place while you tighten the socket cap screws.

8. Torque tighten the socket cap screws to 10 Nm and remove the alignment screw.

9 Changing the Snubber


The main function of the snubber is to provide dorsiflexion resistance. Harder snubbers provide more resistance and softer snubbers provide less resistance.

1. Lubricate the snubber with a suitable rubber lubricant.
2. Press the cone tool onto the housing.
3. Place the dome over the ankle stem.

Note: Make sure the 'A' on top of the dome faces anteriorly.

10 Assembling and Finishing the Foot Shell

1. Place the keel inside the glide sock and slide them both into the foot shell.

Note: Make sure that the heel of the keel is fully engaged in the location feature in the heel of the foot shell.

Note: After assembling to the limb, make sure that the glide sock is not trapped under the dome.

2. Pull the glide sock under slight tension.
3. Use a cable tie to secure the glide sock to the pylon just above the adapter.

4. Remove the end of the cable tie, then fold the glide sock over.


5. Abrade the top surface of the foot shell.
6. Bond a 12 mm pelite interface to the foot shell and shape as appropriate. Use Evostik 528 (926213) or an equivalent.
7. Bond a foam cosmesis to the top surface of the pelite interface and shape as appropriate. Use Thixofix (926204), Evostik 528 (926213) or an equivalent.

11 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (size 26):	565 g (20 oz)
Activity Level:	1-3
Maximum User Weight:	125 kg (275 lb)
Proximal Alignment Attachment:	Male pyramid (Blatchford)
Range of Adjustment	±7 ° angular
Build Height:	95 mm
Heel Height:	10 mm

Fitting length

12 Ordering Information

Order example

This device is available in sizes 22–30.

Note: For dark tone foot shells, add a D to the end of the order number.

NAV 25L W O D NAV25LWOD

Size & Side Ankle Ball Snubber Optional: Dark-Tone Foot Shell

Ankle and Keel Assemblies (USA only)

Size	Type	Part Number
Small	White 60°	409156
Medium	White 60°	409157
Large	White 60°	409158
Extra Large	White 60°	409159
Small	Orange 70°	409256
Medium	Orange 70°	409257
Large	Orange 70°	409258
Extra Large	Orange 70°	409259
Small	Blue 80°	409356
Medium	Blue 80°	409357
Large	Blue 80°	409358
Extra Large	Blue 80°	409359

Snubbers

Type	Part Number
White 60°	405605
Orange 70°	405606
Blue 80°	405607

Glide Sock 405815

Foot Shells (Add a D for dark tone)

Size (cm)	Left Part Number	Right Part Number
22	539020	539021
23	539022	539023
24	539024	539025
25	539026	539027
26	539028	539029
27	539030	539031
28	539032	539033
29	539034	539035
30	539036	539037

Tools for Assembly

Part	Part Number
Alignment Screw	910414
Cone Tool	940083
Cable Tie	929029

Ankle Ball and Pyramids

Type	Part Number
White 60°	409456
Orange 70°	409556
Blue 80°	409656

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CE conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class I Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

Warranty

This device is warranted for 24 months, the foot shell is warranted for 12 months, and the glide sock is warranted for 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

Manufacturer's Registered Address

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Blatchford SAS
Endolite India Ltd