# Silcare Breathe Locking Liner

## Instructions for Use

SBTTLP22L-SBTTLP40L

EN	Instructions for Use	2
FR	Instructions D'Utilisation	10
DE	Gebrauchsanweisung	18
IT	Istruzioni per L'Uso	26
ES	Instrucciones de Uso	34
NO	Bruksanvisning	42
RU	Инструкция протезиста	50
TR	Kullanım Talimatları	58

**Blatchford:**

# 1 Description and Purpose

EN

## Application

These instructions are for the practitioner. The term *device* is used throughout these instructions for use to refer to Silcare Breathe Locking Liner.

This device is an interface component for use only as a component of a lower limb prosthesis.

Please give users the user guide and make sure that they understand all the information in these instructions for use, particularly the maintenance and safety information.

This device is a biocompatible, perforated locking liner that allows air and moisture to escape through its perforations during stance. During the gait cycle, a silicone valve provides vacuum suspension by controlling the airflow through the perforations on the distal cap.

## Activity Level

Although this device is suitable for users of all Activity Levels, it is best suited to users who have the potential to achieve Activity Level 3–4.

## Contra-indications

Deep scars on the distal end of residual limbs can cause failure of the vacuum suspension. These deep scars reduce skin contact with the silicone inside the liner and make a less effective seal.

Conical residual limbs can cause failure of the vacuum suspension. Conical residual limbs reduce skin contact with the silicone inside the liner and make a less effective seal.

Short residual limbs can cause failure of the silicone side of the liner. After donning, if the unperforated fabric is on or above the patella tendon, the silicone is subject to more stress during knee flexion.

# 2 Safety Information

 This warning symbol highlights important safety information.

 Unwanted material can stick to silicone and irritate skin.

 The fabric can collect static electricity.

 Users with sensitive skin, diabetes, and problems of the vascular system must examine the condition of their residual limb and apply lubricant to sensitive areas.

 To avoid danger of suffocation, keep the device away from babies and children.

 Keep the device away from direct heat sources.

 Enlarged perforations can trap the skin and cause blisters. If the perforations enlarge, stop using the device and contact a Blatchford sales representative.

 **Do not** overtighten the locking pin.

 If the residual limb has tissue damage, stop using the device.

 **Do not** pull or stretch the fabric.

⚠ Fingernails, sharp jewelry, and the locking pin can tear the fabric. If the fabric is torn, stop using it and contact a Blatchford sales representative.

⚠ Only use the device in combination with corrosion resistant components.

⚠ Sockets with sharp proximal edges can tear the device.


⚠ **Do not** use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

---

## 3 Construction

Principal parts:

- Perforated Fabric (polyamide and lycra)
- Unperforated Fabric (polyamide and lycra)
- Valve (silicone)
- Washer (nylon)
- Spacer\*\* (silicone)
- Retaining Screw\* (nylon)
- Casting Dummy (silicone)

\* For transportation only. **Do not** use for fitting.

\*\* For use with some pin locks. (See [section 9](#).)

---

## 4 Function

During stance, the valve opens to allow moisture and air to escape through the perforations in the distal cap. During swing, the valve closes and stops air from re-entering the device, thus creating a vacuum, which both improves proprioception and reduces pistoning.

---

## 5 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric
- Tissue damage on the residual limb
- Enlarged perforations
- Changes in either body weight or activity level
- Changes in the performance of the device, for example:
  - Moisture in the distal end of the device. (See [section 12.1.](#))
  - Loss of vacuum. (See [section 12.2.](#))

### Notes:

- Perspiration discolors some socket materials.

### 5.1 Cleaning the device

Suitable for machine washing at 30 °C.

Clean daily.

Dry thoroughly before use.

**Do not** tumble dry.

 Unwanted material can stick to silicone and irritate skin.

 The fabric can collect static electricity.

 **Do not** use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

1. Invert the liner so that its silicone side points externally.
2. Clean the silicone, the valve, and the distal cap with a solution of water and unperfumed, pH balanced soap.
3. Invert the liner so that its silicone side points internally.
4. Fill the liner with warm water while both holding the proximal end closed and squeezing the distal end to flush the water through the perforations.
5. Rinse the device with clean water to remove all residues.
6. Either dry the device with a lint-free cloth, or leave it to air dry.

**Note:** Always dry the device with the silicone side of the liner pointing internally. Otherwise, the fabric can stretch and become distorted.

## 5.2 Cleaning the residual limb


 Users with sensitive skin, diabetes, and problems of the vascular system must examine the condition of their residual limb and apply lubricant to sensitive areas.

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Apply lotion to dry skin, as recommended by a practitioner.

---

## 6 Choosing the Correct Size

1. Measure the circumference of the residual limb 4 cm up from its distal end.
2. If the circumference of the residual limb matches one of the sizes in the table, choose the next size down.
3. If the circumference of the residual limb is between one of the sizes in the table, choose the smaller size.

		Size (cm)									
		22	23.5	25	26.5	28	30	32	34	36	40
Valve Size		Small			Medium		Large			Extra Large	

---

## 7 Casting the Socket with the Casting Dummy

### Before you begin

Allow the user to wear the device for 10 minutes. (See [section 11.](#))

938413


---

## 8 Fitting the Valve

 Do not overtighten the locking pin.

We recommend using a locking pin that has a shoulder of between 13.5 mm and 19 mm in diameter.


1. Put the valve on the distal cap.
2. Put the washer between the locking pin and the valve.
3. Apply Loctite 243 to the threads of the locking pin.
4. Torque the locking pin to 3 Nm, or tighten it with your fingers by either  $\frac{1}{8}$  to  $\frac{1}{4}$  of a turn.

---

## 9 Fitting the Spacer

Fit the spacer between the pin lock and the valve only if the pin lock prevents the valve from opening. (See [section 12](#).)

---

## 10 Trimming the Liner

### Before you begin

We recommend using a rotating disc cutter to achieve a bevelled edge.

Leave enough material on the liner so that its proximal edge can cover the edge of the socket.

Avoid trimming through the perforations.

- Trim the proximal edge of the liner to suit the shape of the user's socket.

If necessary, trim the proximal edge of the liner in a wave pattern to reduce shear stress.

---

## 11 Donning the Device

⚠ Do not pull or stretch the fabric.

⚠ Fingernails, sharp jewelry, and the locking pin can tear the device. Stop using the device if it is torn or damaged.

⚠ Unwanted material can stick to silicone and irritate skin.

⚠ The fabric can collect static electricity.

### Before you begin


Fit the valve. (See [section 8](#).)

1. Invert the liner so that its silicone side points externally.
2. Put the distal cap on the distal end of the residual limb.
3. Align the locking pin with the long axis of the residual limb.
4. Roll the liner on the residual limb while releasing any trapped air.

**Note:** The unperforated fabric should be below the patella tendon.

After donning the device, allow the user to wear it for 10 minutes. If the user feels numbness within 10 minutes, do the following:

1. Do off the device, and wait until normal sensation returns.
2. Don the device again.
3. If the user feels numbness again, stop using the device.

---

## 12 Fitting Advice

### 12.1 Moisture in the distal end of the device

If moisture collects inside the distal end of the device, refer to the following table:

Cause	Solution
The perforations are clogged.	Clean the device. (See <a href="#">section 5.1.</a> )
The socket does not have enough space in its distal end for the valve to open.	Cast another socket. (See <a href="#">section 7.</a> )
The pin lock is stopping the valve from opening.	Put the spacer between the pin lock and the valve. (See <a href="#">section 9.</a> )

### 12.2 Loss of vacuum

If the device loses vacuum, refer to the following table:

Cause	Solution
The valve is damaged.	Stop using the device, and contact a Blatchford sales representative.
The locking pin is not aligned with the long axis of the residual limb.	Re-align the locking pin with the long axis of the residual limb.

### 12.3 Separation of the valve from the distal cap

If the valve falls off the distal cap, refer to the following table:


Cause	Solution
The valve is not fitted correctly.	Make sure that the washer is between the pin and the valve. (See <a href="#">section 8.</a> )

---

# 13 Ordering Information

## Order Example

Part	Part Number
Small Valve Kit	559015
Medium Valve Kit	559016
Large Valve Kit	559017
Extra Large Valve Kit	559018
Casting Dummy	559019

\* SBTTLP23L is for size 23.5. SBTTLP26L is for size 26.5.

---

### Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

### CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

### Warranty

This device warranted for 6 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions.

See the Blatchford website for the current full warranty statement.

### Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH, UK.

Silcare Breathe and Blatchford are registered trademarks of Blatchford Products Limited.

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© Blatchford Products Limited 2019. All rights reserved.
Endolite India Ltd