Silcare Breathe Locking Liner

User Guide

SBTTLP22L–SBTTLP40L

EN	User Guide	2
FR	Guide de l'Utilisateur	9
DE	Bedienungsanleitung für Anwender	16
IT	Guida per il Paziente	23
ES	Guia de Usuario	30
NO	Brukerhåndbok	37
RU	Руководство пользователя	44
TR	Kullanım Kılavuzu	51
AR	دليل المستخدم	58

Blatchford:

1 Description and Purpose

EN

Application

These instructions are for the user.

The term *device* is used throughout these instructions for use to refer to Silcare Breathe Locking Liner.


This device is an interface component for use only as a component of a lower limb prosthesis.

This device is a biocompatible, perforated locking liner that allows air and moisture to escape through its perforations. During the gait cycle, a silicone valve provides vacuum suspension by controlling airflow through the perforations on the distal cap.

2 Safety Information

-  This warning symbol highlights important safety information.
-  Unwanted material can stick to silicone and irritate skin.
-  The fabric can collect static electricity.
-  Users with sensitive skin, diabetes, and problems of the vascular system must examine the condition of their residual limb and apply lubricant to sensitive areas.
-  To avoid danger of suffocation, keep the device away from babies and children.
-  Keep the device away from direct heat sources.
-  Enlarged perforations can trap the skin and cause blisters. If the perforations enlarge, stop using the device and contact your practitioner.
-  **Do not** overtighten the locking pin.
-  If the residual limb has tissue damage, stop using the device, and contact your practitioner.
-  **Do not** pull or stretch the fabric.
-  Fingernails, sharp jewelry, and the locking pin can tear the device. Stop using the device if it is torn or damaged.
-  Only use the device in combination with corrosion resisitant components.
-  **Do not** use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

3 Construction

* For use with some locks. (See [section Z.](#))

4 Function

During stance, the valve opens to allow moisture and air to escape through the perforations in the distal cap. During swing, the valve closes and stops air from re-entering the device, thus creating a vacuum, which both improves proprioception and reduces pistoning.

5 Maintenance

Report the following to your practitioner:

- Tears in the fabric
- Tissue damage on the residual limb
- Enlarged perforations of the perforated fabric
- Changes in either body weight or activity level
- Changes in the performance of the device, for example:
 - Moisture in the distal end of the device. (See [section 9.1.](#))
 - Loss of vacuum. (See [section 9.2.](#))

Notes:

- Perspiration discolours some socket materials.

5.1 Cleaning the device

Suitable for machine washing at 30 °C.

Clean daily.

Dry thoroughly before use.

Do not tumble dry.

Unwanted material can stick to silicone and irritate skin.

Do not use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

1. Invert the liner so that its silicone side points externally.
2. Clean the silicone, the valve, and the distal cap with a solution of water and unperfumed, pH balanced soap.
3. Invert the liner so that its silicone side points internally.
4. Fill the liner with warm water while both holding the proximal end closed and squeezing the distal end to flush the water through the perforations.
5. Rinse the device with clean water to remove all residues.
6. Either dry the device with a lint-free cloth, or leave it to air dry.

Note: Always dry the device with the silicone side of the liner pointing internally. Otherwise, the fabric can stretch and become distorted.

5.2 Cleaning the residual limb

 Users with sensitive skin, diabetes, and problems of the vascular system must examine the condition of their residual limb and apply lubricant to sensitive areas.


1. Clean the skin daily with a solution of water and unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Apply lotion to dry skin as recommended by your practitioner.

6 Fitting the Valve

 Do not overtighten the locking pin.


We recommend using a locking pin that has a shoulder of between 13.5 mm and 19 mm in diameter.

1. Put the valve on the distal cap.
2. Put the washer between the locking pin and the valve.
3. Apply Loctite 243 to the threads of the locking pin.
4. Torque the locking pin to 3 Nm, or tighten it with your fingers by either $\frac{1}{8}$ or $\frac{1}{4}$ of a turn.

7 Fitting the Spacer

Fit the spacer between the lock and the valve only if the lock prevents the valve from opening. (See [section 9.](#))


8 Donning the Device

- ⚠ Do not pull or stretch the fabric.
- ⚠ Fingernails, sharp jewelry, and the locking pin can tear the device. Stop using the device if it is torn or damaged.
- ⚠ Unwanted material can stick to silicone and irritate skin.
- ⚠ The fabric can collect static electricity.

Before you begin

Fit the valve. (See [section 6.](#))

1. Invert the liner so that its silicone side points externally.
2. Put the distal cap on the distal end of the residual limb.
3. Align the locking pin with the long axis of the residual limb.
4. Roll the liner on the residual limb while releasing any trapped air.

Note: The distal end of the perforated fabric should be just below the patella tendon.

After donning the device, wear it for 10 minutes. If you feels numbness within 10 minutes, do the following:

1. Doff the device, and wait until normal feeling returns.
2. Don the device again.
3. If you feel numbness again, stop using the device.

9 Fitting Advice

9.1 Moisture in the distal end of the device

If moisture collects inside the distal end of the device, refer to the following table:

Cause	Solution
The perforations are clogged.	Clean the device. (See section 5.1.)
The lock is stopping the valve from opening.	Put the spacer between the lock and the valve. (See section 7.)

If moisture continues to collect inside the device, contact your practitioner.

9.2 Loss of vacuum

If the device loses vacuum, refer to the following table:

Cause	Solution
The valve is damaged.	Stop using the device, and contact your practitioner.
The locking pin is not aligned with the long axis of the residual limb.	Re-align the locking pin with the long axis of the residual limb.

If the device continues to lose vacuum, contact your practitioner.

12.3 Separation of the valve from the distal cap

If the valve falls of the distal cap, refer to the following table:

Cause	Solution
The valve is not fitted correctly.	Make sure that the washer is between the pin and the valve. (See section 6.)

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

Warranty

This device warranted for 6 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions.

See the Blatchford website for the current full warranty statement.

Manufacturer's Registered Address

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Blatchford SAS
Endolite India Ltd