



THERAPEUTIC SHOES

HOMECARE

MAY 2016



Dear Customers,

We are happy to present you our therapeutic shoes in a completely new light.

The catalogue at hand contains all FIOR & GENTZ therapeutic shoes especially designed for a stay in a retirement home, (rehab) clinic or for simply wearing them at home.

The combination of protection and high comfort is what distinguishes our HOMECARE shoes. All shoes have flexible soles and are well suited for indoors as well as outdoors.

All the shoes' specific details can be found on each shoe's particular product page next to the order information. You can find information on the measurement methods, inner shoe dimensions and available insoles from page 16 on.

Our goal is to provide you with products enabling you to treat your patients the best possible way. You will find further therapeutic shoes in our **product catalogues POST-OP and DIABETES**.

With best regards from Lüneburg
Jörg Fior and Ralf Gentz

All our products meet our high standard of quality. This quality standard has been controlled and certified by an independent certification organisation according to the international standards ISO 9001 and ISO 13485.



SALES TEAM IN-HOUSE

In order to provide you with up to date contact details of our distributors in foreign countries, we have listed all our current distributors online at: www.fior-gentz.com/contact/export

Ana María von Corvin
English/Spanish



Michael Krüger
English



Anna-Lena Kühn
English/French



SALES TEAM BENELUX

Luc Schols
Dutch/English/French



IN-HOUSE SALES TEAM

For any questions, orders or customer data changes just call +49 4131 24445-0. We are available on workdays from 8.00am to 5.00pm (Central European Time).

If you place your order by 4.15pm, we will dispatch the goods on that same day. Express deliveries on Saturday are also possible.

FIT FOR DAILY (REHAB) ROUTINE

Indoors & outdoors

All HOMECARE shoes are particularly well suited for indoor as well as outdoor activities.

Breathable and washable

All shoes provide a breathable upper and can be washed at 40°C.



Fastener system

Due to hook and loop fasteners, all HOMECARE shoes are easy to open and close. While wearing the shoe, large loop areas ensure individual comfort. Besides, the hook and loop fasteners have additional seams so that the fabric can be neatly cut, if necessary.

Adaptable inner volume

Due to the hook and loop fasteners, which can be neatly cut, the inner volume can be adapted to the patient's individual need for space.



Outsoles

The flexible soles are 30mm thick, made of EVA material and have a slightly wedged heel. This sole design supports a natural heel to toe movement while walking.

Insoles

Inside the shoe you will find three insoles made of EVA material. The top insole has an antimicrobial coating. Due to reworkable insoles, the inner volume and the hardness of the cushioning can be adapted to individual needs for space and wear. The insoles can be replaced by individual foot orthotics, if applicable.

Order information on page 17

Reinforced heel cap

All HOMECARE shoes have a reinforced heel cap which provides support, cushions the heel and ensures comfortable stability.



PRODUCT INFORMATION

INDICATIONS

KONSTANZ is particularly suited for patients with wounds (e.g. traumatic injuries) or swellings (caused by other indications) who need an open shoe with flexible sole and adaptable inner volume.

SPECIAL FEATURES

KONSTANZ is a light shoe which provides much air and space. Due to the open design and by using the comfortable hook and loop fasteners, this shoe can easily be put on and taken off.

ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TH4562-35	TH4562-35/L	TH4562-35/R	43	TH4562-43	TH4562-43/L	TH4562-43/R
36	TH4562-36	TH4562-36/L	TH4562-36/R	44	TH4562-44	TH4562-44/L	TH4562-44/R
37	TH4562-37	TH4562-37/L	TH4562-37/R	45	TH4562-45	TH4562-45/L	TH4562-45/R
38	TH4562-38	TH4562-38/L	TH4562-38/R	46	TH4562-46	TH4562-46/L	TH4562-46/R
39	TH4562-39	TH4562-39/L	TH4562-39/R	47	TH4562-47	TH4562-47/L	TH4562-47/R
40	TH4562-40	TH4562-40/L	TH4562-40/R	48	TH4562-48	TH4562-48/L	TH4562-48/R
41	TH4562-41	TH4562-41/L	TH4562-41/R	49	TH4562-49	TH4562-49/L	TH4562-49/R
42	TH4562-42	TH4562-42/L	TH4562-42/R	50	TH4562-50	TH4562-50/L	TH4562-50/R



PRODUCT INFORMATION

INDICATIONS

WEIMAR is particularly suited for patients with wounds (e.g. traumatic injuries) or swellings (caused by other indications) who need a half-open shoe with flexible sole and adaptable inner volume.

SPECIAL FEATURES

WEIMAR is a light shoe which - despite its sandal shape - has a toe cap to protect the toe area (e.g. against bumps).

ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS4502-35	TS4502-35/L	TS4502-35/R	43	TS4502-43	TS4502-43/L	TS4502-43/R
36	TS4502-36	TS4502-36/L	TS4502-36/R	44	TS4502-44	TS4502-44/L	TS4502-44/R
37	TS4502-37	TS4502-37/L	TS4502-37/R	45	TS4502-45	TS4502-45/L	TS4502-45/R
38	TS4502-38	TS4502-38/L	TS4502-38/R	46	TS4502-46	TS4502-46/L	TS4502-46/R
39	TS4502-39	TS4502-39/L	TS4502-39/R	47	TS4502-47	TS4502-47/L	TS4502-47/R
40	TS4502-40	TS4502-40/L	TS4502-40/R	48	TS4502-48	TS4502-48/L	TS4502-48/R
41	TS4502-41	TS4502-41/L	TS4502-41/R	49	TS4502-49	TS4502-49/L	TS4502-49/R
42	TS4502-42	TS4502-42/L	TS4502-42/R	50	TS4502-50	TS4502-50/L	TS4502-50/R





PRODUCT INFORMATION

INDICATIONS

HEIDELBERG is particularly suited for patients with wounds (e.g. traumatic injuries) or swellings (caused by other indications) who need a closed shoe with flexible sole and adaptable inner volume.

SPECIAL FEATURES

HEIDELBERG provides a variable volume. The partially flexible upper fabric of this shoe is elastic and adapts to different demands on space (e.g. due to varying or thick bandages).

ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS4552-35	TS4552-35/L	TS4552-35/R	43	TS4552-43	TS4552-43/L	TS4552-43/R
36	TS4552-36	TS4552-36/L	TS4552-36/R	44	TS4552-44	TS4552-44/L	TS4552-44/R
37	TS4552-37	TS4552-37/L	TS4552-37/R	45	TS4552-45	TS4552-45/L	TS4552-45/R
38	TS4552-38	TS4552-38/L	TS4552-38/R	46	TS4552-46	TS4552-46/L	TS4552-46/R
39	TS4552-39	TS4552-39/L	TS4552-39/R	47	TS4552-47	TS4552-47/L	TS4552-47/R
40	TS4552-40	TS4552-40/L	TS4552-40/R	48	TS4552-48	TS4552-48/L	TS4552-48/R
41	TS4552-41	TS4552-41/L	TS4552-41/R	49	TS4552-49	TS4552-49/L	TS4552-49/R
42	TS4552-42	TS4552-42/L	TS4552-42/R	50	TS4552-50	TS4552-50/L	TS4552-50/R



PRODUCT INFORMATION

INDICATIONS

LEIPZIG is particularly suited for patients with wounds (e.g. traumatic injuries) or swellings (caused by other indications) who need a closed shoe with flexible sole and a huge, adaptable inner volume.

SPECIAL FEATURES

LEIPZIG provides much space for specific demands. Due to the generous hook and loop fasteners and the dorsal zip, this shoe can be opened very wide and thus easily put on and taken off. In addition, this shoe has an integrated toe cap for more stability and protection in the toe area.

ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5512-35	TS5512-35/L	TS5512-35/R	43	TS5512-43	TS5512-43/L	TS5512-43/R
36	TS5512-36	TS5512-36/L	TS5512-36/R	44	TS5512-44	TS5512-44/L	TS5512-44/R
37	TS5512-37	TS5512-37/L	TS5512-37/R	45	TS5512-45	TS5512-45/L	TS5512-45/R
38	TS5512-38	TS5512-38/L	TS5512-38/R	46	TS5512-46	TS5512-46/L	TS5512-46/R
39	TS5512-39	TS5512-39/L	TS5512-39/R	47	TS5512-47	TS5512-47/L	TS5512-47/R
40	TS5512-40	TS5512-40/L	TS5512-40/R	48	TS5512-48	TS5512-48/L	TS5512-48/R
41	TS5512-41	TS5512-41/L	TS5512-41/R	49	TS5512-49	TS5512-49/L	TS5512-49/R
42	TS5512-42	TS5512-42/L	TS5512-42/R	50	TS5512-50	TS5512-50/L	TS5512-50/R



INNER SHOE DIMENSIONS

Dear Customer,

We measured our therapeutic shoes according to the Mondopoint system. Afterwards, the shoe sizes were recalculated by means of the Paris Point.

Mondopoint System

The Mondopoint system is an international standardised system to indicate shoe sizes. In contrast to other systems, not only the inner shoe length but also the inner shoe width is measured.


Paris Point

To determine the shoe sizes we used the Paris Point which is calculated as follows:

shoe size = inner shoe length in cm x 1.5
Length and width are indicated in mm.

ATTENTION!

In this catalogue, we state outdoor shoe sizes. These shoe sizes do not cover additional space for bandages or the like.




ASYMMETRIC SHOES

Inner shoe length

The inner shoe length (c) was measured from the middle of the heel area to the middle of the great toe with a standard measurement device.

Inner shoe widths (ball and heel area)

The inner shoe widths were measured at the insole. We measured from the toe-cap and drew a crossline at 1/3 of the insole for the inner shoe width in ball area. At this right-angled straight line a 15° angle was measured. The resulting line was where we measured the inner shoe width in ball area (a). The inner shoe width in heel area (b) was measured at the level of the middle of the heel.



INNER SHOE DIMENSIONS STANDARD (F00053)

Shoe size	Inner shoe dimensions [mm]			Shoe size	Inner shoe dimensions [mm]		
	a	b	c		a	b	c
35	93	66	231	43	106	75	288
36	94	67	237	44	109	76	295
37	95	69	245	45	109	79	299
38	98	69	250	46	110	79	305
39	99	71	258	47	111	81	313
40	101	73	266	48	113	81	318
41	104	74	271	49	115	81	325
42	105	75	277	50	115	81	332

INSOLE SET STANDARD

Shoe size	Art. no. set	Shoe size	Art. no. set
35	F00053-35	43	F00053-43
36	F00053-36	44	F00053-44
37	F00053-37	45	F00053-45
38	F00053-38	46	F00053-46
39	F00053-39	47	F00053-47
40	F00053-40	48	F00053-48
41	F00053-41	49	F00053-49
42	F00053-42	50	F00053-50



The top insole is antimicrobial. The insoles are made of EVA material and provide a different shore hardness: grey = 46°, white = 23°, antimicrobial = 28°

THERAPEUTIC SHOES POST-OP



You will find our therapeutic shoe catalogue
POST-OP on our website in the download area:
www.fior-gentz.com/download

THERAPEUTIC SHOES DIABETES



You will find our therapeutic shoe catalogue
DIABETES on our website in the download area:
www.fior-gentz.com/download

The information in this catalogue is valid at the date of printing. The contained product information serves as a guideline. Subject to technical modifications. All rights, particularly the distribution, copy and translation of this manual or any part of it, in paper or as electronic document, must be authorised in writing by FIOR & GENTZ Gesellschaft für Entwicklung und Vertrieb von orthopädiotechnischen Systemen mbH. Reprints, copies and any other electronic reproduction, even partial, must be authorised in writing by FIOR & GENTZ Gesellschaft für Entwicklung und Vertrieb von orthopädiotechnischen Systemen mbH.



FIOR & GENTZ

Gesellschaft für Entwicklung und Vertrieb
von orthopädietechnischen Systemen mbH

Dorette-von-Stern-Straße 5
21337 Lüneburg (Germany)

☎ +49 4131 24445-0
☎ +49 4131 24445-57

✉ info@fior-gentz.de
🏠 www.fior-gentz.com

FIOR & GENTZ
ORTHOPÄDIETECHNIK MIT SYSTEM