

# S400

Polycentric locking knee, which folds away for a low profile cosmetic line when sitting. The knee can be set to operate as either a SAKL for maximum security or an MKL where the user might be able to utilise the swing function in appropriate settings.



# S200

Secure knee joint for use where locked knee extension is preferred. The lock is released for sitting and then locks automatically when the knee is fully extended.



# S500

Versatile knee with pneumatic swing phase control in a polycentric design which offers good sitting cosmesis. The knee has a training lock feature to assist users during rehabilitation in the early stages.



# endolite S200



## FEATURES

- Rotatable pyramid for ease of alignment
- Low-profile semi-automatic lock
- Will unlock without toe load
- Audible lock engagement for safety

## SPECIFICATIONS

|                   |       |
|-------------------|-------|
| Max. user weight: | 100kg |
| Activity level:   | 1-2   |
| Component weight: | 358g  |
| Build height:     | 85mm  |

## PART NUMBER

|                           |         |
|---------------------------|---------|
| Knee                      | ENS200  |
| Thigh release upgrade kit | P239443 |

# endolite S400



## FEATURES

- Good sitting cosmesis
- Can be set as a SAKL, MKL or free knee
- Socket adapters provide range of set up options
- Geometric lock at heel strike (when set in free knee)
- Adjustable extension assist spring
- Adjustable geometry to reposition ICR

## SPECIFICATIONS

|                   |                    |
|-------------------|--------------------|
| Max. user weight: | 125kg/2<br>100kg/3 |
| Activity level:   | 2-3                |
| Component weight: | 678g               |
| Build height:     | 130mm              |

## PART NUMBER

|                        |         |
|------------------------|---------|
| Knee                   | ENS400  |
| SAKL Thigh Release Kit | P239443 |
| MKL Thigh Release Kit  | 239640  |

## SOCKET ADAPTERS

|                         |                  |
|-------------------------|------------------|
| Lotus Adapter (KD)      | 2-01-S400L=LA-01 |
| Socket Attachment Plate | 2-01-S400=SP-01  |



# endolite S500



## FEATURES

- Slim line polycentric with pneumatic swing control
- Geometric lock at heel strike and stance
- Independent flexion extension control adjustable
- Socket adapters provide range of set up options
- Geometric lock at heel strike (when set in free knee)
- Adjustable spring extension assist

## SPECIFICATIONS

|                   |       |
|-------------------|-------|
| Max. user weight: | 125kg |
| Activity level:   | 2-3   |
| Component weight: | 882g  |
| Build height:     | 220mm |

## PART NUMBER

|      |        |
|------|--------|
| Knee | ENS500 |
|------|--------|

## SOCKET ADAPTERS

|                         |                 |
|-------------------------|-----------------|
| Lotus Adapter (KD)      | 2-01-S500=LA-01 |
| Socket Attachment Plate | 2-01-S500=SP-01 |



Follow Blatchford: | blatchford.co.uk  
 @BlatchfordGrp | 0114 263 7900  
 /BlatchfordGrp | sales@blatchford.co.uk



Endolite is part of the Blatchford Group. Registered address, Chas A Blatchford & Sons Ltd, Unit D, Antura, Kingsland Business Park, Basingstoke, RG24 8PZ, United Kingdom. t: +44 (0) 1256 316600 f: +44 (0) 1256 316710



An annual visual inspection is recommended. Check for visual defects that may affect proper function. Maintenance must be carried out by competent personnel.