



endolite

get busy living


mercury

Hi-ACTIVITY

**Blatchford
Products Ltd**


Lister Road, Basingstoke, Hampshire,
RG22 4AH, United Kingdom
Telephone: +44 (0)1256 316600
Facsimile: +44 (0)1256 479705
Email: sales@blatchford.co.uk
www.endolite.co.uk

**MHPYR (239352)
MHDPYRS
MHDPYRM
239040, 239041**



endolite



Load Line → 0-5mm



endolite




To maximize the function of the knee, distal pyramid based components should be fixed in neutral alignment with the shin components vertical.



Note: Hi-Activity Shin and Kneeball not shown




**MH30 (239252)
MHD30S
MHD30M
239037, 239038**



endolite

Load Line → 0-5mm



endolite