

KX06



DESCRIPTION

The KX06 is a unique polycentric knee with hydraulic swing and stance control. It is designed for the most demanding physical activity of level 3-4 users. The geometric and hydraulic combination draws on the best functionality of each to enhance the overall gait smoothness on all terrain.

FEATURES

- Polycentric knee with hydraulic swing and stance or swing only
- 160° knee flexion
- Water resistant design
- Variable cadence
- Cylinder flexion lock
- Mode switch for cycling
- Multiple proximal adapter choices
- Distal pyramid

SPECIFICATIONS

Max. user weight/activity level:	150kg/3 145kg/4
Component weight:	1.23kg
Build height:	see diagram
Control unit:	endolite hydraulic

PART NUMBERS

KX06
KX06 - SWING


ACCESSORIES

Mercury discontinuous cosmesis kit see page 77



Cylinder adjuster tool 940091

SOCKET ADAPTERS



BUILD HEIGHT

