

DESCRIPTION

A robust, single axis knee for activity level 3-4 users. The Mercury's tough frame and hydraulic control meet the demands of the most active lifestyles. Available with swing or swing and stance control.

FEATURES

- Tough carbon fibre frame
- 125° knee flexion
- Long life needle roller bearings
- Reinforced hydraulic unit
- Variable cadence
- Cylinder flexion lock
- Mode switch for cycling
- Durable urethane kneeler pad
- Compact design

SPECIFICATIONS

Max. user weight/activity level:	150kg/3 145kg/4
Component weight:	1.13kg
Build height:	see diagram
Control unit:	endolite hydraulic



PART NUMBERS

	MH30 MH30 SWING MHPYR MHPYR SWING
Discontinuous:	MHDPYRS - small kneeball 100mm MHDPYRM - medium kneeball 110mm MHD30S - small kneeball 100mm MHD30M - medium kneeball 110mm

ACCESSORIES

	Mercury discontinuous cosmesis kit see page 77
Cylinder adjuster tool	940091

SOCKET ADAPTERS



BUILD HEIGHT

