

# endolite ESK+



Pyramid option shown

## DESCRIPTION

The ESK+ with weight activated stance control and Stanceflex offers a safe and smooth gait for moderately active Level 2-3 users. The pneumatic swing control may be set from regular to dynamic walking pace.

## FEATURES

- Weight activated adjustable stance control
- Pneumatic swing phase control
- Stanceflex option improves comfort
- 140° of knee flexion
- 4-Bolt option reduces build height

## SPECIFICATIONS

Max. user weight:	125kg
Activity level:	2-3
Component weight:	1.13kg Stanceflex 1.15kg 4-Bolt
Build height:	see diagram
Control unit:	Pneumatic

## PART NUMBER

Stanceflex options	ESKSFP30 ESKSFP35 ESKSFPYR
4-Bolt options	ESK4BP30 ESK4BP35 ESK4BPPYR

## STANCEFLEX RUBBER HARDNESS

25°	35°	45°	70°
231608	231606	231607	231609
Soft		Fitted as standard	Hard

Careful selection of Stanceflex rubber is essential as it influences the action of the unit (see fitting information)

## SOCKET ADAPTERS

Stanceflex options

409041	409042	189806	189123 (rotation)	189127

4-Bolt options

EP-MSR	EP-FSR	P189310	239017	239032

189128

## BUILD HEIGHT

