

endolite S200

DESCRIPTION

Secure knee joint for use where locked knee extension is preferred. The lock is released for sitting and then locks automatically when the knee is fully extended.

FEATURES

- Rotatable pyramid for ease of alignment
- Low-profile semi-automatic lock
- Will unlock without toe load
- Audible lock engagement for safety

SPECIFICATIONS

Max. user weight:	100kg
Activity level:	1-2
Component weight:	358g
Build height:	85mm



PART NUMBER

Knee	ENS200
Thigh release upgrade kit	P239443
Service/Repair Kit	2-01-4S40J-SPK

