

endolite S400



DESCRIPTION

Polycentric locking knee, which folds away for a low profile cosmetic line when sitting. The knee can be set to operate as either a SAKL for maximum security or an MKL where the user might be able to utilise the swing function in appropriate settings.

FEATURES

- Compact low-profile polycentric locking knee
- Good sitting cosmesis
- Can be set as a SAKL, MKL or free knee
- Socket Adapters provide range of set up options
- Geometric lock at heel strike (when set in free knee)
- Adjustable extension assist spring
- Adjustable geometry to reposition ICR

SPECIFICATIONS

Max. user weight:	125kg/2 100kg/3
Activity level:	2-3
Component weight:	678g
Build height:	130mm

PART NUMBER

Knee **ENS400**



SAKL Thigh Release Kit	P239443
MKL Thigh Release Kit	239640
Service/Repair Kit	2-01-S400L-SPK

SOCKET ADAPTERS

Lotus Adapter (KD) **2-01-S400L=LA-01**



Socket Attachment Plate **2-01-S400=SP-01**

