


navigator



BUILD HEIGHT



* sizes 22-24 75mm
 25/26 80mm
 27/28 85mm
 29/30 90mm

DESCRIPTION

The stability and comfort provided by this multi-axial foot with flexible keel ensures that activity level 1-3 users are safer on uneven ground. The ankle center is positioned anatomically to promote a natural gait from heel strike to toe off.

FEATURES

- Multi-axial flexible keel
- Integrated ankle joint
- Choice of ankle stiffness
- Sandal toe

Footshell and glide sock included

SPECIFICATIONS

Max. user weight:	125kg
Activity level:	1-3
Size range:	22cm-30cm
Component weight:	565g†
Build height:	95mm
Heel height:	10mm

ORDER EXAMPLE

NAV **25L** **W** **O**
 Size Ball & snubber
 Side

For dark tone add suffix D.
 Foot example: navigator, size 25 left, ball rating white, snubber rating orange.

Ball and Snubber
 W=white, soft
 O=orange, medium
 B=blue, hard

Activity	User Weight							kg lbs
	44-52 100-115	53-59 116-130	60-68 131-150	69-77 151-170	78-88 171-195	89-100 196-220	101-125 221-275	
1	WW	WW	WW	WO	OO	OB	BB	Ball and snubber
2	WW	WW	WO	OO	OB	BB	BB	
3	WW	WO	OO	OB	BB	BB	BB	

†Component weight shown is for a size 26cm