




## BUILD HEIGHT



### Multiflex Ankle Standard 30mm and 35mm

	30mm	35mm
Ankle Kits:	379342 (W)	409125 (W)
	379343 (O)	409126 (O)
	379344 (B)	409127 (B)
Snubber Kits:	409007 (W)	
	409008 (O)	
	409009 (B)	



## DESCRIPTION

All the multi-axial benefits of the Multiflex system within a standard footshell profile.

## FEATURES

- Multi-axial function
- Ankle level footshell
- Sandal toe

## SPECIFICATIONS

Max. user weight:	125kg
Activity level:	1-2
Size range:	25cm-28cm
Component weight:	425g <sup>†</sup>
Build height:	40mm
Heel height:	10 ± 15mm

## PART NUMBERS

Size	Left	Right
25cm	519525	519526
26cm	519527	519528
27cm	519529	519530
28cm	519531	519532

### Multiflex Ankle Pyramid

Ankle Kits:	409353 (W)
	409354 (O)
	409355 (B)
Snubber Kits:	409007 (W)
	409008 (O)
	409009 (B)



Multiflex ankle and snubber selection.

Activity	User Weight								kg
	44-52	53-59	60-68	69-77	78-88	89-100	101-125	126-150	
1-2	WW	WW	WW	WO	OO	OB	BB		Ball and snubber

<sup>†</sup>Component weight shown is for a size 26cm