

DESCRIPTION

A cushioned foot with flexible keel for Level 1-2 activity. An adjustable heel screw allows heel cushion adjustment for an individual fit.

FEATURES

- Flexible keel
- Adjustable heel cushion
- Sandal toe
- Integrated pyramid

SPECIFICATIONS

Max. user weight:	125kg
Activity level:	1-2
Size range:	22cm-28cm
Component weight:	449g†
Build height:	70mm size 22-24 75mm size 25-26 80mm size 27-28
Heel height:	19mm



PART NUMBERS

Size	Left	Right
22cm	P529101	P529102
23cm	P529103	P529104
24cm	P529105	P529106
25cm	P529107	P529108
26cm	P529109	P529110
27cm	P529111	P529112
28cm	P529113	P529114

BUILD HEIGHT

Size	A
22cm	70mm
23cm	
24cm	
25cm	75mm
26cm	
27cm	80mm
28cm	

†Component weight shown is for a size 26cm