

senior



BUILD HEIGHT

Size	A
22cm	70mm
23cm	
24cm	
25cm	75mm
26cm	
27cm	80mm
28cm	

DESCRIPTION

An ultra lightweight foot for level 1-2 activity with excellent heel cushion for walking comfort.

FEATURES

- Replaceable foot shell
- Sandal toe
- Flexible keel
- Integrated pyramid

SPECIFICATIONS

Max. user weight:	125kg
Activity level:	1-2
Size range:	22cm-28cm
Component weight:	395g†
Build height:	70mm size 22-24 75mm size 25-26 80mm size 27-28
Heel height:	19mm

ORDER EXAMPLE

SNR

25
Size

L
Side

For dark tone add suffix D.
Foot example: senior, size
25 left.