

## DESCRIPTION

A waterproof prosthesis with integral shin designed for barefoot use in the shower and wet environments. The moulded components allow fast assembly.

## FEATURES

- Integral shin and cosmesis
- Rapid assembly
- Full alignment device included
- Anti-slip tread pattern on sole

Please call Customer Services for a lead time on 80kg Aqualimb

## SPECIFICATIONS

Max. user weight:	80kg system / 100kg system
Activity level:	1-2
Size range:	22cm-27cm
Component weight:	1.39kg†100kg
Build height:	280 to 415mm 100kg 100 to 415mm 80kg
Heel height:	5mm



## PART NUMBERS

### 100kg:


Size	Left	Right
22cm	629705L	629706L
23cm	629707L	629708L
24cm	629709L	629710L
25cm	629711L	629712L
26cm	629713L	629714L
27cm	629715L	629716L

### 80kg:


Size	Left	Right
22cm	629405LF	629406LF
23cm	629407LF	629408LF
24cm	629409LF	629410LF
25cm	629411LF	629412LF
26cm	629413LF	629414LF
27cm	629415LF	629416LF

## BUILD HEIGHT

80kg:



100kg:



†Component weight shown is for a size 26cm