

DESCRIPTION

The SmartIP employs intelligent programming technology that dynamically programs the knee to optimum swing settings for the individual Level 2-3 activity user, reducing the amount of energy expended by trans-femoral users.

FEATURES

- Comfortable gait from very slow to very fast walking speeds
- Smoother speed progression because of enhanced program
- Weight activated stance control
- Stanceflex option improves comfort
- Battery life expectancy 9-12 months under normal usage
- 140° knee flexion
- 4-Bolt option reduces build height

SPECIFICATIONS

Max. user weight:	125kg
Activity level:	2-3
Component weight:	1.27kg Stanceflex 1.15kg 4-Bolt
Build height:	see diagram
Control unit:	Swing: MPK/Pneumatic Stance: mechanical



Pyramid option shown

PART NUMBERS

Stanceflex options	SMARTIPSF30 SMARTIPSF30PYR
4-Bolt options	SMARTIP4B30 SMARTIP4BPYR

STANCEFLEX RUBBER HARDNESS

25°	35°	45°	70°
231608	231606	231607	231609
Soft		Fitted as standard	Hard

Careful selection of Stanceflex rubber is essential as it influences the action of the unit (see fitting information)

SOCKET ADAPTERS

Stanceflex options



4-Bolt options



BUILD HEIGHT

