

endolite S500

DESCRIPTION

Versatile knee with pneumatic swing phase control in a polycentric design which offers good sitting cosmesis. The knee has a training lock feature to assist users during rehabilitation in the early stages.

FEATURES

- Slim line polycentric with pneumatic swing control
- Geometric lock at heel strike and stance
- Independent flexion/extension control adjustable
- Socket adapters provide range of set up options
- Geometric lock at heel strike (when set in free knee)
- Adjustable spring extension assist

SPECIFICATIONS

Max. user weight:	125kg
Activity level:	2-3
Component weight:	882g
Build height:	220mm

PART NUMBER

Knee	ENS500
Service/Repair Kit	2-01-S500-SPK

SOCKET ADAPTERS

Lotus Adapter (KD)	2-01-S500=LA-01
Socket Attachment Plate	2-01-S500=SP-01

