

STEALTH360™, PYR-SL-TI™

Pyramid: Sliding or Rotating Product Information

The **STEALTH360™** is a 360° rotating pyramid made of titanium with a base of hard-anodized aluminum. It is compatible with the European 4-hole pattern and the OB harmony system. The rotation is locked into place by tightening the 4 adjustment set screws on the interfacing connecting component, such as PDI's PR-SL. If connecting to a PDI lock, torque the M6 bolts to 12 Nm (9 ft lb). **US Pat No. 6,033,440 and 6,231,618B.**

Suggested L-Code:

L5910 – BK
L5920 – AK



STEALTH360

SPECIFICATIONS

Added Height:	7.53 mm (0.3 in)
Weight:	59 g (1.9 oz)
Load Rating:	
K1-K2 Activity Level	265 lbs (120 Kg)
K3-K4 Activity Level	220 lbs (100 Kg)

The **PYR-SL-TI™** features a 12.7 mm (1/2") slide off center, with closed ends on the base and can be oriented for either AP or ML slide. This component is manufactured in titanium and is compatible with the European 4-hole pattern. The slide is locked into place by tightening the 4 adjustment set screws on the interfacing connection component. If attach to a PDI lock, torque the M6 bolts to 12Nm (9ft/lb). **US Pat No. 6,033,440 and 6,231,618.**

Suggested L-Code:

L5910 – BK
L5920 – AK



PYR-SL-TI

SPECIFICATIONS

Added Height:	13.48 mm (0.53 in)
Weight:	123 g (4.34 oz)
Load Rating:	
K1-K2 Activity Level	265 lbs (120 Kg)
K3-K4 Activity Level	220 lbs (100 Kg)

Warranty Information: PDI warrants its components to purchasers of the prosthetic components only. No warranties are made to consumers. All components are warranted for 30 days from original purchase date. This warranty applies only to manufacturing defects and does not cover: (1) any damage due to misuse, altering of the component, improper service or accidents; (2) any secondary costs that might be incurred in the exchange of components; (3) wear; normal, or due to incorrect selection of components. During the warranty period, defective parts will be replaced or repaired at no charge. Missing bolts will be replaced at retail prices. This warranty is in lieu of all other warranties, expressed or implied, including any implied warranty of merchantability or fitness for particular purpose. The remedies provided above are exclusive; in no event shall Prosthetic Design, Inc. be liable for any special, indirect, incidental, or consequential damages.

Returns: A restocking fee totaling 2% of retail price will be charged to customer should the customer wish to return any PDI product for credit. The part must be unused, and in its original, unopened packaging. In addition, the return must be initiated within 30 days of its purchase in order to qualify for the credit. If returned after 30 days of purchase date a 15% restocking fee will be charged.



Prosthetic Design, Inc., 700 Harco Drive, Clayton, OH 45315 USA
www.prostheticdesign.com Telephone: 800-459-0177
International Telephone: 937-836-1464 Fax: 937-832-5361

