FC15-TD8CFI, FC15-TD10CFI, FC15-TD12CFI Flexcon™ Product Information

IMPROVING LIVES THROUGH INNOVATION

Introduction: PDI's lightweight *carbon fiber* FLEXCONTM is a modular alternative to labor–intensive laminations for patients with flexion contractures. The FLEXCONTM flexes slightly under weight-bearing forces, enhancing gait during toe-off. Simply attach FLEXCONTM between PDI's pyramid and PDI's Modular Lock. The FLEXCONTM compensates for the flexion contracture angle and creates an offset based on the length of the residual limb and the angle of the flexion contracture. **US Patent** # 6,991,658

Choose the correct FLEXCONTM for your application:



Carbon Fiber
Limb Length from Trochanter to Distal

15° Flexion Contracture Angle

Suggested L-Code: L5845

Offset	PDI Part Number
2.07"	FC15-TD8CFI
2.59"	FC15-TD10CFI
3.11"	FC15-TD12CFI



SPECIFICATIONS For FC15-TD8CF

Thickness: 9.5 mm (3/8")
Weight: 98 g (3.5 oz)

Load Rating:

High Activity 220 lbs (100 Kg)
Moderate Activity 265 lbs (120 Kg)

SPECIFICATIONS For FC15-TD10CF

Thickness: 9.5 mm (3/8")
Weight: 98 q (3.5 oz)

Load Rating:

High Activity 220 lbs (100 Kg)
Moderate Activity 265 lbs (120 Kg)

SPECIFICATIONS For FC15-TD12CF

Thickness: 9.5 mm (3/8")
Weight: 98 g (3.5 oz)

Load Rating:

High Activity 220 lbs (100 Kg)
Moderate Activity 265 lbs (120 Kg)

Over for more information

