

FIOR & GENTZ
ORTHOPÄDIETECHNIK MIT SYSTEM



THERAPEUTIC SHOES
DIABETES
MAY 2016



Dear Customers,

We are happy to present you our therapeutic shoes in a completely new light.

The catalogue at hand contains all FIOR & GENTZ therapeutic shoes especially designed for diabetic feet.

To ensure a consistent communication, we summarised a diabetic foot syndrome classification on pages 8-9. You will find references to this classification on each product page.

The table of contents offers an overview of the varied shaft and sole designs. However, all the shoes' details can be found on each shoe's particular product page. Matching insoles are also described on these pages.

You can find information on shoe sizes and measurement methods from page 34 and matching accessory parts, such as protective caps, from page 36 on.

Our goal is to provide you with products enabling you to treat your patients the best possible way. You will find further therapeutic shoes in our **product catalogues POST-OP and HOMECARE.**

With best regards from Lüneburg
Jörg Fior and Ralf Gentz

All our products meet our high standard of quality. This quality standard has been controlled and certified by an independent certification organisation according to the international standards ISO 9001 and ISO 13485.



SALES TEAM IN-HOUSE

In order to provide you with up to date contact details of our distributors in foreign countries, we have listed all our current distributors online at: www.fior-gentz.com/contact/export

Ana María von Corvin
English/Spanish



Michael Krüger
English



Anna-Lena Kühn
English/French



SALES TEAM BENELUX

Luc Schols
Dutch/English/French



IN-HOUSE SALES TEAM

For any questions, orders or customer data changes just call +49 4131 24445-0. We are available on workdays from 8.00am to 5.00pm (Central European Time).

If you place your order by 4.15pm, we will dispatch the goods on that same day. Express deliveries on Saturday are also possible.

INTRODUCTION

WELCOME.....	2-3
SERVICE AND SUPPORT.....	4-5
DIABETIC FOOT SYNDROME.....	8-9
SHOE CHARACTERISTICS.....	10-11

SHOES FOR BANDAGES

ASYMMETRIC SHOES WITH TALL SHAFT

Flexible sole with pitch, 30mm thick

LEIPZIG.....	12-13
--------------	-------

Rigid, extra strong rocker bottom sole with pitch, 30mm thick

FLENSBURG.....	14-15
AUGSBURG.....	16-17
DRESDEN.....	18-19
KASSEL.....	20-21

ASYMMETRIC SHOES WITH LOW SHAFT

Flexible sole with pitch, 30mm thick

KONSTANZ.....	22-23
WEIMAR.....	24-25
HEIDELBERG.....	26-27

Rigid, extra strong rocker bottom sole with pitch, 30mm thick

NÜRNBERG.....	28-29
POTSDAM.....	30-31

HEEL RELIEF SHOES

ASYMMETRIC SHOES WITH TALL SHAFT

Rigid, extra strong rocker bottom sole with pitch, 30mm thick

WETZLAR.....	32-33
--------------	-------

ADDITIONAL INFORMATION

INNER SHOE DIMENSIONS.....	34-35
SPARE AND ACCESSORY PARTS.....	36-38
MORE SHOE CATALOGUES.....	39



CLASSIFICATION OF THE DIABETIC FOOT SYNDROME

The Diabetic Foot Syndrome (DFS) results from diabetes mellitus and involves a whole series of various symptoms. Either a neuropathic or ischaemic foot or a combination of both is causing a diabetic foot syndrome.

If it is the former, the nerves are damaged because of a permanently heightened blood glucose level which is weakening the muscles and causing the ligamentous apparatus to collapse.

Then, diabetic patients experience symptoms such as, for example, numbness, pressure points on the foot's sole and/or the toes as well as dry and vulnerable skin. Manifesting ulcers can become chronic wounds which only heal poorly, if at all.

If it is an ischaemic foot, the arterial blood flow is interrupted. This may lead to the death of entire tissue areas.

Generally, wound treatment is an essential part of DFS therapy as the wounds' ability to heal is often impaired. And when, in addition, weight had been put on the foot the wrong way, this may result in amputations of the extremities.

The cause of DFS must definitely be considered when treating diabetic feet as each therapy differs from the other in many ways or is even contraindicated.

Moreover, a correct assessment of the existent risk group is decisive for avoiding (further) injuries and to ensure the proper treatment of each individual patient.

This catalogue provides diabetic shoes for risk group III and VII. You can see which risk groups and wound stages each diabetic shoe is especially recommended for in the paragraph Product Information at Indications on each shoe page. Our specifications are based on the classifications on the next page.

Risk Groups of the Diabetic Foot Syndrome*

Risk group 0	Diabetes mellitus without peripheral neuropathy (PNP)/peripheral artery disease (PAD)
Risk group I	Same as risk group 0, with foot deformity
Risk group II	Diabetes mellitus with sensory neuropathy caused by PNP/relevant PAD
Risk group III	Condition after plantar ulcer
Risk group IV	Same as risk group II with deformities or disproportions
Risk group V	Diabetic neuro-osteoarthropathy (DNOAP) (LEVIN III)
Risk group VI	Same as risk group II with partial foot amputation
Risk group VII	Severe lesion

* according to "AG Fuß: Arbeitsgemeinschaft Diabetischer Fuß der Deutschen Diabetes Gesellschaft" (Working group on the diabetic foot in the German Diabetic Society) (as of 25th February 2006)

Wound Stages of Diabetic Feet**

	Grade 0	Grade 1	Grade 2	Grade 3	Grade 4	Grade 5
A	Pre- or post-ulcerative lesion	Superficial wound	Wound up to level of tendons or capsule	Wound up to level of bones and joints	Necrosis of parts of foot	Necrosis of the entire foot
B	With infection	With infection	With infection	With infection	With infection	With infection
C	With ischaemia	With ischaemia	With ischaemia	With ischaemia	With ischaemia	With ischaemia
D	With infection and ischaemia	With infection and ischaemia	With infection and ischaemia	With infection and ischaemia	With infection and ischaemia	With infection and ischaemia

** according to Wagner (grade 0-5) and Armstrong (stage A-D)

Classification of the Diabetic Foot Syndrome According to Wagner

Grade 0	Pre- and post-ulcerative foot
Grade 1	Superficial wound, but affected skin and subcutaneous tissue
Grade 2	Wound up to the level of tendons or capsule
Grade 3	Wound up to level of bones and joints
Grade 4	Necrosis of parts of foot
Grade 5	Necrosis of the entire foot

Armstrong's Stages of Diabetic Foot Syndrome

Stage A	Neuropathic foot without infection
Stage B	Neuropathic foot with infection
Stage C	With ischaemia
Stage D	With infection and ischaemia

DIABETIC SHOES

Risk group III

Usually, a diabetic shoe with either a diabetic foot orthotic or an orthopaedic shoe modification is used for any condition after a plantar ulcer. That is why we offer an insert blank/an insole set which can be reworked to relieve plantar pressure points and used in all our diabetic shoes. The insert/insole can be worn as an interim solution until the diabetic foot orthotic customised to the patient's needs is available. Furthermore, our diabetic shoes offer sufficient space so that a thick insert/insole/diabetic foot orthotic can easily be inserted without having to squeeze the foot into the shoe. If even more space should be necessary, because of bandages, foot malpositions/deformities or the like, the shoes' inner volume can still be adapted, the (inner) fabric reworked or you could choose not a STANDARD but a WIDE shoe.

Risk group VII

In general, a relief/bandage/interim or orthotic shoe with a diabetic foot orthotic or an orthopaedic shoe modification is used for any lesions. Our shoes provided for such an indication - especially our WIDE shoes - offer sufficient space for individual foot orthotics as well as for additionally worn bandages and/or orthoses. Due to its hook and loop fasteners, WIDE versions and the reworkable (inner) fabric, our shoes can also relieve pressure to the point on non-plantar shoe parts. Additionally, the rigid sole design and the tall shaft ensure stability and security while walking and standing.




WASHABLE

All diabetic shoes can be washed at 40°C.

BREATHABLE

The upper fabric of all shoes is breathable.



OUTSOLES

The soles are 30mm thick and made of solid EVA material. The rigid soles can be adapted to the patient's needs by grinding and/or adhering material.



REINFORCED HEEL CAP

All diabetic shoes have a reinforced heel cap which provides support, cushions the heel and ensures comfortable stability.



FASTENER SYSTEM

Due to hook and loop fasteners, all diabetic shoes are easy to open and close. While wearing the shoe, large loop areas ensure individual comfort. Besides, the hook and loop fasteners have additional seams so that the fabric can be neatly cut, if necessary.



ADAPTABLE INNER VOLUME

Due to hook and loop fasteners - which can be neatly cut - and removable insoles, the inner volume can be adapted to the patient's individual need for space.



PRODUCT INFORMATION

INDICATIONS

LEIPZIG is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

Risk group III: condition after plantar ulcer
Grade 0: condition after lesion, vulnerable foot

SOLE DESIGN

LEIPZIG has a 30mm thick, flexible sole with pitch. This sole design supports a natural heel to toe movement while walking.

INSOLES


Inside the shoe you will find three 5mm thick insoles made of EVA material (F00053-XX). The top insole has an antimicrobial coating. You can use these insoles as an interim solution for pressure relief until replacing them by a diabetic foot orthotic.

SPECIAL FEATURES

LEIPZIG provides much space for specific demands. Due to the hook and loop fasteners and the dorsal zip, this shoe can be opened very wide and thus easily be put on and taken off. In addition, this shoe has an integrated toe cap for more stability and protection in the toe area.

ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5512-35	TS5512-35/L	TS5512-35/R	43	TS5512-43	TS5512-43/L	TS5512-43/R
36	TS5512-36	TS5512-36/L	TS5512-36/R	44	TS5512-44	TS5512-44/L	TS5512-44/R
37	TS5512-37	TS5512-37/L	TS5512-37/R	45	TS5512-45	TS5512-45/L	TS5512-45/R
38	TS5512-38	TS5512-38/L	TS5512-38/R	46	TS5512-46	TS5512-46/L	TS5512-46/R
39	TS5512-39	TS5512-39/L	TS5512-39/R	47	TS5512-47	TS5512-47/L	TS5512-47/R
40	TS5512-40	TS5512-40/L	TS5512-40/R	48	TS5512-48	TS5512-48/L	TS5512-48/R
41	TS5512-41	TS5512-41/L	TS5512-41/R	49	TS5512-49	TS5512-49/L	TS5512-49/R
42	TS5512-42	TS5512-42/L	TS5512-42/R	50	TS5512-50	TS5512-50/L	TS5512-50/R





PRODUCT INFORMATION

INDICATIONS

FLENSBURG is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

- Risk group III:** condition after plantar ulcer
- Grade 0:** condition after lesion
- Risk group VII:** severe lesion
- Grade 1:** superficial wound
- Grade 2:** wounds up to tendon/capsule, plantar ulcer at the forefoot area
- Grade 4:** endemic necrosis at the forefoot area, condition after toe amputation
- Grade 4-5:** partial forefoot amputation, e.g. condition after single toe removals but preservation of the metatarsophalangeal joint of the big toe



SHOE TYPES

The asymmetric diabetic shoe FLENSBURG is available in a STANDARD and WIDE version.

SOLE DESIGN

FLENSBURG has a 30mm thick, rigid and extra strong rocker bottom sole with pitch. The rigid sole prevents flexion of the metatarsophalangeal joints and ensures the patient's foot being stabilised and immobilised. That is how an existent wound is relieved and the wound's healing supported in an optimal way.



INSOLES

Inside the shoe you will find three 5mm thick insoles made of EVA material (STANDARD F00053-XX and WIDE F00054-XX). The top insole has an antimicrobial coating. You can use these insoles as an interim solution for pressure relief until replacing them by a diabetic foot orthotic.



SPECIAL FEATURES


FLENSBURG, especially the WIDE shoe, provides much space for specific demands. By using the dorsal zip, this shoe can be put on and taken off without even using the hook and loop fasteners. In addition, this shoe has an integrated toe cap for more stability and protection in the toe area.

FLENSBURG STANDARD

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5311-35	TS5311-35/L	TS5311-35/R
36	TS5311-36	TS5311-36/L	TS5311-36/R
37	TS5311-37	TS5311-37/L	TS5311-37/R
38	TS5311-38	TS5311-38/L	TS5311-38/R
39	TS5311-39	TS5311-39/L	TS5311-39/R
40	TS5311-40	TS5311-40/L	TS5311-40/R
41	TS5311-41	TS5311-41/L	TS5311-41/R
42	TS5311-42	TS5311-42/L	TS5311-42/R
43	TS5311-43	TS5311-43/L	TS5311-43/R
44	TS5311-44	TS5311-44/L	TS5311-44/R
45	TS5311-45	TS5311-45/L	TS5311-45/R
46	TS5311-46	TS5311-46/L	TS5311-46/R
47	TS5311-47	TS5311-47/L	TS5311-47/R
48	TS5311-48	TS5311-48/L	TS5311-48/R
49	TS5311-49	TS5311-49/L	TS5311-49/R
50	TS5311-50	TS5311-50/L	TS5311-50/R

FLENSBURG WIDE

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5911-35	TS5911-35/L	TS5911-35/R
36	TS5911-36	TS5911-36/L	TS5911-36/R
37	TS5911-37	TS5911-37/L	TS5911-37/R
38	TS5911-38	TS5911-38/L	TS5911-38/R
39	TS5911-39	TS5911-39/L	TS5911-39/R
40	TS5911-40	TS5911-40/L	TS5911-40/R
41	TS5911-41	TS5911-41/L	TS5911-41/R
42	TS5911-42	TS5911-42/L	TS5911-42/R
43	TS5911-43	TS5911-43/L	TS5911-43/R
44	TS5911-44	TS5911-44/L	TS5911-44/R
45	TS5911-45	TS5911-45/L	TS5911-45/R
46	TS5911-46	TS5911-46/L	TS5911-46/R
47	TS5911-47	TS5911-47/L	TS5911-47/R
48	TS5911-48	TS5911-48/L	TS5911-48/R
49	TS5911-49	TS5911-49/L	TS5911-49/R
50	TS5911-50	TS5911-50/L	TS5911-50/R



PRODUCT INFORMATION

INDICATIONS

AUGSBURG is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

- Risk group III:** condition after plantar ulcer
- Grade 0:** condition after lesion
- Risk group VII:** severe lesion
- Grade 1:** superficial wound
- Grade 2:** wounds up to tendon/capsule, plantar ulcer at the forefoot area
- Grade 4:** endemic necrosis at the forefoot area, condition after toe amputation
- Grade 4-5:** partial forefoot amputation, e.g. condition after single toe removals but preservation of the metatarsophalangeal joint of the big toe



SHOE TYPES

The asymmetric diabetic shoe AUGSBURG is available in a STANDARD and WIDE version.

SOLE DESIGN

AUGSBURG has a 30mm thick, rigid and extra strong rocker bottom sole with pitch. The rigid sole prevents flexion of the metatarsophalangeal joints and ensures the patient's foot being stabilised and immobilised. That is how an existent wound is relieved and the wound's healing supported in an optimal way.



INSOLES

Inside the shoe you will find three 5mm thick insoles made of EVA material (STANDARD F00053-XX and WIDE 00054-XX). The top insole has an antimicrobial coating. You can use these insoles as an interim solution for pressure relief until replacing them by a diabetic foot orthotic.



SPECIAL FEATURES

AUGSBURG, especially the WIDE shoe, provides much space for bandages and other specific demands. Besides, and if applicable, a protective forefoot cap can be attached to the toe area that keeps the foot dry and warm.

AUGSBURG STANDARD

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5381-35	TS5381-35/L	TS5381-35/R
36	TS5381-36	TS5381-36/L	TS5381-36/R
37	TS5381-37	TS5381-37/L	TS5381-37/R
38	TS5381-38	TS5381-38/L	TS5381-38/R
39	TS5381-39	TS5381-39/L	TS5381-39/R
40	TS5381-40	TS5381-40/L	TS5381-40/R
41	TS5381-41	TS5381-41/L	TS5381-41/R
42	TS5381-42	TS5381-42/L	TS5381-42/R
43	TS5381-43	TS5381-43/L	TS5381-43/R
44	TS5381-44	TS5381-44/L	TS5381-44/R
45	TS5381-45	TS5381-45/L	TS5381-45/R
46	TS5381-46	TS5381-46/L	TS5381-46/R
47	TS5381-47	TS5381-47/L	TS5381-47/R
48	TS5381-48	TS5381-48/L	TS5381-48/R
49	TS5381-49	TS5381-49/L	TS5381-49/R
50	TS5381-50	TS5381-50/L	TS5381-50/R

AUGSBURG WIDE

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5481-35	TS5481-35/L	TS5481-35/R
36	TS5481-36	TS5481-36/L	TS5481-36/R
37	TS5481-37	TS5481-37/L	TS5481-37/R
38	TS5481-38	TS5481-38/L	TS5481-38/R
39	TS5481-39	TS5481-39/L	TS5481-39/R
40	TS5481-40	TS5481-40/L	TS5481-40/R
41	TS5481-41	TS5481-41/L	TS5481-41/R
42	TS5481-42	TS5481-42/L	TS5481-42/R
43	TS5481-43	TS5481-43/L	TS5481-43/R
44	TS5481-44	TS5481-44/L	TS5481-44/R
45	TS5481-45	TS5481-45/L	TS5481-45/R
46	TS5481-46	TS5481-46/L	TS5481-46/R
47	TS5481-47	TS5481-47/L	TS5481-47/R
48	TS5481-48	TS5481-48/L	TS5481-48/R
49	TS5481-49	TS5481-49/L	TS5481-49/R
50	TS5481-50	TS5481-50/L	TS5481-50/R



DRESDEN STANDARD

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5811-35	TS5811-35/L	TS5811-35/R
36	TS5811-36	TS5811-36/L	TS5811-36/R
37	TS5811-37	TS5811-37/L	TS5811-37/R
38	TS5811-38	TS5811-38/L	TS5811-38/R
39	TS5811-39	TS5811-39/L	TS5811-39/R
40	TS5811-40	TS5811-40/L	TS5811-40/R
41	TS5811-41	TS5811-41/L	TS5811-41/R
42	TS5811-42	TS5811-42/L	TS5811-42/R
43	TS5811-43	TS5811-43/L	TS5811-43/R
44	TS5811-44	TS5811-44/L	TS5811-44/R
45	TS5811-45	TS5811-45/L	TS5811-45/R
46	TS5811-46	TS5811-46/L	TS5811-46/R
47	TS5811-47	TS5811-47/L	TS5811-47/R
48	TS5811-48	TS5811-48/L	TS5811-48/R

DRESDEN WIDE

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5411-35	TS5411-35/L	TS5411-35/R
36	TS5411-36	TS5411-36/L	TS5411-36/R
37	TS5411-37	TS5411-37/L	TS5411-37/R
38	TS5411-38	TS5411-38/L	TS5411-38/R
39	TS5411-39	TS5411-39/L	TS5411-39/R
40	TS5411-40	TS5411-40/L	TS5411-40/R
41	TS5411-41	TS5411-41/L	TS5411-41/R
42	TS5411-42	TS5411-42/L	TS5411-42/R
43	TS5411-43	TS5411-43/L	TS5411-43/R
44	TS5411-44	TS5411-44/L	TS5411-44/R
45	TS5411-45	TS5411-45/L	TS5411-45/R
46	TS5411-46	TS5411-46/L	TS5411-46/R
47	TS5411-47	TS5411-47/L	TS5411-47/R
48	TS5411-48	TS5411-48/L	TS5411-48/R

PRODUCT INFORMATION

INDICATIONS

DRESDEN is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

- Risk group III:** condition after plantar ulcer
- Grade 0:** condition after lesion
- Risk group VII:** severe lesion
- Grade 1:** superficial wound
- Grade 2:** wounds up to tendon/capsule, plantar ulcer at the forefoot area
- Grade 4:** endemic necrosis at the forefoot area, condition after toe amputation
- Grade 4-5:** partial forefoot amputation, e.g. condition after single toe removals but preservation of the metatarsophalangeal joint of the big toe

SHOE TYPES

The asymmetric diabetic shoe DRESDEN is available in a STANDARD and WIDE version.

SOLE DESIGN


DRESDEN has a 30mm thick, rigid and extra strong rocker bottom sole with pitch. The rigid sole prevents flexion of the metatarsophalangeal joints and ensures the patient's foot being stabilised and immobilised. That is how an existent wound is relieved and the wound's healing supported in an optimal way.

INSOLES

The shoe is delivered without insoles. As interim solution, you could use an insole set with three insoles with different shore hardness (STANDARD F00033-XX and WIDE F00034-XX) or an insert blank (STANDARD F00021-XX and WIDE F00015-XX). For a successful treatment, in the long run, the shoe should be used with a diabetic foot orthotic.

SPECIAL FEATURES

DRESDEN, especially the WIDE shoe, provides much space for specific demands. The firm outer fabric can be separated from the soft inner fabric in order to avoid possible pressure points and prevent the manifestation of oedema. Furthermore, DRESDEN is an appropriate **TwinShoe** and thus height compensation for the diabetic shoe WETZLAR.





KASSEL STANDARD

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5421-35	TS5421-35/L	TS5421-35/R
36	TS5421-36	TS5421-36/L	TS5421-36/R
37	TS5421-37	TS5421-37/L	TS5421-37/R
38	TS5421-38	TS5421-38/L	TS5421-38/R
39	TS5421-39	TS5421-39/L	TS5421-39/R
40	TS5421-40	TS5421-40/L	TS5421-40/R
41	TS5421-41	TS5421-41/L	TS5421-41/R
42	TS5421-42	TS5421-42/L	TS5421-42/R
43	TS5421-43	TS5421-43/L	TS5421-43/R
44	TS5421-44	TS5421-44/L	TS5421-44/R
45	TS5421-45	TS5421-45/L	TS5421-45/R
46	TS5421-46	TS5421-46/L	TS5421-46/R
47	TS5421-47	TS5421-47/L	TS5421-47/R
48	TS5421-48	TS5421-48/L	TS5421-48/R

KASSEL WIDE

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5321-35	TS5321-35/L	TS5321-35/R
36	TS5321-36	TS5321-36/L	TS5321-36/R
37	TS5321-37	TS5321-37/L	TS5321-37/R
38	TS5321-38	TS5321-38/L	TS5321-38/R
39	TS5321-39	TS5321-39/L	TS5321-39/R
40	TS5321-40	TS5321-40/L	TS5321-40/R
41	TS5321-41	TS5321-41/L	TS5321-41/R
42	TS5321-42	TS5321-42/L	TS5321-42/R
43	TS5321-43	TS5321-43/L	TS5321-43/R
44	TS5321-44	TS5321-44/L	TS5321-44/R
45	TS5321-45	TS5321-45/L	TS5321-45/R
46	TS5321-46	TS5321-46/L	TS5321-46/R
47	TS5321-47	TS5321-47/L	TS5321-47/R
48	TS5321-48	TS5321-48/L	TS5321-48/R

PRODUCT INFORMATION

INDICATIONS

KASSEL is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

- Risk group III:** condition after plantar ulcer
- Grade 0:** condition after lesion
- Risk group VII:** severe lesion
- Grade 1:** superficial wound
- Grade 2:** wounds up to tendon/capsule, plantar ulcer at the forefoot area
- Grade 4:** endemic necrosis at the forefoot area, condition after toe amputation
- Grade 4-5:** partial forefoot amputation, e.g. condition after single toe removals but preservation of the metatarsophalangeal joint of the big toe

SHOE TYPES

The asymmetric diabetic shoe KASSEL is available in a STANDARD and WIDE version.

SOLE DESIGN


KASSEL has a 30mm thick, rigid and extra strong rocker bottom sole with pitch. The rigid sole prevents flexion of the metatarsophalangeal joints and ensures the patient's foot being stabilised and immobilised. That is how an existent wound is relieved and the wound's healing supported in an optimal way.

INSOLES

The shoe is delivered without insoles. As interim solution, you could use an insole set with three insoles with different shore hardness (STANDARD F00033-XX and WIDE F00034-XX) or an insert blank (STANDARD F00021-XX and WIDE F00015-XX). For a successful treatment, in the long run, the shoe should be used with a diabetic foot orthotic.

SPECIAL FEATURES

KASSEL, especially the WIDE shoe, provides much space for specific demands. The firm outer fabric can be separated from the soft inner fabric in order to avoid possible pressure points and prevent the manifestation of oedema. If applicable, a protective forefoot cap can be attached to the toe area that keeps the foot dry and warm.





PRODUCT INFORMATION

INDICATIONS

KONSTANZ is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

Risk group III: condition after plantar ulcer
Grade 0: condition after lesion, vulnerable foot

SOLE DESIGN

KONSTANZ has a 30mm thick, flexible sole with pitch. This sole design supports a natural heel to toe movement while walking.

INSOLES

Inside the shoe you will find three 5mm thick insoles made of EVA material (F00053-XX). The top insole has an antimicrobial coating. You can use these insoles as an interim solution for pressure relief until replacing them by a diabetic foot orthotic.

SPECIAL FEATURES

KONSTANZ is a light shoe which provides much air and space. Due to the open design and by using the comfortable hook and loop fasteners, this shoe can easily be put on and taken off.



ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TH4562-35	TH4562-35/L	TH4562-35/R	43	TH4562-43	TH4562-43/L	TH4562-43/R
36	TH4562-36	TH4562-36/L	TH4562-36/R	44	TH4562-44	TH4562-44/L	TH4562-44/R
37	TH4562-37	TH4562-37/L	TH4562-37/R	45	TH4562-45	TH4562-45/L	TH4562-45/R
38	TH4562-38	TH4562-38/L	TH4562-38/R	46	TH4562-46	TH4562-46/L	TH4562-46/R
39	TH4562-39	TH4562-39/L	TH4562-39/R	47	TH4562-47	TH4562-47/L	TH4562-47/R
40	TH4562-40	TH4562-40/L	TH4562-40/R	48	TH4562-48	TH4562-48/L	TH4562-48/R
41	TH4562-41	TH4562-41/L	TH4562-41/R	49	TH4562-49	TH4562-49/L	TH4562-49/R
42	TH4562-42	TH4562-42/L	TH4562-42/R	50	TH4562-50	TH4562-50/L	TH4562-50/R



PRODUCT INFORMATION

INDICATIONS

WEIMAR is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

- Risk group III: condition after plantar ulcer
- Grade 0: condition after lesion, vulnerable foot

SOLE DESIGN

WEIMAR has a 30mm thick, flexible sole with pitch. This sole design supports a natural heel to toe movement while walking.

INSOLES

Inside the shoe you will find three 5mm thick insoles made of EVA material (F00053-XX). The top insole has an antimicrobial coating. You can use these insoles as an interim solution for pressure relief until replacing them by a diabetic foot orthotic.


SPECIAL FEATURES

WEIMAR is a light shoe which - despite its sandal shape - has a toe cap to protect the toe area (e.g. against bumps).

ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS4502-35	TS4502-35/L	TS4502-35/R	43	TS4502-43	TS4502-43/L	TS4502-43/R
36	TS4502-36	TS4502-36/L	TS4502-36/R	44	TS4502-44	TS4502-44/L	TS4502-44/R
37	TS4502-37	TS4502-37/L	TS4502-37/R	45	TS4502-45	TS4502-45/L	TS4502-45/R
38	TS4502-38	TS4502-38/L	TS4502-38/R	46	TS4502-46	TS4502-46/L	TS4502-46/R
39	TS4502-39	TS4502-39/L	TS4502-39/R	47	TS4502-47	TS4502-47/L	TS4502-47/R
40	TS4502-40	TS4502-40/L	TS4502-40/R	48	TS4502-48	TS4502-48/L	TS4502-48/R
41	TS4502-41	TS4502-41/L	TS4502-41/R	49	TS4502-49	TS4502-49/L	TS4502-49/R
42	TS4502-42	TS4502-42/L	TS4502-42/R	50	TS4502-50	TS4502-50/L	TS4502-50/R





PRODUCT INFORMATION

INDICATIONS

HEIDELBERG is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

Risk group III: condition after plantar ulcer
Grade 0: condition after lesion, vulnerable foot

SOLE DESIGN

HEIDELBERG has a 30mm thick, flexible sole with pitch. This sole design supports a natural heel to toe movement while walking.

INSOLES

Inside the shoe you will find three 5mm thick insoles made of EVA material (F00053-XX). The top insole has an antimicrobial coating. You can use these insoles as an interim solution for pressure relief until replacing them by a diabetic foot orthotic.

SPECIAL FEATURES

HEIDELBERG provides a variable volume. The partially flexible upper fabric of this shoe is elastic and adapts to different demands on space (e.g. due to varying or thick bandages).

ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS4552-35	TS4552-35/L	TS4552-35/R	43	TS4552-43	TS4552-43/L	TS4552-43/R
36	TS4552-36	TS4552-36/L	TS4552-36/R	44	TS4552-44	TS4552-44/L	TS4552-44/R
37	TS4552-37	TS4552-37/L	TS4552-37/R	45	TS4552-45	TS4552-45/L	TS4552-45/R
38	TS4552-38	TS4552-38/L	TS4552-38/R	46	TS4552-46	TS4552-46/L	TS4552-46/R
39	TS4552-39	TS4552-39/L	TS4552-39/R	47	TS4552-47	TS4552-47/L	TS4552-47/R
40	TS4552-40	TS4552-40/L	TS4552-40/R	48	TS4552-48	TS4552-48/L	TS4552-48/R
41	TS4552-41	TS4552-41/L	TS4552-41/R	49	TS4552-49	TS4552-49/L	TS4552-49/R
42	TS4552-42	TS4552-42/L	TS4552-42/R	50	TS4552-50	TS4552-50/L	TS4552-50/R



NÜRNBERG STANDARD

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS4311-35	TS4311-35/L	TS4311-35/R
36	TS4311-36	TS4311-36/L	TS4311-36/R
37	TS4311-37	TS4311-37/L	TS4311-37/R
38	TS4311-38	TS4311-38/L	TS4311-38/R
39	TS4311-39	TS4311-39/L	TS4311-39/R
40	TS4311-40	TS4311-40/L	TS4311-40/R
41	TS4311-41	TS4311-41/L	TS4311-41/R
42	TS4311-42	TS4311-42/L	TS4311-42/R
43	TS4311-43	TS4311-43/L	TS4311-43/R
44	TS4311-44	TS4311-44/L	TS4311-44/R
45	TS4311-45	TS4311-45/L	TS4311-45/R
46	TS4311-46	TS4311-46/L	TS4311-46/R
47	TS4311-47	TS4311-47/L	TS4311-47/R
48	TS4311-48	TS4311-48/L	TS4311-48/R

NÜRNBERG WIDE

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS4411-35	TS4411-35/L	TS4411-35/R
36	TS4411-36	TS4411-36/L	TS4411-36/R
37	TS4411-37	TS4411-37/L	TS4411-37/R
38	TS4411-38	TS4411-38/L	TS4411-38/R
39	TS4411-39	TS4411-39/L	TS4411-39/R
40	TS4411-40	TS4411-40/L	TS4411-40/R
41	TS4411-41	TS4411-41/L	TS4411-41/R
42	TS4411-42	TS4411-42/L	TS4411-42/R
43	TS4411-43	TS4411-43/L	TS4411-43/R
44	TS4411-44	TS4411-44/L	TS4411-44/R
45	TS4411-45	TS4411-45/L	TS4411-45/R
46	TS4411-46	TS4411-46/L	TS4411-46/R
47	TS4411-47	TS4411-47/L	TS4411-47/R
48	TS4411-48	TS4411-48/L	TS4411-48/R

PRODUCT INFORMATION

INDICATIONS

NÜRNBERG is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

Risk group III: condition after plantar ulcer
Grade 0: condition after lesion (vulnerable foot)

SHOE TYPES

The asymmetric diabetic shoe NÜRNBERG is available in a STANDARD and WIDE version.

SOLE DESIGN


NÜRNBERG has a 30mm thick, rigid and extra strong rocker bottom sole with pitch. The rigid sole prevents flexion of the metatarsophalangeal joints and ensures the patient's foot being stabilised and immobilised. That is how an existent wound is relieved and the wound's healing supported in an optimal way.

INSOLES

The shoe is delivered without insoles. As interim solution, you could use an insole set with three insoles with different shore hardness (STANDARD F00033-XX and WIDE F00034-XX) or an insert blank (STANDARD F00021-XX and WIDE F00015-XX). For a successful treatment, in the long run, the shoe should be used with a diabetic foot orthotic.

SPECIAL FEATURES

NÜRNBERG, especially the WIDE shoe, provides much space for specific demands. Besides, the firm outer fabric can be separated from the soft inner fabric in order to avoid possible pressure points and prevent the manifestation of oedema.





PRODUCT INFORMATION

INDICATIONS

POTSDAM is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

Risk group III: condition after plantar ulcer
Grade 0: condition after lesion (vulnerable foot)

SOLE DESIGN


POTSDAM has a 30mm thick, rigid and extra strong rocker bottom sole with pitch. The rigid sole prevents flexion of the metatarsophalangeal joints and ensures the patient's foot being stabilised and immobilised. That is how an existent wound is relieved and the wound's healing supported in an optimal way.

INSOLES

Inside the shoe you will find three 5mm thick insoles made of EVA material (F00053-XX). The top insole has an antimicrobial coating. You can use these insoles as an interim solution for pressure relief until replacing them by a diabetic foot orthotic.


SPECIAL FEATURES

POTSDAM is a therapeutic slipper and the ideal shoe to easily slip into for short distances.



ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS4391-35	TS4391-35/L	TS4391-35/R	43	TS4391-43	TS4391-43/L	TS4391-43/R
36	TS4391-36	TS4391-36/L	TS4391-36/R	44	TS4391-44	TS4391-44/L	TS4391-44/R
37	TS4391-37	TS4391-37/L	TS4391-37/R	45	TS4391-45	TS4391-45/L	TS4391-45/R
38	TS4391-38	TS4391-38/L	TS4391-38/R	46	TS4391-46	TS4391-46/L	TS4391-46/R
39	TS4391-39	TS4391-39/L	TS4391-39/R	47	TS4391-47	TS4391-47/L	TS4391-47/R
40	TS4391-40	TS4391-40/L	TS4391-40/R	48	TS4391-48	TS4391-48/L	TS4391-48/R
41	TS4391-41	TS4391-41/L	TS4391-41/R	49	TS4391-49	TS4391-49/L	TS4391-49/R
42	TS4391-42	TS4391-42/L	TS4391-42/R	50	TS4391-50	TS4391-50/L	TS4391-50/R



PRODUCT INFORMATION

INDICATIONS

WETZLAR is particularly suited for the following risk groups and wound stages of the diabetic foot syndrome:

- Risk group VII: severe lesion
- Grade 1: superficial wound
- Grade 2: wounds up to tendon/capsule, plantar ulcer at the forefoot area

SHOE TYPES

The asymmetric diabetic shoe WETZLAR is available in a left or right version.

SOLE DESIGN


WETZLAR has a 30mm thick, rigid and extra strong rocker bottom sole with pitch. The rigid sole prevents flexion of the metatarsophalangeal joints and ensures the patient's foot being stabilised and immobilised. That is how an existent wound is relieved and the wound's healing supported in an optimal way.

INSOLES

The shoe is delivered without insoles. As interim solution, you could use an insole set with three insoles with different shore hardness (F00034-XX) or an insert blank (F00015-XX). For a successful treatment, in the long run, the shoe should be used with a diabetic foot orthotic.

SPECIAL FEATURES

WETZLAR provides an open heel area. You can attach a protective heel cap to this area which is included in the scope of delivery. The shoe has no midsole in the heel area, thus, the sole can directly be ground from the inside. The low framing brim keeps a subsequently inserted (diabetic) insole in position and provides additional seams where the fabric can be neatly cut if pressure relief should be necessary.



ORDER INFORMATION

Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe	Shoe size	Art. no. pair	Art. no. left shoe	Art. no. right shoe
35	TS5431-35	TS5431-35/L	TS5431-35/R	42	TS5431-42	TS5431-42/L	TS5431-42/R
36	TS5431-36	TS5431-36/L	TS5431-36/R	43	TS5431-43	TS5431-43/L	TS5431-43/R
37	TS5431-37	TS5431-37/L	TS5431-37/R	44	TS5431-44	TS5431-44/L	TS5431-44/R
38	TS5431-38	TS5431-38/L	TS5431-38/R	45	TS5431-45	TS5431-45/L	TS5431-45/R
39	TS5431-39	TS5431-39/L	TS5431-39/R	46	TS5431-46	TS5431-46/L	TS5431-46/R
40	TS5431-40	TS5431-40/L	TS5431-40/R	47	TS5431-47	TS5431-47/L	TS5431-47/R
41	TS5431-41	TS5431-41/L	TS5431-41/R	48	TS5431-48	TS5431-48/L	TS5431-48/R

INNER SHOE DIMENSIONS

Dear Customer,

We measured our therapeutic shoes according to the Mondopoint system. Afterwards, the shoe sizes were recalculated by means of the Paris Point.

Mondopoint System

The Mondopoint system is an international standardised system to indicate shoe sizes. In contrast to other systems, not only the inner shoe length but also the inner shoe width is measured.

Paris Point


To determine the shoe sizes we used the Paris Point which is calculated as follows:

$$\text{shoe size} = \text{inner shoe length in cm} \times 1.5$$

Length and width are indicated in mm.

ATTENTION!

In this catalogue, we state outdoor shoe sizes. These shoe sizes do not cover additional space for bandages or the like.




ASYMMETRIC SHOES

Inner shoe length

The inner shoe length (c) was measured from the middle of the heel area to the middle of the great toe with a standard measurement device.

Inner shoe widths (ball and heel area)

The inner shoe widths were measured at the insole. We measured from the toe-cap and drew a crossline at 1/3 of the insole for the inner shoe width in ball area. At this right-angled straight line a 15° angle was measured. The resulting line was where we measured the inner shoe width in ball area (a). The inner shoe width in heel area (b) was measured at the level of the middle of the heel.



INNER SHOE DIMENSIONS STANDARD (F00033 and F00053)

Shoe size	Inner shoe dimensions [mm]			Shoe size	Inner shoe dimensions [mm]		
	a	b	c		a	b	c
35	93	66	231	43	106	75	288
36	94	67	237	44	109	76	295
37	95	69	245	45	109	79	299
38	98	69	250	46	110	79	305
39	99	71	258	47	111	81	313
40	101	73	266	48	113	81	318
41	104	74	271	49	115	81	325
42	105	75	277	50	115	81	332

AUGSBURG STANDARD, DRESDEN STANDARD, FLENSBURG STANDARD, HEIDELBERG, KASSEL STANDARD, KONSTANZ, LEIPZIG, NÜRNBERG STANDARD, POTSDAM, WEIMAR

INNER SHOE DIMENSIONS WIDE (F00034 and F00054)

Shoe size	Inner shoe dimensions [mm]			Shoe size	Inner shoe dimensions [mm]		
	a	b	c		a	b	c
35	102	67	233	43	114	76	282
36	103	69	240	44	116	79	287
37	104	71	247	45	118	80	295
38	107	73	249	46	121	81	302
39	109	73	255	47	123	82	308
40	111	74	263	48	126	83	319
41	111	75	269	49	126	83	325
42	112	76	277	50	127	87	333

AUGSBURG WIDE, DRESDEN WIDE, FLENSBURG WIDE, KASSEL WIDE, NÜRNBERG WIDE, WETZLAR

SPARE AND ACCESSORY PARTS



The top insole is antimicrobial. The insoles are made of EVA material and provide a different shore hardness: grey = 46°, white = 23°, antimicrobial = 28°

INSOLE SET STANDARD

Shoe size	Art. no. set	Shoe size	Art. no. set
35	F00053-35	43	F00053-43
36	F00053-36	44	F00053-44
37	F00053-37	45	F00053-45
38	F00053-38	46	F00053-46
39	F00053-39	47	F00053-47
40	F00053-40	48	F00053-48
41	F00053-41	49	F00053-49
42	F00053-42	50	F00053-50

AUGSBURG STANDARD, DRESDEN STANDARD, FLENSBURG STANDARD, HEIDELBERG, KASSEL STANDARD, KONSTANZ, LEIPZIG, NÜRNBERG STANDARD, POTSDAM, WEIMAR

INSOLE SET WIDE

Shoe size	Art. no. set	Shoe size	Art. no. set
35	F00054-35	43	F00054-43
36	F00054-36	44	F00054-44
37	F00054-37	45	F00054-45
38	F00054-38	46	F00054-46
39	F00054-39	47	F00054-47
40	F00054-40	48	F00054-48
41	F00054-41	49	F00054-49
42	F00054-42	50	F00054-50

AUGSBURG WIDE, DRESDEN WIDE, FLENSBURG WIDE, KASSEL WIDE, NÜRNBERG WIDE, WETZLAR

INSOLE SET STANDARD


Shoe size	Art. no. set	Shoe size	Art. no. set
35	F00033-35	43	F00033-43
36	F00033-36	44	F00033-44
37	F00033-37	45	F00033-45
38	F00033-38	46	F00033-46
39	F00033-39	47	F00033-47
40	F00033-40	48	F00033-48
41	F00033-41	49	F00033-49
42	F00033-42	50	F00033-50

AUGSBURG STANDARD, DRESDEN STANDARD, FLENSBURG STANDARD, HEIDELBERG, KASSEL STANDARD, KONSTANZ, LEIPZIG, NÜRNBERG STANDARD, POTSDAM

INSOLE SET WIDE

Shoe size	Art. no. set	Shoe size	Art. no. set
35	F00034-35	43	F00034-43
36	F00034-36	44	F00034-44
37	F00034-37	45	F00034-45
38	F00034-38	46	F00034-46
39	F00034-39	47	F00034-47
40	F00034-40	48	F00034-48
41	F00034-41	49	F00034-49
42	F00034-42	50	F00034-50

AUGSBURG WIDE, DRESDEN WIDE, FLENSBURG WIDE, KASSEL WIDE, NÜRNBERG WIDE, WETZLAR



The insert blank is made of EVA material and provides the following shore hardness = 60°

INSERT BLANK STANDARD

Shoe size	Art. no. left shoe	Art. no. right shoe	Unit
35	F00021-35/L	F00021-35/R	Piece
36	F00021-36/L	F00021-36/R	Piece
37	F00021-37/L	F00021-37/R	Piece
38	F00021-38/L	F00021-38/R	Piece
39	F00021-39/L	F00021-39/R	Piece
40	F00021-40/L	F00021-40/R	Piece
41	F00021-41/L	F00021-41/R	Piece
42	F00021-42/L	F00021-42/R	Piece
43	F00021-43/L	F00021-43/R	Piece
44	F00021-44/L	F00021-44/R	Piece
45	F00021-45/L	F00021-45/R	Piece
46	F00021-46/L	F00021-46/R	Piece
47	F00021-47/L	F00021-47/R	Piece
48	F00021-48/L	F00021-48/R	Piece

DRESDEN STANDARD, KASSEL STANDARD, NÜRNBERG STANDARD

INSERT BLANK WIDE

Shoe size	Art. no. left shoe	Art. no. right shoe	Unit
35	F00015-35/L	F00015-35/R	Piece
36	F00015-36/L	F00015-36/R	Piece
37	F00015-37/L	F00015-37/R	Piece
38	F00015-38/L	F00015-38/R	Piece
39	F00015-39/L	F00015-39/R	Piece
40	F00015-40/L	F00015-40/R	Piece
41	F00015-41/L	F00015-41/R	Piece
42	F00015-42/L	F00015-42/R	Piece
43	F00015-43/L	F00015-43/R	Piece
44	F00015-44/L	F00015-44/R	Piece
45	F00015-45/L	F00015-45/R	Piece
46	F00015-46/L	F00015-46/R	Piece
47	F00015-47/L	F00015-47/R	Piece
48	F00015-48/L	F00015-48/R	Piece

DRESDEN WIDE, KASSEL WIDE, NÜRNBERG WIDE, WETZLAR



The insoles are made of EVA material and provide a different shore hardness: black = 40°, grey = 25°, white = 50°



OUTSOLE MATERIAL, RUBBER

Art. no.	Unit	Dimensions [mm]
F00092-S	Sheet	220 x 110 x 1.5
F00092-M	Sheet	295 x 145 x 1.5
F00092-L	Sheet	357 x 158 x 1.5

AUGSBURG, FLENSBURG, POTSDAM



PROTECTIVE FOREFOOT CAP

Shoe size	Art. no.
35-36	F09520-XXS
37-38	F09520-XS
39-40	F09520-S
41-42	F09520-M
43-44	F09520-L
45-46	F09520-XL
47-48	F09520-XXL

The protective forefoot cap can easily be attached to the KASSEL STANDARD and WIDE shoe and keeps the foot dry and warm.



PROTECTIVE FOREFOOT CAP

Shoe size	Art. no.
35-36	F09580-XXS
37-38	F09580-XS
39-40	F09580-S
41-42	F09580-M
43-44	F09580-L
45-46	F09580-XL
47-48	F09580-XXL
49-50	F09580-XXXL

The protective forefoot cap can easily be attached to the AUGSBURG STANDARD and WIDE shoe and keeps the foot dry and warm.



PROTECTIVE HEEL CAP

Shoe size	Art. no.
35-36	F09543-XXS
37-38	F09543-XS
39-40	F09543-S
41-42	F09543-M
43-44	F09543-L
45-46	F09543-XL
47-48	F09543-XXL

The protective heel cap can easily be attached to the WETZLAR shoe and keeps the foot dry and warm.

THERAPEUTIC SHOES
POST-OP



You will find our therapeutic shoe catalogue POST-OP on our website in the download area:
www.fior-gentz.com/download

THERAPEUTIC SHOES
HOMECARE



You will find our therapeutic shoe catalogue HOMECARE on our website in the download area:
www.fior-gentz.com/download

The information in this catalogue is valid at the date of printing. The contained product information serves as a guideline. Subject to technical modifications. All rights, particularly the distribution, copy and translation of this manual or any part of it, in paper or as electronic document, must be authorised in writing by FIOR & GENTZ Gesellschaft für Entwicklung und Vertrieb von orthopädiotechnischen Systemen mbH. Reprints, copies and any other electronic reproduction, even partial, must be authorised in writing by FIOR & GENTZ Gesellschaft für Entwicklung und Vertrieb von orthopädiotechnischen Systemen mbH.



FIOR & GENTZ

Gesellschaft für Entwicklung und Vertrieb
von orthopädietechnischen Systemen mbH

Dorette-von-Stern-Straße 5
21337 Lüneburg (Germany)

☎ +49 4131 24445-0
☎ +49 4131 24445-57

✉ info@fior-gentz.de
🌐 www.fior-gentz.com

FIOR & GENTZ
ORTHOPÄDIETECHNIK MIT SYSTEM