

Uniprox. Quality and Function

F90 Uniprox Multiaxial Carbonfiber Foot



For a softly and flexible movement

The special combination of carbon and plastic spring elements as well as elastomer bumpers of the **F90 Uniprox Multi-axial Carbonfiber Foot** offer the user a smooth rollover. The excellent flexibility of the multiaxial ankle joint allows the foot to adapt to any kind of ground conditions. The elastomer bumpers can be individually adapted to the user's needs. A complete bumper set is included.

Features:

- Available Size: 24 - 29 cm
- Heel Height: 10 mm
- Weight Limit: < 125 kg
- Mobility Level: 2, 3, 4
- Installation Height: 135 mm
- Weighth: 810 g (Size 27)

Selection of bumpers

Weight Activity Level	41-68 kg			69 - 76 kg			77 - 95 kg			96 -112 kg			113 - 125 kg		
	Anterior	Axial	Posterior	Anterior	Axial	Posterior	Anterior	Axial	Posterior	Anterior	Axial	Posterior	Anterior	Axial	Posterior
AL 2	85	85	55	95	95	65	95	95	75	95	95	75	95	95	75
AL 3	85	85	65	95	95	75	95	95	75	95	95	75	95	95	85
AL 4	85	85	65	95	95	75	95	95	75	105	95	85	105	95	85

Bumper Anterior	Bumper Axial	Bumper Posterior
105 Shore A	95 Shore A	85 Shore A
95 Shore A	85 Shore A	75 Shore A
85 Shore A		65 Shore A
		55 Shore A



Accessories and Spare Parts for F90 Uniprox Multiaxial Carbonfiber-Foot

Bumper

Order No.	Description	Article No.
E-F80-PV2	Bumper anterior 85 shore	4 602 029 00 02 000
E-F80-PV3	Bumper anterior 95 shore	4 602 029 00 03 000
E-F80-PV4	Bumper anterior 105 shore	4 602 029 00 04 000
E-F80-PM1	Bumper axial 85 shore	4 602 029 00 05 000
E-F80-PM2	Bumper axial 95 shore	4 602 029 00 06 000
E-F80-PH1	Bumper posterior 55 shore	4 602 029 00 16 000
E-F80-PH2	Bumper posterior 65 shore	4 602 029 00 17 000
E-F80-PH3	Bumper posterior 75 shore	4 602 029 00 18 000
E-F80-PH4	Bumper posterior 85 shore	4 602 029 00 19 000
E-F80-PS	Bumper Set for F80 and F90	4 602 029 00 28 000

Accessories

Order No.	Description	Article No.
E-F80-OR	O-Ring	4 602 029 00 00 001
E-F80-SS	Spectra-Sock for F80	4 602 029 00 39 000
MM202	Lithium grease	4 690 240 00 89 100

Foot shell F90 Uniprox Carbon Foot and connection plate



Order-No.	Size / Side	Description	Available Sizes	Article No.
E-F80-	<input type="checkbox"/> <input type="checkbox"/>	Foot shell	24-29	4 603 029 0X 00 0ZZ
E-F80-AP1		Connection plate, small	24-26	4 602 029 00 00 002
E-F80-AP2		Connection plate, big	27-29	4 602 029 00 00 003