

Sparkle

STOCKINETTES & BRAIDS

by **LimbTex**
PROSTHETIC & ORTHOTIC



FIRST EDITION



FIESTA RED

Code: SPARK_ _FRED



GOLD

Code: SPARK_ _GLD



GUN METAL

Code: SPARK_ _GUNM



JET BLUE

Code: SPARK_ _JBLU



PINK

Code: SPARK_ _PINK



RAINBOW

Code: SPARK_ _RAINB



SILVER

Code: SPARK_ _SIL



TURQUOISE

Code: SPARK_ _TURQ



RED MULTI

Code: SPARK_ _REDMULTI

Sparkle

STOCKINETTES & BRAIDS

by **LimbTex**
PROSTHETIC & ORTHOTIC

This is our unique range of stockinettes and braids made with the addition of Lurex® yarns and film for the “Sparkle effect”. The stockinettes have been burst tested by an independent test house and have been proven to have comparable strength to Perlon.

What this means is that unlike other fabric finishes, the inner or outer double Perlon on a prosthetic socket can be replaced with a double Sparkle layer which adds no weight, no extra thickness and the manufacturing process remains the same. These are an excellent addition to other available finishes that we offer which are easy to use and cost effective.

CODE (Insert required width e.g. SPARK10FRED)	Available Widths in cm	Description	Length of 6, 8, 10cm Width Stockinette	Length of 12. 15. 20cm Width Stockinette
SPARK__FRED	6-8-10-12-15-20	Nylon/Lurex® Fiesta Red Stockinette	2m	2.5m
SPARK__GLD	6-8-10-12-15-20	Nylon/Lurex® Gold Stockinette	2m	2.5m
SPARK__GUNM	6-8-10-12-15-20	Nylon/Lurex® Gun Metal Stockinette	2m	2.5m
SPARK__JBLU	6-8-10-12-15-20	Nylon/Lurex® Jet Blue Stockinette	2m	2.5m
SPARK__PINK	6-8-10-12-15-20	Nylon/Lurex® Pink Stockinette	2m	2.5m
SPARK__RAINB	6-8-10-12-15-20	Nylon/Lurex® Rainbow Stockinette	2m	2.5m
SPARK__SIL	6-8-10-12-15-20	Nylon/Lurex® Silver Stockinette	2m	2.5m
SPARK__TURQ	6-8-10-12-15-20	Nylon/Lurex® Turquoise Stockinette	2m	2.5m
SPARK__REDMULTI	6-8-10-12-15-20	Nylon/Lurex® Red Multi Stockinette	2m	2.5m



LimbTex

PROSTHETIC & ORTHOTIC

BASALT FIBRE & LUREX® BRAIDS

This is a decorative structural material made from Basalt fibres braided with bright and shiny Lurex® film. It can be used to create both inside and outside socket finishes. The braids have been developed to give a consistent diamond pattern which keep their shine once laminated giving a pleasingly different look. It is recommended to use a backing layer of either standard Basalt braid or NYLONBLK to keep the Lurex® defined.



Available in Braids only

Code	Description	Weight g/m
BFBSL130/5	130mm Ø +/- 45°Basalt Fibre & Lurex® Braid, 5m Roll	185
BFBSL150/5	150mm Ø +/- 45°Basalt Fibre & Lurex® Braid, 5m Roll	370
BFBSL225/5	225mm Ø +/- 45°Basalt Fibre & Lurex® Braid, 5m Roll	461

LimbTex Ltd

Units 1-3, Elizabeth Business Park,
Tigers Close, South Wigston, Leicester LE18 4AE

Tel: +44 (0) 116 2785440 | Fax: +44 (0) 116 2780227 | Email: info@limbtex.com | Web: www.limbtex.com