

Lightweight blade, ideal for moderate impact sports and activities.

Features

- Small sizes available
- e-carbon foot with independent springs
- Dynamic pylon enhances shock absorption
- Easy to finish cosmetically
- Sandal Toe Footshell

Foot shell and glide sock included.

Specifications

Max. User Weight:	166kg
Activity Level:	(2), 3, (4*)
Size Range:	22cm-30cm
Component Weight:	444g [†]
Build Height:	200mm - 365mm (min - max)
Heel Height:	10mm



Suitable for submersion

Order Example




Product Code	Size	Side	Width*	Spring Set	Sandal Toe	+ Top Housing
JAV	25	L	N	3	S	Adapter

*Narrow (N) and Wide (W) available for sizes 25-27 only.

For dark tone add suffix D.

Example: foot size 25, left, narrow, spring rating 3, sandal toe, adapter.


Accessories

Foot Shell Cover	539087	
Top housing male	539601	
Top housing female	539602	

Javelin



Build Height



Activity Level 2 users who would benefit from this foot will require softer springs as appropriate for the individual. Refer to Instructions for Use for further information.

Spring set recommendations are for trans-tibial users. For trans-femoral we suggest selecting a spring set one level lower.

Activity	User Weight									
	44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166 kg
	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365 lbs
3	1	1	2	3	4	5	6	7	8	9



*Maximum user weight 100 kg and always use one higher spring rate category than shown in Spring selection table.

[†]Component weight shown is for a size 26cm without foot shell.