



The ALPS Skin Reliever eliminates shear forces, friction, and prevents abrasions on the skin of the stump. It also accommodates volume fluctuations over the life of the prosthesis. There is no need to make a new socket, ALPS skin Reliever can be worn under silicone, gel, and pelite liners.

Features and Applications:

- Reduces shearing and abrasions
- High elasticity
- Replaces pelite and socks
- Softens and soothes the skin

PRODUCT NO.


Description

ENCP-XX

Skin Reliever Tapered from 3mm to 1.5mm
Size: 16-20-28

XX indicates size. For details please see the sizing chart in the back page.


Skin Reliever ENCP. Cross section
3 mm - 1.5mm Tapered - No Fabric or Distal Attachment



Skin Reliever ENCP. Sizing chart
Measurements taken 6 cm from the distal end

Measured size	3/1,5 mm
10-19 cm	ENCP16
20-27 cm	ENCP20
28-37 cm	ENCP28

How to size
Measurements taken 6 cm from the distal end



ALPS SOUTH



www.easyliner.eu

+47 55 91 88 60
post@ortopro.no

ortopro.no
fOrtopro

ORTOPRO

3250 01108.1E3 B