



FIXQUICK POWERED BY
BOA FIT SYSTEM

Semi-rigid orthosis line with **BOA**® closure

PRESS, TURN AND READY!





ORLIMAN®

FIXQUICK POWERED BY
BOA FIT SYSTEM

The **BOA®** closure and adjustment system is reinventing the performance of adjustment systems for medical devices, sports equipment, footwear, etc. as an innovation over traditional closure systems.

More than **25 million** products around the world already use this innovative system.



Orliman wanted to join this revolution with

FIXQUICK POWERED BY BOA FIT SYSTEM

a whole new concept in orthopaedics aimed at providing quick, comfortable and effective solutions through design and technology with a closure system that offers a perfect fit just by turning a wheel.



RELIABILITY

The three components of the **BOA**® system have been designed to ensure Orliman Fixquick orthoses are multi-purpose and can be used in any environment

- 1 THE DIAL**
Provides a strong hold and a close fit, down to the millimetre.
- 2 THE LACES**
With stainless steel lining, they are ultra-resistant and light.
- 3 THE GUIDES**
They ensure there is no friction between the laces and offer a completely balanced fit without pressure points.

PRECISION, SIMPLICITY, LIGHTNESS AND COMFORT.

A whole new line of cutting-edge orthopaedic products designed above all for the click generation who want fast and effective results and for those looking for a difference without losing any of the essence and functionality of orthopaedic and sports products.

PRESS, TURN AND READY!

Fast and easy to put on, **Fixquick** meets these requirements: press, turn and ready!

YOUR DAY TO DAY, WITH ORLIMAN FIXQUICK

POWERED BY

BOA® FIT SYSTEM

Whether due to overuse, repetitive actions, deformations, sports injuries or accidents, there are a number of pathologies that can affect the joints we use the most (wrist, knee, ankle) and, therefore, our daily lives.

A SYSTEM THAT ADAPTS TO YOUR PACE OF LIFE



This new line of products has been designed to treat the pathologies we suffer in our day to day such as overstrain, contusions, sprains, tendinitis, etc.





Main characteristics

SPEED SIMPLICITY

The Boa® closure system makes it easier to adjust Orliman® Fixquick orthoses so they can be put on quickly with a simplified hold that can be readjusted throughout the day. All you have to do is activate the system by pressing down on the dial and then turning it towards the largest [] arrows to get the best fit or in the direction of the smallest ones [] for a looser fit. To unlock it, just pull the disc up!



PRESS



TIGHTEN



LOOSEN



RELEASE

COMFORT

The laces on the Boa® closure system have a 100% stainless steel inner structure for immense flexibility as well as 100% polyamide lining so they may be slid along the guides without any friction for adjustment and without any pressure points.

Orliman® Fixquick products have high-quality finishes and are specially designed to ensure user wellbeing. The inner suede and elastic lining have been tested and are guaranteed not to irritate the skin.

PRECISION

The gradual cog adjustment mechanism allows micro-adjustments with millimetric precision:

1 click = 1 millimetre more of adjustment or looseness.

This precision allows total control when adjusting Orliman® Fixquick orthoses to get the desired pressure and comfort.

Ref.: **BCS50D** | RIGHT

Ref.: **BCS50I** | LEFT



SEMI-RIGID WRIST BRACE WITH A PALM SPLINT

Features

Wrist braces made of semi-rigid breathable fabric with the Boa® closure system for quick and easy adjustment. The gradual adjustment enables micro-adjustments with millimetric precision. It includes a malleable aluminium palm and/or thumb splint in accordance with the chosen model to perfectly secure the wrist or thumb based on the desired angles in addition to two rigid plastic dorsal reinforcements.

Recommended for

- Wrist trauma.
- Sprains.
- Contusions.
- Tendinitis.
- Carpal tunnel syndrome.
- Post-operative immobilisation.

Colour: Black.



SIZES	MEASUREMENTS
1	13.5 - 14.5 cm
2	14.5 - 16 cm
3	16 - 17.5 cm
4	17.5 - 19 cm

Ref.: **BCS51D** | RIGHT

Ref.: **BCS51I** | LEFT



SEMI-RIGID WRIST BRACE WITH A PALM/THUMB SPLINT

Padded lining for comfort and perfect fit. The dorsal side is reinforced with a lined tongue with a thermal plastic core. The elastic fabric makes it easy to put on.

Recommended for

- Wrist and thumb trauma.
- Wrist/thumb sprains.
- Scaphoid injuries.
- De Quervain syndrome.
- Acute rhizarthrosis.
- Post-operative immobilisation.

Colour: Black.



SIZES	MEASUREMENTS
1	13.5 - 14.5 cm
2	14.5 - 16 cm
3	16 - 17.5 cm
4	17.5 - 19 cm

Ref.: **BCS400**



TENNIS ELBOW BRACE


Features

A brace made of semi-rigid breathable fabric with the Boa® closure system for quick and easy adjustment. The gradual adjustment enables micro-adjustments with millimetric precision. Includes an adjustable and removable pressure pad. Padded lining for comfort and perfect adaptation.

Recommended for

- Epicondylitis (tennis elbow).
- Medial epicondylitis (golfer's elbow).
- The prevention of sports and occupational injuries.

Colour: Black.



3 cm from the elbow

SIZES	MEASUREMENTS
1	22 - 27 cm
2	27 - 32 cm

Ref.: **BCS900**



STABILISING ANKLE BRACE

Features

An ankle brace made of semi-rigid breathable fabric with the Boa® closure system for quick and easy fit. The gradual adjustment provides adequate protection and ankle stabilisation. Includes anatomically designed medial and lateral reinforcements for excellent ankle/foot eversion/inversion control. Padded lining for comfort and perfect adaptation.

Recommended for

- Ankle trauma.
- Sprains.
- Ligament instability.
- Post-operative recovery.
- Resumption of activities and prevention.

Colour: Black.



SIZES	MEASUREMENTS
1	19 - 22 cm
2	22 - 26 cm
3	26 - 29 cm
4	29 - 32 cm

PRESS, TURN AND READY!

Advertising Material



Catalogue

Ref.: **Cat#BOAI**

28x21 cm



Display

Ref.: **Displays-BOAI**

20x30 cm

