



***SOFY-PLANT***®

**MAIN FEATURES**




1 2 3

**MAXIMUM ELASTICITY silicone insoles**

The SOFY-PLANT silicone insoles unload the whole foot, especially the support areas of the metatarsals and the heel. Its two densities (one especially soft for the unloading areas) absorb shocks and impacts from the ground while walking, providing relief and comfort by avoiding pressure on the area of the heel and the metatarsals.

Especially indicated for tired, weak and delicate feet and in post-surgery situations. Also indicated for the practice of any sport that requires an extended effort.



Comparing figures A and B a significant reduction can be seen in the maximum pressure in the area of the heel and the top of the 5th metatarsal, moving from pressure maximums of 600Kpa to maximums of about 250 Kpa. SOFY-PLANT insoles and heelpieces are made of a notably pleasant material to touch and are made of 100% silicone, which makes them perfectly tolerable and easily washed with cold water and without any type of soap.



**Ref.: TL-601**

**SILICONE HEEL CUP**

SIZES	FOOTWEAR SIZE
1	35-38
2	39-42
3	43-46



**Units for packaging:** 1 pair.

**Characteristics:**

Due to their visco-elastic properties, these silicone heel cups cushion the load points both during walking and in static position.

**Indications:**

- Relieve talalgia.
- Bursitis.
- Tendonitis of the Achilles tendon.

**Colour:**  Transparent.



**Ref.: TL-611**

**SILICONE HEEL CUP WITH CENTRAL SPUR**

SIZES	FOOTWEAR SIZE
1	35-38
2	39-42
3	43-46



**Units for packaging:** 1 pair.

**Characteristics:**

Due to their viscoelastic properties, these silicone heel supports cushion the load points both during walking and in static position. They have a circular area of softer consistency in the heel (blue area in).

**Indications:**

- Treatment of centralised calcaneus spur.
- Bursitis of the heel.
- Absorption of impacts.

**Colour:**  Transparent.



**Ref.: TL-612**

**SILICONE HEEL CUP WITH LATERAL SPUR**

SIZES	FOOTWEAR SIZE
1	35-38
2	39-42
3	43-46



**Units for packaging:** 1 pair.

**Characteristics:**

Due to their viscoelastic properties, these silicone heel supports cushion the load points both during walking and in static position. They have an oval area of softer consistency in the heel (blue area).

**Indications:**

- Lateralised spur.
- Bursitis of the heel.

**Colour:**  Transparent.



**Ref.: TL-613**

**ANATOMICAL SILICONE HEEL CUP FOR CENTRAL SPUR**

SIZES	FOOTWEAR SIZE
1	35-38
2	39-42
3	43-46



**Units for packaging:** 1 pair.

**Characteristics:**

Heel support made of silicone which, due to its viscoelastic properties, cushions the load points both during walking and in static position. It has an asymmetric area of softer consistency in the heel (blue zone). Base with anatomic form and greater width in the front part favouring the distribution of the loads between retroflex and half-flex, the control of the calcaneus and improves stability.

**Indications:**

- Treatment of calcaneus spur and bursitis of the heel.
- Tendonitis of the Achilles tendon.
- Plantar fasciitis and talalgia.

**Colour:**  Transparent.



Ref.: TL-617-10 (10 mm)

Ref.: TL-617-05 (5 mm)

**FLAT SILICONE HEEL CUP WITH CENTRAL SPUR**

SIZES	FOOTWEAR SIZE
1	35-38
2	39-42
3	43-46
4	47-50



Units for packaging: 1 pair.

**Characteristics:**

Due to their viscoelastic properties, these silicone heel supports cushion the load points both during walking and in static position. They have a circular area of softer consistency in the heel (blue area). Without lateral edges and with a height of 10 mm (TL-617-10) and 5 mm (TL-617-05).

**Indications:**

- Treatment of centralised calcaneus spur.
- Bursitis of the heel and absorption of impacts.
- In lower member dysmetrias of 1 cm (TL-617-10) and 0.5 cm (TL-617-05).

Colour:  Transparent.



Ref.: TL-619

**SILICONE HEEL CUP FOR EXCESSIVE PRONATION OR SUPINATION**

SIZES	FOOTWEAR SIZE
1	35-38
2	39-42
3	43-46



Units for packaging: 1 pair.

**Characteristics:**

Silicone heel cup for excessive pronation or supination. It corrects heel deviation to improve posture and prevent disorders caused by poor positioning when walking.

**Indications:**

- For patients with excessive pronation or supination of the heels.

Colour:  Transparent.



Ref.: TL-620

**SILICONE PRONOSUPINATION HEEL CUP**

SIZES	FOOTWEAR SIZE
1	27-33
2	33-39
3	39-45



Units for packaging: 1 pair.

**Characteristics:**

Silicone heel cup for aligning the heel in pronation or supination. It corrects heel deviation to improve posture and prevent disorders caused by poor positioning when walking. Marked slant (18:7 mm).

**Indications:**

- Patients with excessive pronation or supination of the heel.

Colour:  Transparent.




Ref.: TL-618F

**LINED SILICONE HEEL CUP**

SIZES	FOOTWEAR SIZE
1	35-38
2	39-42
3	43-46



Units for packaging: 1 pair.



**Characteristics:**

Perfect housing for the heel due to its bowl shape and the hypoallergenic polyester fabric which improves comfort and breathability and impedes bacterial growth. Easy to wash. New attachment system that uses Velcro to stick to the footwear so it does not move. Ideal height which allows it to be adapted to any type of footwear. It is a long term product due to its greater silicone resistance. Made from 2 density silicone.

**Indications:**

- Provides swift relief for ligaments and muscles especially in excess load conditions.
- It is also suitable for bursitis and tendonitis of the Achilles Heel (impact absorption).

Colour:  Transparent and grey lining.



Ref.: PL-755

**LONG SILICONE INSOLE WITHOUT RETROCAPITAL SUPPORT**

SIZES	FOOTWEAR SIZE
0	35-36
1	37-38
2	39-40
3	41-42
4	43-44
5	45-46



Units for packaging: 1 pair.

**Characteristics:**

Insoles made with viscoelastic silicone with metatarsal and calcaneus supports in different density and with bridge.

**Indications:**

- Talalgia.
- Bursitis.
- Metatarsalgia.
- Calcaneus spur.
- Plantar fasciitis.
- Diabetic foot.
- Painful, tired or post-surgical feet. In sporting practice.

Colour:  Transparent.



**Ref.: PL-755F**

**LINED LONG SILICONE INSOLE WITHOUT RETROCAPITAL SUPPORT**

SIZES	FOOTWEAR SIZE
0	35-36
1	37-38
2	39-40
3	41-42
4	43-44
5	45-46



Units for packaging: 1 pair.

**Characteristics:**

Insoles made with viscoelastic silicone with metatarsal and calcaneus supports in different density and bridge. Lined with hypoallergenic fabric in polyester that prevents bacterial growth and is easily washable.

**Indications:**

- Talalgia and bursitis.
- Metatarsalgia.
- Calcaneus spur.
- Plantar fasciitis.
- Diabetic foot.
- Painful, tired or post-surgical feet. In sporting practice.

**Colour:**  Transparent and grey lining.



**Ref.: PL-750**

**SILICONE INSOLE**

SIZES	FOOTWEAR SIZE
0	35-36
1	37-38
2	39-40
3	41-42
4	43-44
5	45-46



Units for packaging: 1 pair.

**Characteristics:**

Insoles made with viscoelastic silicone with metatarsal and calcaneus supports in different density and with retrocapital support and bridge.

**Indications:**

- Talalgia and bursitis.
- Metatarsalgia and calcaneus spur.
- Plantar fasciitis and diabetic foot.
- Painful, tired or post-surgical feet. In sports.

**Colour:**  Transparent.



**Ref.: PL-750F**

**LINED SILICONE INSOLE**

SIZES	FOOTWEAR SIZE
0	35-36
1	37-38
2	39-40
3	41-42
4	43-44
5	45-46



Units for packaging: 1 pair.

**Characteristics:**

Insoles made with viscoelastic silicone with metatarsal and calcaneus supports in different density and with retrocapital support and bridge. Lined with hypoallergenic fabric in polyester that prevents bacterial growth and is easily washable.

**Indications:**

- Talalgia and bursitis.
- Metatarsalgia and calcaneus spur.
- Plantar fasciitis and diabetic foot.
- Painful, tired or post-surgical feet. In sporting practice.

**Colour:**  Transparent and grey lining.



**Ref.: PL-760**

**3/4 SILICONE INSOLE**

SIZES	FOOTWEAR SIZE
1	37-38
2	39-40
3	41-42
4	43-44



Units for packaging: 1 pair.

**Characteristics:**

Insoles made with viscoelastic silicone with metatarsal and calcaneus supports in different density and with bridge. The point of the foot is free.

**Indications:**

- Talalgia and bursitis.
- Metatarsalgia and calcaneus spur.
- Plantar fasciitis and diabetic foot.
- Painful, tired or post-surgical feet. In sporting practice.

**Colour:**  Transparent.



**Ref.: PL-760F**

**3/4 LINED SILICONE INSOLE**

SIZES	FOOTWEAR SIZE
1	37-38
2	39-40
3	41-42
4	43-44



Units for packaging: 1 pair.

**Characteristics:**

Insoles made with viscoelastic silicone with metatarsal and calcaneus supports in different density and with bridge. Lined with hypoallergenic fabric in polyester that prevents bacterial growth and is easily washable.

**Indications:**

- Talalgia and bursitis.
- Metatarsalgia and calcaneus spur.
- Plantar fasciitis and diabetic foot.
- Painful, tired or post-surgical feet. In sporting practice.

**Colour:**  Transparent and grey lining.



Ref.: PL-702

**INSOLE WITH METATARSAL ARCH**

SIZES	FOOTWEAR SIZE
1	34-36
2	37-39
3	40-42



Units for packaging: 1 pair.

**Characteristics:**

An insole made from superfine flexible viscoelastic silicone, which can adapt to most footwear due to its reduced thickness, and featuring a metatarsal arch that helps to align the metatarsal heads, providing pain relief to the forefoot and preventing the formation of calluses in people who spend considerable time standing or walking. The insole can be cut to fit the footwear.

**Indications:**


- Metatarsalgia and forefoot pain.

**Colour:**  Transparent.

Ref.: PL-701F

**LINED EXTRA-FINE INSOLE**

SIZES	FOOTWEAR SIZE
0	35-36
1	37-38
2	39-40
3	41-42
4	43-44
5	45-46



Units for packaging: 1 pair.

**Characteristics:**

Insoles made with flexible extra-fine viscoelastic silicone, which allows great adaptation to the majority of footwear due to thinness, with metatarsal and calcaneus supports in different density and with retrocapital support and bridge. Lined with hypoallergenic fabric in polyester that prevents bacterial growth and is easily washable.

**Indications:**

- Talalgia and bursitis.
- Metatarsalgia and calcaneus spur.
- Plantar fasciitis and diabetic foot.
- Painful, tired or post-surgical feet. In sporting practice.

**Colour:**  Transparent and grey lining.



**Exhibitor**

Ref.: Exp-FINE-Cartón

Units: 10  
Sizes: 34 x 25 cm




Ref.: PL-700F

**EXTRA-THIN LINED SILICONE INSOLES WITH METATARSAL PAD FOR LADIES**

SIZES	FOOTWEAR SIZE
1	34-36
2	37-39
3	40-42

Units for packaging: 1 pair.



**Characteristics:**

Flexible extra-thin silicone insoles, with metatarsal pad which adapt perfectly to ladies shoes (high-heels or flats), due to their minimum thickness, ideal for all-day use. Lined with hypoallergenic polyester fabric that prevents bacterial growth, easy to clean. At the heel they incorporate an anti-slip system, preventing them from sliding out of place.

**Indications:**

- At the front of the foot they include a pad that provides a continuous damping effect in the metatarsal area, making them extremely comfortable.
- They relieve pain from the sole of the foot, preventing the formation of hard patches, they are shock-absorbing and reduce callosities.
- Pre-cut insoles.

**Colour:**  Transparent and grey lining.

**Exhibitor**

Ref.: Exp-FINE-Cartón

Units: 14  
Sizes: 34 x 25 cm



Ref.: PS-17

**EXTRA-FINE SILICONE FOREFOOT CUSHION**

SIZES	ONE SIZE

Units for packaging: 1 pair.

**Characteristics:**


The forefoot cushions are made with flexible extra-fine silicone, which allows perfect adaptation to the shoe, throughout the medial area of the foot. Its design incorporates raised circles on the external side.

**Indications:**

- Provides pain relief of the metatarsal and forefoot area, preventing formation of calluses and hardness.

**Colour:**  Transparent.





**5.3%**  
LESS VIBRATION IN THE LUMBAR REGION

**17.2%**  
LESS VIBRATION IN THE KNEE

**12.6%**  
LESS VIBRATION IN THE ANKLE

**44%**  
GREATER FOOT TREAD STABILITY

FOOT ORLIMAN FEETPAD



**Ref.: PL-BION**

**CUSTOM BIOMECHANICAL INSOLE TO PREVENT JOINT PAIN**

SIZES	36	37	38	39	40	41	42	43	44	45	46
-------	----	----	----	----	----	----	----	----	----	----	----


**Characteristics:**

A custom insole that can adapt to any type of foot. It can adapt to the anatomy of the foot to provide greater comfort when walking. 100% Total Dry. Composed of On Steam fabric - a microfibre that provides maximum breathability. It keeps the feet at the necessary temperature. Antibacterial treatment through the use of carbon.

**Indications:**

- Mild disorders in the biomechanics of the foot.
- Reduces vibration and joint pain.

**Colour:**  Grey.



Before

After



**Exhibitor**

**Ref.: EXP-BION**  
Exhibitor Bion-Tech  
32x45 cm / 8 units



FITTING VIDEO



TECHNOLOGY AND DESIGN FOR THE WELLBEING OF THE FEET



Ref.: PL-710

DAILY ACTIVITY GEL INSOLES

SIZES	FOOTWEAR SIZE
1	36-39
2	40-43
3	44-47

Units for packaging: 1 pair.

Characteristics:

- ① Heel support: gel with anatomical design adapted to the heel for greater cushioning and improved momentum when walking.
- ② Plantar arch: thermoplastic elastomer structure with greater rigidity and an anatomical design.
- ③ Metatarsal support: gel with anatomical design adapted to the metatarsus or greater cushioning.
- ④ Contact insole: soft and breathable fabric with anti-bacterial treatment containing: Dimethyloctadecyl[3-(trimethoxysilyl)propyl]ammonium chloride, to prevent bad smells.

Indications:

Help to absorb impacts effectively and reducing pressure on the feet caused by everyday footwear and feeling comfortable and rested throughout the day. In addition to helping to provide greater stability and support to the plantar arch and cushioning to the metatarsus and heel. The sizing and adjustment system enables each size to be anatomically adapted to the appropriate shoe size. They can be used with everyday footwear.

Colour: ■ Blue.

Ref.: PL-720

INTENSE ACTIVITY GEL INSOLES

SIZES	FOOTWEAR SIZE
1	36-39
2	40-43
3	44-47

Units for packaging: 1 pair.

Characteristics:

- ① Heel support: thermoplastic elastomer gel with anatomical design for greater cushioning.
- ② Plantar arch: special thermoplastic PU structure for prolonged periods of standing with an anatomical design for increased stability and support in the arch and heel
- ③ Metatarsal support: perforated SBR pad with an anatomical design for a high degree of energy absorption, impact protection and improved momentum when walking.
- ④ Contact insole: High strength and breathable fabric with anti-bacterial treatment containing: Dimethyloctadecyl[3-(trimethoxysilyl) propyl] ammonium chloride, to prevent bad smells.

Indications:

Help to absorb impacts effectively and reduce the excessive pressure that the feet have to withstand during prolonged periods standing. They help provide reinforced stability and support to the plantar arch and cushioning to the metatarsus and heel. The sizing and adjustment system enables each size to be anatomically adapted to the appropriate shoe size. They can be used with everyday footwear.

Colour: ■ Grey.

**MAIN FEATURES**



**Polymer gel:** This is the name of a special mixture of polymers which are not toxic, to which a medicinal mineral oil has been incorporated.

► **POLYMER GEL TECHNOLOGY**

Thanks to the medicinal mineral oil which it contains, it hydrates, softens, protects and absorbs the affected area. Orliman Gels are transparent, not toxic, hypoallergenic, dermatologically tested and do not favor the growth of bacteria.

► **NUTRITIOUS ENVIRONMENT**

The specific gels based on the latest advances of Dermal Application Systems Technology provide a nutritious environment which favours healing and healthy skin.

► **WASHABLE AND REUSABLE**

All Orliman products are washable and reusable yet however they retain all their advantages. In fact, adequate cleaning and care contribute to a long life of the product.

► **CUSHIONS AND PROTECTORS**

Given the viscoelastic properties similar to those of fluids, our gels "flow" providing relief in situations of impact, sudden force, pressure, shock or vibration.



**Ref.: GL-100**

**TOE SPREADERS**

SIZES	UNITS	
S	2	
M		
L		

► **Characteristics:**

Toe spreader made with non-toxic, hypoallergenic, dermatologically-tested, viscoelastic polymer gel that does not favour bacterial growth

► **Indications:**

- Helps to preserve the alignment of the metatarsal phalangeal joint (ideal for the deformity of the bunion and deviation of the toe in post-operative care).
- Reduces friction, abrasion and irritation between adjacent or superimposed toes.

► **Colour:**  Transparent.



**Ref.: GL-123**

**TOE SPREADER WITH TOE LOOP**

SIZES	UNITS	
M	2	
L		

► **Characteristics:**

The toe separator with ring for the big toe separates the toes to help maintain correct alignment, prevent chafing and relieve pressure on the big toe joint and painful bunions. Anatomically moulded for comfortable insertion between the toes. The flexible ring on the big toe holds the toe separator in place. The gel slowly releases medical-grade mineral oil and vitamins to protect the skin.

► **Indications:**

- Aligns deviating toes.
- Relieves pressure on the metacarpophalangeal joint.
- Separates and protects the toes.
- Washable and reusable.
- Patented design.

► **Colour:**  Transparent.



**Ref.: GL-101**

**TOE SEPARATORS**

SIZES	UNITS	
S	2	
M		
L		
XL		

► **Characteristics:**

Half-moon separator made with nontoxic, hypoallergenic, dermatologically-tested, viscoelastic polymer gel that does not favour bacterial growth.

► **Indications:**

- Helps to protect against interdigit friction, abrasion and irritation caused by painful formations, such as spots, verrucas or ungual problems.
- Maintains alignment during healing for prolonged and post-operative use.

► **Colour:**  Transparent.

**NEW DISPLAY PODOLOGY**

**Exhibitor**

Ref.: EXP-SG

Units: 60

Sizes: 50x32 cm







Ref.: GL-102

**PURE GEL BUNION SHIELD**



SIZES	UNITS	POLYMER GEL
S	2	
L		

**Characteristics:**

Protector shaped around the bunion made with non-toxic, hypoallergenic, dermatologically tested, viscoelastic polymer gel that does not favour bacterial growth.

**Indications:**

- The bunion protector cushions and protects the bony protuberance.
- Absorbs pressure and friction.

**Colour:**  Transparent.



Ref.: GL-124

**BUNION TOE SPREADER**



SIZES	UNITS	POLYMER GEL
ONE SIZE	1	

**Characteristics:**

The bunion toe spreader relieves pain caused by bunions and properly aligns the big toe. It absorbs stress placed on the metacarpophalangeal joint. Provides protection between big toe and shoe. The product slides easily over the big toe.

**Indications:**

- Reduces pressure and friction.
- Relieves bunion pain.
- Helps big toe realignment.
- Fits both feet.
- Patented design.

**Colour:**  Transparent.



Ref.: GL-122

**PURE GEL SIMPLE BUNION SHIELD**



SIZES	UNITS	POLYMER GEL
ONE SIZE	2	

**Colour:**  Transparent.



Ref.: GL-121

**BUNION RELIEF SLEEVE**



SIZES	UNITS	GEL INSIDE
S	1	
L		

**Characteristics:**

Absorbs shock and shear forces. The 3mm thin gel pad fits easily into shoes and other footwear. Soft, mineral oil gel softens and comforts painful and sensitive bunions (hallux valgus). The 3mm gel pad is in direct contact with the skin, moisturising and softening the skin. The comfort fabric secures comfortably over the big toe and around the foot, preventing sliding.

**Colour:**  Beige.



Ref.: GL-125

**DOUBLE LOOP BUNION SPREADER**



SIZES	UNITS	POLYMER GEL
ONE SIZE	1	

**Characteristics:**

Double loop tailor bunion spreader eases bunion pain and provides proper toe alignment. With two soft gel loops, the spreader stays in perfect position while reducing pressure on the MP Joint. Product easily slips over the 4th and 5th toes.

**Indications:**



- Tailor's bunion.
- Feet callus.

**Colour:**  Transparent.



Ref.: GL-103

**PROTECTIVE BUNION SHIELD IN GEL WITH FABRIC**

SIZES	UNITS	 
S	1	
L		

**Characteristics:**

Protector that covers the area of the bunion, being placed on the toe. Made in viscoelastic polymer gel covered in elastic surgical fabric and with straps secured to the forefoot.

**Indications:**


- Cushions and protects the bony protuberance.
- Absorbs pressure and friction.

**Colour:**  Beige.



Ref.: GL-104

**DIGITAL PADS**

SIZES	UNITS	
S	2	
M		
L		

**Characteristics:**

Ring-shaped cut tubes and with non-toxic, hypoallergenic, dermatologically tested viscoelastic polymer gel, which does not favour bacterial growth.

**Indications:**



- Relieves the pressure and friction of the affected toe.
- Contributes to reducing scar tissue.
- For ampoules, hammer toes, corns and toe traumas.

**Colour:**  Beige.



Ref.: GL-105

**DIGITAL CAPS LINED**

SIZES	UNITS	 
S	1	
M		
L		

**Characteristics:**

Protector that covers the area of the bunion, being placed on the toe. Made in viscoelastic polymer gel covered in elastic surgical fabric and with straps secured to the forefoot.

**Indications:**



- Cushions and protects the bony protuberance.
- Absorbs pressure and friction.

**Colour:**  Beige.



Ref.: GL-117

**PURE GEL DIGITAL CAP**

SIZES	UNITS	 
S	1	
M		
L		

**Characteristics:**

Tubular, toe-shaped, with non-toxic, hypoallergenic, dermatologically tested, viscoelastic polymer gel, which does not favour bacterial growth.

**Indications:**

- Helps relieve pressure and friction on the toes.
- Hammer toes, superimposed or infraposed toes.
- Keratotic injury or partially amputated toes.
- Skin or ungual problems due to dryness.

**Colour:**  Transparent.



Ref.: GL-106

**LINED GEL TUBING**

SIZES	UNITS	LENGTH	  
XS	1	15	
S			
M			
L			

**Characteristics:**

15 cm tubes and with non-toxic, hypoallergenic, dermatologically tested, viscoelastic polymer gel that does not favour bacterial growth. They can be cut according to length.

**Indications:**

- Helps relieve pressure and friction.
- Protection and hydration of toes.
- It helps reduce scar tissue.

**Colour:**  Beige.








Ref.: GL-107

**CORN PADS**

SIZES	UNITS	LENGTH
S	1	10
M		
L		

**Characteristics:**

10 cm tubes and with one side of non-toxic, hypoallergenic, dermatologically tested, viscoelastic polymer gel, which does not favour bacterial growth. They can be cut according to length.

**Indications:**

- Helps relieve pressure and friction.
- Protection and hydration of toes.
- It helps reduce scar tissue.



**Colour:**  Beige.



Ref.: GL-116

**PURE GEL TUBING**

SIZES	UNITS	LENGTH
S	1	15
M		
L		

**Characteristics:**

15cm tube of non-toxic, hypoallergenic, dermatologically-tested, viscoelastic polymer gel, which does not favour bacterial growth.

**Indications:**

- Helps relieve pressure and friction.
- Protection and hydration of toes.
- Helps to reduce scar tissue.




**Colour:**  Transparent.



Ref.: GL-118

**GEL FULLY COATED DIGITAL MESH TUBE**

SIZES	UNITS	LENGTH
S	1	15
L		

**Characteristics:**

The fabric is fully lined with vitamin enriched gel. Helps relieve pain from corns, calluses, nail problems, dry cuticles and abrasions. The flexible mesh fabric stretches easily for ease of application. The gel will not flatten out, unlike foam products. The gel tubes can easily be cut to size. Washable and reusable.

**Colour:**  Beige.




Ref.: GL-115D RIGHT

Ref.: GL-115I LEFT

**PURE GEL TOE CREST**

SIZES	UNITS
S	1
M	
L	



FOOT	
R	L

**Characteristics:**

Small pads of non-toxic, viscoelastic polymer gel with ring secured to 2nd toe and with protuberance that raises the head of the first metatarsal.

**Indications:**

- They provide relief for people with claw or hammer toes.
- They offer maximum comfort and reduce pressure on the toes and metatarsal heads.


**Colour:**  Transparent.



Ref.: GL-200

**METATARSAL CUSHION LINED**

SIZES	UNITS
S	1 Pair
L	



**Characteristics:**

Pad in non-toxic, viscoelastic polymer gel, with ring secured to 2nd toe, covered in elastic surgical fabric.

**Indications:**

- Protects the metatarsal area of the foot against friction and abrasion.
- Indicated in metatarsalgia.
- Helps to relieve the burning sensation and pain under the metatarsal heads.
- Relieves the rubbing of corns, small neuromas, and tarsal adipose cushions.
- It can be used as complement of soft cushion for metatarsalgia.

**Colour:**  Beige.





**Ref.: GL-201**

**PURE GEL METATARSAL CUSHION**

SIZES	UNITS
S	1 Pair
L	



**Characteristics:**

Pad with ring secured to 2nd toe, in non-toxic, viscoelastic polymer gel.

**Indications:**

- Protects the metatarsal area of the foot against friction and abrasion
- Indicated in metatarsalgia which require an important unload.
- Helps to relieve the burning sensation and pain under the metatarsal heads.
- Relieves the rubbing of corns, small neuromas, and tarsal adipose cushions.
- It can be used as complement of soft cushion for metatarsalgia.

**Colour:**  Transparent.



**Ref.: GL-206**

**METATARSAL ELASTIC BAND WITH ELEVATION**

SIZES	SHOE SIZE	PERIMETER CM	UNITS
S	35-37	16-19	2
M	38-40	19-22	
L	41-43	22-26	

**Characteristics:**

Manufactured in stretch fabric, it provides support for the 2nd, 3rd and 4th metatarsals in the plantar area.

**Indications:**

Metatarsalgia, collapsed anterior arch of the foot (metatarsus latus)

**Colour:**  Beige.



**Ref.: GL-207**

**ELASTIC BAND**

SIZES	SHOE SIZE	PERIMETER CM	UNITS
S	35-37	16-19	2
M	38-40	19-22	
L	41-43	22-26	

**Characteristics:**

Manufactured in stretch fabric.

**Indications:**

Metatarsalgia, collapsed anterior arch of the foot (metatarsus latus)

**Colour:**  Beige.



**Ref.: GL-208**

**DRESS BUNION RELIEF SLEEVE WITH GEL METATARSAL PAD**

SIZES	UNITS
S	1 Pair
L	



**Characteristics:**

The sleeve is made from a fine stretch fabric and contains a polymer gel pad that cushions the metatarsal and bunion area. Discreet and non-slip. Can fit either foot.

**Indications:**

- Reduces pressure and friction.
- Eases bunion pain.
- Absorbs shock.
- Prevents the formation of corns and reduces calluses.
- Relieves pain in the metatarsal area.

**Colour:**  Beige.



**Ref.: GL-202D RIGHT**

**Ref.: GL-202I LEFT**

**METATARSAL ELASTIC BAND WITH GEL CUSHION**

SIZES	UNITS
S	1
L	



**Characteristics:**

Tubular elastic band to be placed on the forefoot. With separation for the first toe and with plantar cushion in non-toxic, viscoelastic polymer gel.

**Indications:**

- Designed to cushion the plantar area and the metatarsal head.
- Relieves plantar pain.
- Shock absorber. Reduces calluses.

**Colour:**  Beige.

**Ref.: FR-GL-202 D/I Metatarsal elastic band with cushion 2cm more longer than Ref: GL-202 D/I**



**Ref.: GL-203**

**GEL FOOT COVER**



SIZES	UNITS	
S	1 Pair	
M		
L		

**Characteristics:**

Forefoot protector in non-toxic, viscoelastic polymer gel coated with elastic surgical fabric.

**Indications:**

- Helps to retain the body temperature in the toes.
- Comforts and relieves symptoms in individuals suffering from arthritis or cold toes due to poor circulation.

**Colour:** ■ Beige.



**Ref.: GL-300**

**COMPLETE FOOT PROTECTOR IN GEL PINKY WITH FABRIC**



SIZES	UNITS	
S 36-38	1 Pair	
M 39-41		
L 42-44		

**Characteristics:**

Non-toxic, viscoelastic polymer gel sock covered with elastic surgical fabric.

**Indications:**

- Ideal for tired feet.
- Ideal for diabetic feet.
- Helps to retain body temperature in the feet.
- Comforts or alleviates symptoms in those suffering from arthritis or cold toes due to poor circulation.

**Colour:** ■ Beige.



**Ref.: GL-204** *Conforgel*

**MINI-PLANTAR BANDS WITH GEL**



SIZES	UNITS	
S <41	1 Pair	
L ≥41		

**Characteristics:**

Protect and relieve the pain in the metatarsal area, preventing the formation of hardness of the skin, shock absorber and reduces calluses.

**Indications:**

- The Mini plantars Conforgel bands, are manufactured with a fine elastic material and in the interior part they have a polymer gel cushion that provides a rapid relief to the metatarsal area.
- The minibands don't move.
- They are discreet and of sparse thickness, so you can use them in all kind of shoe.

**Colour:** ■ Beige.



**Ref.: GL-205**

**TOE STRAIGHTENER DOUBLE TOE**

SIZES	UNITS
ONE SIZE	1

**Characteristics:**

The toe splint is the ideal solution to relieve pain caused by axial deviation, fractures or hammer toe. The splint helps to align the toes by means of adjustable rings while the soft double-layered anti-slip pad provides comfort for the metatarsal area of the foot.

**Indications:**

- Separates and aligns deformed, fractured and/or hammer toes.

**Colour:** ■ Beige.



**Ref.: PS-18**

**INVISIBLE FOREFOOT GEL MINIPADS.(SOFY-PLANT MINI)**

SIZES	UNITS
ONE SIZE	1 Pair

**Characteristics:**

Maximum relaxation and total discretion; forget foot sole pain all day

**Indications:**

- SOFY-PLANT MINI Invisible Forefoot Gel Minipads are designed to afford comfort and relief in the anterior region of the sole of the foot.
- They prevent the formation of corns and hardened areas, and offer anti-slipping characteristics and comfort.

**Colour:**  Transparent.

**SELF ADHESIVE**



**Ref.: PS-19**

**INVISIBLE FOREFOOT RELIEF PAD/BALL OF FOOT GEL CUSHION**

SIZES	UNITS
ONE SIZE	1 Pair

**Characteristics:**

Relief of burning feet pain. Designed to absorb shock and pressure. Recommended for high heel shoes, sandals and boots. Provide all-day cushioning. Ultra slim to fit comfortably into shoes, sandals and boots. Durable. Transparent to match any shoe. Self adhesive. Non slip. Washable and re-usable. The gel pads will not compress like foams and felt products.

**Colour:**  Transparent.

**SELF ADHESIVE**



**Ref.: PS-20**

**SANDAL SPREADER WITH TEARDROP**

SIZES	UNITS
ONE SIZE	1 Pair

**Characteristics:**

Helps to reduce friction and irritation between the toes, ensuring maximum comfort. The pad provides comfort for the middle metatarsals. Absorbs shock and pressure. Slim design for maximum comfort. Self adhesive for easy fit. Stay in place. Transparent. Washable and reuseable. One size.

**Colour:**  Transparent.





**Exhibitor**

**Ref.: Exp-SF-Carton**  
Units: 10  
Sizes: 30 x 20 x 20cm



**SELF ADHESIVE**



**Ref.: PS-21**

**FLEXIBLE SELF ADHESIVE HEEL GRIPS**

SIZES	UNITS
ONE SIZE	1 Pair

**Characteristics:**

Washable, durable and hypoallergenic. The gel will not bottom out like foams. Ultra slim for easy fit. Self adhesive for easy fit. Transparent. Suitable for boots, shoes, and sandals. One size.

**Colour:**  Transparent.

**SELF ADHESIVE**



**Ref.: PS-22**

**3/4 LENGTH ADHESIVE GEL INSOLES WITH PRESSURE-RELIEF AREAS ON HEEL & FOREFOOT**

SIZES	UNITS
ONE SIZE	1 Pair

**Characteristics:**

3/4-length adhesive insoles with antislip invisible gel that absorb pressure and provide comfort for the feet when wear your favourite shoes (high heels or flats) all day long, whether you are at work or during your leisure time. The ultra-thin insoles adapt easily inside the shoes preventing the burning feeling on the feet. The strategically located pressure relief areas provide additional cushioning under the heel and the forefoot. They absorb impacts and pressure. The adhesive prevents the insole from slipping and moving. Transparent. Washable and reusable. One size.

**Colour:**  Transparent.

**SELF ADHESIVE**



**Ref.: PS-23**

**ADHESIVE GEL HEEL CUSHIONS**

SIZES	UNITS
ONE SIZE	1 Pair

**Characteristics:**

Ultra-thin anti-slip transparent gel pads that easily adapt to the back part of the inside of your shoes to prevent that burning feeling under the heels, the pressure-relief area absorbs impacts and pressure. The adhesive prevents them slipping and moving.

**Colour:**  Transparent.



Ref.: HV-30 RIGHT | Ref.: HV-31 LEFT

## NIGHTTIME HALLUX-VALGUS CORRECTOR

SIZES	ONE SIZE	
-------	----------	--

### Characteristics:

Made in internal cotton bouclé in contact with the Beige. Comprises an abductor splint of moldable aluminium to regulate the position of the big toe. It covers the entire foot, and the perimeter can be regulated.

### Indications:

- Correction of the Hallux-Valgus posture.
- Post-surgical.

**Colour:** Beige.



Ref.: HV-32

## DAYTIME HALLUX-VALGUS CORRECTOR

SIZES	FOOTWEAR SIZE
1	35-38
2	39-44



### Characteristics:

Comprises an elastic band around the instep and a cap for the big toe that is extended by an elastic band behind the Achilles tendon and fastens with Velcro for traction up to optimum position.

### Indications:

- Hallux- valgus correction.

**Colour:** Beige.



Ref.: HV-33D RIGHT | Ref.: HV-33I LEFT

## HALLUX-VALGUS IN THERMOPLASTIC

SIZES	FOOTWEAR SIZE
1	34-36
2	37-40
3	41-46

### Characteristics:

Night-time corrector for bunions, manufactured in thermoplastic and velour, with micro-fastener seal. It includes a cushioning pad on the big toe which, together with the pressure strap, at the height of the bunion, maintains the first metatarsus aligned and in correct position during nighttime use, reducing pain and preventing its possible evolution.

### Indications:

- Hallux- valgus correction.

**Colour:** Blue.



Ref.: HV-34

## conex<sup>®</sup> HALLUX VALGUS CORRECTOR

SIZES	ONE SIZE
Maximum half-foot perimeter: 30 cm.	

### Characteristics:

Orthosis for dynamic correction of the position of the big toe. Comprising an elastic band that wraps around the midfoot and an elastic strap with gel that acts on the big toe through controlled pressure that promotes gradual valgus alignment. Its elasticity and malleability enable a dynamic action during walking and provide a perfect fit to the foot. The intensity of the correction is fully adjustable to the specific needs of each patient.

Comfortable and easy to fit, it can be used with any type of footwear. (Patented design).

### Effects:

- Hallux valgus correction.
- Post-surgical maintenance of the correction.
- Preventing the deformity from increasing.

### Indications:

- Correction of big toe valgus deformity;
- Pain relief by reducing capsule and ligament tension.
- Prevents the deformity from increasing by impeding incorrect biomechanical traction of the tendons.

### Colour:

Beige.



