

Limb *Wear*

BODY COMFORT

by LimbTex



1st Edition
PROSTHETIC SOCKS

Ordering information

All of our customers are very important to us, so we try our best to despatch all orders received by 1.30pm the same day.

How to place an order:

When placing an order please provide the following information:

- An order number
- Product codes, descriptions and sizes
- Invoice and delivery addresses

We accept orders in the following ways:

- Telephone: +44 (0) 116 278 5440
- Facsimile: +44 (0) 116 278 0227
- Email: info@limbtex.com
- Post: LimbTex Ltd, Units 1 & 2, Elizabeth Business Park, Tigers Close, South Wigston, Leicestershire, LE18 4AE United Kingdom

Please call +44 (0) 116 2785440 for our current delivery charges

Materials manufactured or ordered to customer requirements:

We are able to offer materials to your specification, please call our customer services team for advice.

Delivery Satisfaction:

Please notify our customer service team of any incorrect or damaged orders within 7 days. After 7 days LimbTex Ltd reserve the right not to accept responsibility for any claims.

Return of Goods:

All materials may be returned within 30 days for a credit less 20% re-stocking charge. All materials received must be in clean, unused condition. Materials manufactured at your request are non-returnable.

Pricing:

Pricing is available separately on request and subject to VAT at the current rate.

Payment:

Our Payment terms are strictly 30 days from date of invoice, Cheques should be made payable to LimbTex Ltd.

Specification:

LimbTex Ltd reserves the right to alter specification or withdraw products from sale.

Ownership of Goods:

All goods supplied remain the property of LimbTex Ltd until paid for in full.

Limb *Wear*
BODY COMFORT
by **LimbTex**

www.limbwear.co.uk

LimbWear by LimbTex

Sheaths, Plain & Terry Upper & Lower Limb Prosthetic Socks

LimbWear Cotton Prosthetic socks are made from the highest quality yarns and all benefit from our "Non-restrictive cuff". Although non-restrictive, the cuff design and the use of Lycra® ensures that they retain their shape and fit perfectly.

Our socks are manufactured on state of the art machinery and do not have any inner seams that may cause discomfort.

As these are class 1 medical products all the materials used conform to the Öko-Tex Standard 100, so meet all necessary requirements. The cotton we use is very soft to give the softest feel possible to our finished products and comfort to the user.

As pin-locking systems for silicone liners are very common, we offer a reinforced distal hole across the whole cotton range.

All of our socks are pre-shrunk, so the size you order retains its size wash after wash.

Our sheer Nylon sheaths are ideal to use as donning aids or directly against the skin to add to the comfort of the user.

If you have particular requirements for other yarns, treatments, sizes etc, please contact us as we are able to custom manufacture including your logo's etc.

Contents

Prosthetic Sock Range: pages 4-5

Prosthetic Sheaths: page 6

Childrens Trans-Tibial & Upper Limb Socks: page 7

Trans-Tibial Socks: pages 8-11

Trans-Femoral Socks: pages 12-16

Chopart Socks: page 17

Sock Selection Chart: page 18

Sock Care Guide: page 19

Limb Wear

BODY COMFORT

Our Sheaths & Socks

Every style of sock/sheath is available in a variety of lengths to suit all users and we can manufacture socks to your specifications if required (min orders will apply). Each sock and sheath is packaged individually and can be ordered individually.



NYLON SHEATHS

Thinnest sock/sheath we offer

Nylon Sheaths are manufactured from fine nylon yarn that offers comfort if worn against the skin or over the liner. The open mesh design allows airflow keeping the residual limb comfortable. These sheaths can also be used for donning prosthetic limbs.



PLAIN KNIT PROSTHETIC SOCKS

Middle thickness sock we offer

Plain knit provides a comfortable, breathable, cushioning interface between the residual limb and the prosthetic limb socket.



TERRY TOWELLING KNIT PROSTHETIC SOCKS

Thickest sock we offer

Terry Towelling knit provides a very comfortable, breathable cushioning interface between the residual limb and the prosthetic limb socket.





Our plain knit and terry socks are available in FIVE widths...

SOCKS WITH A LOCKING PIN HOLE

Available in Plain and Terry knit

A variety of our socks are designed to be used with a Locking Pin Liner so have a reinforced hole in the distal end of the sock for the pin to pass through. These socks are to help with residual limb volume fluctuation or to extend the life of a prosthetic socket if it becomes too loose. The hole is sewn in a specific way so that it reinforces the fibres and prevents the yarns from fraying which could cause problems by fouling the Pin Lock etc.



BI THICKNESS SOCK

Bi Thickness Socks are designed for the Trans-Tibial wearer and are constructed with either Plain knit at the top and Terry knit at the bottom or Terry knit at the top and Plain knit at the bottom. They are used to compensate for volume increase or reduction at the proximal or distal end of the residual limb. This is a much better option than the old fashioned way of using a spare sock cut in half or the more time consuming half leather lining.




LimbWear Nylon Sheaths

Manufactured in 100% Nylon with non-restrictive cuff

TTS:


Trans-Tibial size




Product Code	Trans-Tibial Sheaths	Min Circ of Cuff in cm	Max Circ of Cuff in cm	Length Code Cuff Bottom Colour
TTS25	Trans-Tibial Sheath 25cm Long	30	47	Yellow
TTS30	Trans-Tibial Sheath 30cm Long	30	47	Blue
TTS35	Trans-Tibial Sheath 35cm Long	30	47	Red
TTS45	Trans-Tibial Sheath 45cm Long	30	47	Black
TTS55	Trans-Tibial Sheath 55cm Long	30	47	Green
TTS65	Trans-Tibial Sheath 65cm Long	30	47	Pink

TFS:

Trans-Tibial/Trans-Femoral size



Product Code	Trans-Tibial Sheaths	Min Circ of Cuff in cm	Max Circ of Cuff in cm	Length Code Cuff Bottom Colour
TFS25	Trans-Femoral Sheath 25cm Long	38	52	Yellow
TFS35	Trans-Femoral Sheath 35cm Long	38	52	Red
TFS45	Trans-Femoral Sheath 45cm Long	38	52	Black
TFS55	Trans-Femoral Sheath 55cm Long	38	52	Green
TFS65	Trans-Femoral Sheath 65cm Long	38	52	Pink



LimbWear Childs Trans-Tibial & Arm Mitt Socks

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain and Terry knit

TRANS-TIBIAL
CHILDRENS



CTTP:

Plain Knit


Product Code	Childs Trans-Tibial & Arm Mitt Plain Knit	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
CTTP10	Childs Trans-Tibial & Arm Mitt Plain 10cm Long	11	45		
CTTP15	Childs Trans-Tibial & Arm Mitt Plain 15cm Long	11	45		
CTTP20	Childs Trans-Tibial & Arm Mitt Plain 20cm Long	11	45		
CTTP25	Childs Trans-Tibial & Arm Mitt Plain 25cm Long	11	45		
CTTP30	Childs Trans-Tibial & Arm Mitt Plain 30cm Long	11	45		
CTTP35	Childs Trans-Tibial & Arm Mitt Plain 35cm Long	11	45		



CTTT:

Terry Knit

Product Code	Childs Trans-Tibial & Arm Mitt Terry	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
CTTT10	Childs Trans-Tibial & Arm Mitt Terry 10cm Long	11	45		
CTTT15	Childs Trans-Tibial & Arm Mitt Terry 15cm Long	11	45		
CTTT20	Childs Trans-Tibial & Arm Mitt Terry 20cm Long	11	45		
CTTT25	Childs Trans-Tibial & Arm Mitt Terry 25cm Long	11	45		
CTTT30	Childs Trans-Tibial & Arm Mitt Terry 30cm Long	11	45		
CTTT35	Childs Trans-Tibial & Arm Mitt Terry 35cm Long	11	45		
CTTT40	Childs Trans-Tibial & Arm Mitt Terry 40cm Long	11	45		
CTTT45	Childs Trans-Tibial & Arm Mitt Terry 45cm Long	11	45		



LimbWear Trans-Tibial Socks

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain and Terry knit

TTP:


Plain Knit



Product Code	Trans-Tibial Plain Knit	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTP20	Trans-Tibial Plain 20cm Long	16	55	Teal	Teal
TTP25	Trans-Tibial Plain 25cm Long	16	55	Teal	Olive Green
TTP30	Trans-Tibial Plain 30cm Long	16	55	Teal	Yellow
TTP35	Trans-Tibial Plain 35cm Long	16	55	Teal	Yellow
TTP40	Trans-Tibial Plain 40cm Long	16	55	Teal	Dark Teal
TTP45	Trans-Tibial Plain 45cm Long	16	55	Teal	Light Blue
TTP50	Trans-Tibial Plain 50cm Long	16	55	Teal	Black
TTP60	Trans-Tibial Plain 55cm Long	16	55	Teal	Red

TTT:


Terry Knit



Product Code	Trans-Tibial Terry Knit	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTT20	Trans-Tibial Terry 20cm Long	16	55	Yellow	Teal
TTT25	Trans-Tibial Terry 25cm Long	16	55	Yellow	Olive Green
TTT30	Trans-Tibial Terry 30cm Long	16	55	Yellow	Yellow
TTT35	Trans-Tibial Terry 35cm Long	16	55	Yellow	Yellow
TTT40	Trans-Tibial Terry 40cm Long	16	55	Yellow	Dark Teal
TTT45	Trans-Tibial Terry 45cm Long	16	55	Yellow	Light Blue
TTT50	Trans-Tibial Terry 50cm Long	16	55	Yellow	Black
TTT60	Trans-Tibial Terry 60cm Long	16	55	Yellow	Red

LimbWear Trans-Tibial Socks with Locking Pin Hole


Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain and Terry knit



TTPH:

Plain Knit with Hole

Product Code	Trans-Tibial Plain Knit with Hole	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTPH20	Trans-Tibial Plain with Hole 20cm Long	16	55	Green	Teal
TTPH25	Trans-Tibial Plain with Hole 25cm Long	16	55	Green	Light Green
TTPH30	Trans-Tibial Plain with Hole 30cm Long	16	55	Green	Yellow
TTPH35	Trans-Tibial Plain with Hole 35cm Long	16	55	Green	Yellow-Green
TTPH40	Trans-Tibial Plain with Hole 40cm Long	16	55	Green	Dark Teal
TTPH45	Trans-Tibial Plain with Hole 45cm Long	16	55	Green	Light Blue
TTPH50	Trans-Tibial Plain with Hole 50cm Long	16	55	Green	Black
TTPH60	Trans-Tibial Plain with Hole 60cm Long	16	55	Green	Red




TTTH:

Terry Knit with Hole

Product Code	Trans-Tibial Terry Knit	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTTH20	Trans-Tibial Terry with Hole 20cm Long	16	55	Yellow	Teal
TTTH25	Trans-Tibial Terry with Hole 25cm Long	16	55	Yellow	Light Green
TTTH30	Trans-Tibial Terry with Hole 30cm Long	16	55	Yellow	Yellow
TTTH35	Trans-Tibial Terry with Hole 35cm Long	16	55	Yellow	Yellow-Green
TTTH40	Trans-Tibial Terry with Hole 40cm Long	16	55	Yellow	Dark Teal
TTTH45	Trans-Tibial Terry with Hole 45cm Long	16	55	Yellow	Light Blue
TTTH50	Trans-Tibial Terry with Hole 50cm Long	16	55	Yellow	Black
TTTH60	Trans-Tibial Terry with Hole 60cm Long	16	55	Yellow	Red



TRANS-TIBIAL




LimbWear Trans-Tibial Bi Knit

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain/Terry knit

TTTUPL:


Bi-knit Terry Upper, Plain Lower



Product Code	Trans-Tibial Dual Terry Knit Upper/ Plain Knit Lower	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTTUPL25	Trans-Tibial Bi Thickness Terry/Plain 25cm Long	16	55	Teal	Olive Green
TTTUPL30	Trans-Tibial Bi Thickness Terry/Plain 30cm Long	16	55	Teal	Yellow
TTTUPL35	Trans-Tibial Bi Thickness Terry/Plain 35cm Long	16	55	Teal	Yellow
TTTUPL40	Trans-Tibial Bi Thickness Terry/Plain 40cm Long	16	55	Teal	Teal
TTTUPL45	Trans-Tibial Bi Thickness Terry/Plain 45cm Long	16	55	Teal	Blue

TTPUTL:



Bi-knit Plain Upper, Terry Lower



Product Code	Trans-Tibial Dual Plain Knit Upper/ Terry Knit Lower	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTPUTL25	Trans-Tibial Bi Thickness Plain/Terry 25cm Long	16	55	Black	Olive Green
TTPUTL30	Trans-Tibial Bi Thickness Plain/Terry 30cm Long	16	55	Black	Yellow
TTPUTL35	Trans-Tibial Bi Thickness Plain/Terry 35cm Long	16	55	Black	Yellow
TTPUTL40	Trans-Tibial Bi Thickness Plain/Terry 40cm Long	16	55	Black	Teal
TTPUTL45	Trans-Tibial Bi Thickness Plain/Terry 45cm Long	16	55	Black	Blue

LimbWear Trans-Tibial Bi Knit with Locking Pin Hole

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain/Terry knit



TTUPLH:

Bi-knit Terry Upper, Plain Lower with Hole



Product Code	Trans-Tibial Dual Terry Knit Upper/ Plain Knit Lower	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTUPLH25	Trans-Tibial Bi Thickness Terry/Plain with Hole 25cm Long	16	55	Light Blue	Green
TTUPLH30	Trans-Tibial Bi Thickness Terry/Plain with Hole 30cm Long	16	55	Light Blue	Yellow
TTUPLH35	Trans-Tibial Bi Thickness Terry/Plain with Hole 35cm Long	16	55	Light Blue	Yellow
TTUPLH40	Trans-Tibial Bi Thickness Terry/Plain with Hole 40cm Long	16	55	Light Blue	Dark Teal
TTUPLH45	Trans-Tibial Bi Thickness Terry/Plain with Hole 45cm Long	16	55	Light Blue	Light Blue



TTPUTLH:

Bi-knit Plain Upper, Terry Lower with Hole



Product Code	Trans-Tibial Dual Plain Knit Upper/ Terry Knit Lower	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TTPUTLH25	Trans-Tibial Bi Thickness Plain/Terry with Hole 25cm Long	16	55	Red	Green
TTPUTLH30	Trans-Tibial Bi Thickness Plain/Terry with Hole 30cm Long	16	55	Red	Yellow
TTPUTLH35	Trans-Tibial Bi Thickness Plain/Terry with Hole 35cm Long	16	55	Red	Yellow
TTPUTLH40	Trans-Tibial Bi Thickness Plain/Terry with Hole 40cm Long	16	55	Red	Dark Teal
TTPUTLH45	Trans-Tibial Bi Thickness Plain/Terry with Hole 45cm Long	16	55	Red	Light Blue




LimbWear Trans-Femoral Standard Width Socks

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain/Terry knit

TFP:

Plain Knit



Product Code	Trans-Femoral Plain Knit	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFP20	Trans-Femoral Plain 20cm Long	24	71	Green	Teal
TFP30	Trans-Femoral Plain 30cm Long	24	71	Green	Yellow
TFP40	Trans-Femoral Plain 40cm Long	24	71	Green	Dark Teal
TFP50	Trans-Femoral Plain 50cm Long	24	71	Green	Black
TFP60	Trans-Femoral Plain 60cm Long	24	71	Green	Red

TFT:



Terry Knit



Product Code	Trans-Femoral Terry Knit	Min Circ of Cuff in cm	Max Circ of Cuff in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFT20	Trans-Femoral Terry 20cm Long	24	71		Teal
TFT30	Trans-Femoral Terry 30cm Long	24	71		Yellow
TFT40	Trans-Femoral Terry 40cm Long	24	71		Dark Teal
TFT50	Trans-Femoral Terry 50cm Long	24	71		Black
TFT60	Trans-Femoral Terry 60cm Long	24	71		Red

LimbWear Trans-Femoral Standard Width Socks with Locking Pin Hole

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain/Terry knit Terry knit




TFPH:

Plain Knit with Hole



Product Code	Trans-Femoral Plain Knit with Hole	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFPH20	Trans-Femoral Plain with Hole 20cm Long	24	71	Teal	Yellow
TFPH30	Trans-Femoral Plain with Hole 30cm Long	24	71	Teal	Yellow
TFPH40	Trans-Femoral Plain with Hole 40cm Long	24	71	Teal	Yellow
TFPH50	Trans-Femoral Plain with Hole 50cm Long	24	71	Teal	Black
TFPH60	Trans-Femoral Plain with Hole 60cm Long	24	71	Teal	Red




TFTH:

Terry Knit with Hole



Product Code	Trans-Femoral Terry Knit with Hole	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFTH20	Trans-Femoral Terry with Hole 20cm Long	24	71	Teal	Yellow
TFTH30	Trans-Femoral Terry with Hole 30cm Long	24	71	Teal	Yellow
TFTH40	Trans-Femoral Terry with Hole 40cm Long	24	71	Teal	Yellow
TFTH50	Trans-Femoral Terry with Hole 50cm Long	24	71	Teal	Black
TFTH60	Trans-Femoral Terry with Hole 60cm Long	24	71	Teal	Red




LimbWear Trans-Femoral X Wide Width Socks

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain/Terry knit

TFXWP:

Plain Knit



Product Code	Trans-Femoral X Wide Plain Knit	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFXWP20	Trans-Femoral X Wide Plain 20cm Long	28	89	Orange	Teal
TFXWP30	Trans-Femoral X Wide Plain 30cm Long	28	89	Orange	Yellow
TFXWP40	Trans-Femoral X Wide Plain 40cm Long	28	89	Orange	Dark Teal
TFXWP50	Trans-Femoral X Wide Plain 50cm Long	28	89	Orange	Black
TFXWP60	Trans-Femoral X Wide Plain 60cm Long	28	89	Orange	Red

TFXWT:


Terry Knit



Product Code	Trans-Femoral X Wide Terry	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFXWT20	Trans-Femoral X Wide Terry 20cm Long	28	89	Pink	Teal
TFXWT30	Trans-Femoral X Wide Terry 30cm Long	28	89	Pink	Yellow
TFXWT40	Trans-Femoral X Wide Terry 40cm Long	28	89	Pink	Dark Teal
TFXWT50	Trans-Femoral X Wide Terry 50cm Long	28	89	Pink	Orange
TFXWT60	Trans-Femoral X Wide Terry 60cm Long	28	89	Pink	Red

LimbWear Trans-Femoral X Wide Width Socks with Locking Pin Hole

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain/Terry knit



TFXWPH:

Plain Knit with Hole



Product Code	Trans-Femoral X Wide Plain Knit with Hole	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFXWPH20	Trans-Femoral X Wide Plain with Hole 20cm Long	28	89		
TFXWPH30	Trans-Femoral X Wide Plain with Hole 30cm Long	28	89		
TFXWPH40	Trans-Femoral X Wide Plain with Hole 40cm Long	28	89		
TFXWPH50	Trans-Femoral X Wide Plain with Hole 50cm Long	28	89		
TFXWPH60	Trans-Femoral X Wide Plain with Hole 60cm Long	28	89		



TFXWTH:

Terry Knit with Hole



Product Code	Trans-Femoral X Wide Terry Knit with Hole	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFXWTH20	Trans-Femoral X Wide Terry with Hole 20cm Long	28	89		
TFXWTH30	Trans-Femoral X Wide Terry with Hole 30cm Long	28	89		
TFXWTH40	Trans-Femoral X Wide Terry with Hole 40cm Long	28	89		
TFXWTH50	Trans-Femoral X Wide Terry with Hole 50cm Long	28	89		
TFXWTH60	Trans-Femoral X Wide Terry with Hole 60cm Long	28	89		



LimbWear Trans-Femoral XX Wide Width Socks with and without Locking Pin Hole

Manufactured in soft Cotton/Nylon with non-restrictive cuff in Plain/Terry knit

TFXXWT:


Terry Knit



Product Code	Trans-Femoral XX Wide Terry	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFXXWT30	Trans-Femoral XX Wide Terry 30cm Long	33	95	Red	Yellow
TFXXWT40	Trans-Femoral XX Wide Terry 40cm Long	33	95	Red	Teal
TFXXWT50	Trans-Femoral XX Wide Terry 50cm Long	33	95	Red	Black
TFXXWT60	Trans-Femoral XX Wide Terry 60cm Long	33	95	Red	Orange

TFXXWTH:

Terry Knit with Hole



Product Code	Trans-Femoral XX Wide Terry with Hole	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code Cuff Upper Colour
TFXXWTH30	Trans-Femoral XX Wide Terry with Hole 30cm long	33	95	Olive	Yellow
TFXXWTH40	Trans-Femoral XX Wide Terry with Hole 40cm long	33	95	Olive	Teal
TFXXWTH50	Trans-Femoral XX Wide Terry with Hole 50cm long	33	95	Olive	Orange
TFXXWTH60	Trans-Femoral XX Wide Terry with Hole 60cm long	33	95	Olive	Red

LimbWear Chopart Socks

Manufactured in soft Cotton/Lycra soft ribbed cuff, heel & short foot length

CHOPART:


Plain Knit



Product Code	Chopart Sock	Min Circ in cm	Max Circ in cm	Style Code Cuff Lower Colour	Length Code
					Cuff Upper Colour
CHOPART30	CHOPART SOCK 30cm Long	18	55		Yellow
CHOPART40	CHOPART SOCK 40cm Long	18	55		Teal
CHOPART50	CHOPART SOCK 50cm Long	18	55		Black

notes:

CHOPART SOCKS




LimbWear Sock Selection Chart


Product Code (Add Length)	Description	Min Circ in cm's	Max Circ in cm's	Sock Length in cm's Colour Coding for Upper Part of Cuff										Sock Type Colour Coding Lower Part of cuff			
				10	15	20	25	30	35	40	45	50	60				
CTTP	Childs Trans-Tibial Plain and Arm Mitt	11	44														
CTTT	Childs Trans-Tibial Terry and Arm Mitt	11	44														
TTP	Trans-Tibial Plain	16	54														
TTPH	Trans-Tibial Plain with Hole	16	56														
TTT	Trans-Tibial Terry	16	54														
TTTH	Trans-Tibial Terry with Hole	16	54														
TTUPL	Trans-Tibial Bi Thickness Terry/Plain	16	54														
TTUPLH	Trans-Tibial Bi Thickness Terry/Plain with Hole	16	54														
TTPUTL	Trans-Tibial Bi Thickness Plain/Terry	16	54														
TTPUTLH	Trans-Tibial Bi Thickness Plain/Terry with Hole	16	54														
TFP	Trans-Femoral Plain	24	70														
TFPH	Trans-Femoral Plain with Hole	24	70														
TFT	Trans-Femoral Terry	24	66														
TFTH	Trans-Femoral Terry with Hole	24	64														
TFXWP	Trans-Femoral X Wide Plain	28	88														
TFXWPH	Trans-Femoral X Wide Plain with Hole	28	88														
TFXWT	Trans-Femoral X Wide Terry	28	80														
TFXWTH	Trans-Femoral X Wide Terry with Hole	28	80														
TFXXWT	Trans-Femoral XX Wide Terry	32	94														
TFXXWTH	Trans-Femoral XX Wide Terry with Hole	32	94														

Product Code (Add Length)	Description	Min Circ in cm's	Max Circ in cm's	Sock Length in cm's for Chopart Socks ONLY Colour Coding for Lower Part of Cuff										Sock Type Colour Coding Upper Part of cuff			
				10	15	20	25	30	35	40	45	50	60				
CHOPART	Chopart Sock	18	54														

Product Code	Trans-Tibial Sheaths	Min Circ in cm	Max Circ in cm
TTS25	Trans-Tibial Sheath 25cm Long	30	47
TTS30	Trans-Tibial Sheath 30cm Long	30	47
TTS35	Trans-Tibial Sheath 35cm Long	30	47
TTS45	Trans-Tibial Sheath 45cm Long	30	47
TTS55	Trans-Tibial Sheath 55cm Long	30	47
TTS65	Trans-Tibial Sheath 65cm Long	30	47



Product Code	Trans-Femoral Sheaths	Min Circ in cm	Max Circ in cm
TFS25	Trans-Femoral Sheath 25cm Long	38	52
TFS35	Trans-Femoral Sheath 35cm Long	38	52
TFS45	Trans-Femoral Sheath 45cm Long	38	52
TFS55	Trans-Femoral Sheath 55cm Long	38	52
TFS65	Trans-Femoral Sheath 65cm Long	38	52



Caring for LimbWear prosthetic Socks

For our prosthetic socks to remain comfortable and perform as they should, they will need to be cared for. Please follow these guidelines to get the most use out of them:

- Use a clean sock every day
- Avoid creases
- Do not mend your socks, if socks become uncomfortably thin, hard or have holes in you should discard them and use new ones.
- Our socks are machine washable but using fabric conditioners may cause skin irritation so try to avoid this.

Wash the socks inside out and tumble dry as this will help return them back to original condition. Washing instructions can be seen on each pack.



Limbtex Ltd
Units 1 & 2, Elizabeth Business Park
Tigers Close
South Wigston
Leicestershire
LE18 4AE
United Kingdom

Tel: **+44 (0) 116 278 5440**
Fax: **+44 (0) 116 278 0227**
Email: **info@limbtex.com**
Web: **www.limbwear.co.uk**

