



# Cosmetic Cover

Everyday Protection

Blatchford 

# Cosmetic Cover

The Cosmetic Cover provides optimal protection against wear and tear from everyday activities for Orion3 and Linx prosthesis.



## Protection

Whether it's scratches, impact or wear and tear, the cosmetic cover protects the prosthesis whilst being lightweight in design.

## Natural Look and Comfort

The Cosmetic Cover makes kneeling more comfortable by softening contact with the ground. The design also enhances a more natural shape under clothes.

## Simple to Use

The cover simply clips together so is easy to put on and take off. The design also allows for easy access to charging points and controls.

## Technical Specification

- Component Weight - 425g
- Component Length - 313mm
- Compatible with - Orion3, Linx
- Warranty - 6 months
- Part Number: 561101BLACK

