

Quality and function


NEW

get more information on
www.uniprox.de



The Ultra Stable Breathable Liner: SoftskinAir® Sport

- Perfect sweat transportation through micropores
- Double-layer silicone with a silky inner surface
- Additional distal matrix for more stability
- Shore hardness 40 Shore A
- Wall thickness 3 mm, length 370 mm
- Stainless steel pin holder
- Liner stand included
- Available in sizes 20, 22, 23.5, 25, 26.5, 28, 30



SoftskinAir® Sport

Order No.	Size	Article No.
Sport-20-Air	20	4 164 030 40 00 200
Sport-22-Air	22	4 164 030 40 00 220
Sport-23,5-Air	23,5	4 164 030 40 00 235
Sport-25-Air	25	4 164 030 40 00 250
Sport-26,5-Air	26,5	4 164 030 40 00 265
Sport-28-Air	28	4 164 030 40 00 280
Sport-30-Air	30	4 164 030 40 00 300



unique prosthetic solutions

A company of the Bauerfeind Group

uniprox®