

Patented quality Made in Germany

# Softskin Cushion

Get more information on  
[www.uniprox.de](http://www.uniprox.de)



Silicone Liner with soft,  
lateral tapering silicone padding.


unique prosthetic solutions  
A company of the Bauerfeind Group

 uniprox®

# SoftskinAir® Cushion 40 and Softskin Cushion 40



- Silicone Liner without attachment and without umbrella
- Soft external End-Cap
- Shore 40 medical grade silicone, white color
- Lengths: until size 26,5 appr. 370 mm. from size 28 appr. 400 mm, from size 36 appr. 430 mm



Order Nr.	Available sizes	Description	Article Nr.
SC40- <input type="text" value="Size"/> -Air	23,5 / 25 / 26,5 / 28 / 30 / 32 / 34 / 36 / 38	Air Liner	4 168 020 40 50 XXX
SC40- <input type="text" value="Size"/>	23,5 / 25 / 26,5 / 28 / 30 / 32 / 34 / 36 / 38	Non-Air Liner	4 168 010 40 50 XXX

Example: SoftskinAir® Cushion 40 in Size 26.5 → SC40-26,5-Air → 4 168 020 40 50 265

› The circumference is measured 4 cm above the distal end of the residual limb!

The soft lateral tapering silicone padding is, for example, very suitable for

- Extreme bony trans tibial amputees
- Knee disarticulations
- Trans Femoral amputations with severe scarf tissue



Patented quality Made in Germany

# Softskin Cushion

Get more information on  
[www.uniprox.de](http://www.uniprox.de)

**Uniprox GmbH & Co. KG**

Heinrich-Heine-Straße 4

07937 Zeulenroda

Germany

Phone +49 36628 66 33 00

Fax +49 36628 66 33 55

E-mail [info@uniprox.de](mailto:info@uniprox.de)

Customer-Service +49 36628 66 33 70

Fax +49 36628 66 33 77

A company of the Bauerfeind Group \_ [www.uniprox.de](http://www.uniprox.de)