

Quality and function


NEW

get more information on
www.uniprox.de

It takes two: SoftskinAir® Duo

- Soft on the inside and stable on the outside
- Medical silicone (Shore hardness 30/50)
- Soft padding of residual limb
- Efficient sweat transport through micropores
- Wall thickness 3 mm, length 370 mm
- Stainless steel pin holder
- Liner stand included
- Available in sizes 16 - 45

> Also available as a Softskin without Air function



SoftskinAir® Duo

Order No.	Size	Article No.
S-Duo-20-Air	16	4 164 020 35 00 160
S-Duo-22-Air	18	4 164 020 35 00 180
S-Duo-23,5-Air	20	4 164 020 35 00 200
S-Duo-25-Air	21	4 164 020 35 00 210
S-Duo-26,5-Air	22	4 164 020 35 00 220
S-Duo-28-Air	23,5	4 164 020 35 00 235
S-Duo-30-Air	25	4 164 020 35 00 250
S-Duo-30-Air	26,5	4 164 020 35 00 265
S-Duo-30-Air	28	4 164 020 35 00 280

Order No.	Size	Article No.
S-Duo-30-Air	30	4 164 020 35 00 300
S-Duo-32-Air	32	4 164 020 35 00 320
S-Duo-34-Air	34	4 164 020 35 00 340
S-Duo-36-Air	36	4 164 020 35 00 360
S-Duo-38-Air	38	4 164 020 35 00 380
S-Duo-40-Air	40	4 164 020 35 00 400
S-Duo-42-Air	42	4 164 020 35 00 420
S-Duo-45-Air	45	4 164 020 35 00 450

unique prosthetic solutions

A company of the Bauerfeind Group

 **uniprox**®